IKA

designed for scientists

DRY BLOCK HEATER 1 DRY BLOCK HEATER 2 DRY BLOCK HEATER 3 DRY BLOCK HEATER 4


Betriebsanleitung Ursprungssprache	DE	3
Operating instructions	EN	11
Mode d'emploi	FR	19
Инструкция по эксплуатации	RU	27
使用说明	ZH	35

Indicaciones de seguridad	ES	43
Veiligheidsinstructies	NL	44
Avvertenze per la sicurezza	IT	45
Säkerhetsanvisningar	SV	46
Sikkerhedshenvisninger	DA	47
Sikkerhetsanvisninger	NO	48
Turvallisuusohjeet	FI	49
Indicações de segurança	PT	50
Υποδείξεις ασφάλειας	EL	51
Wskazówki bezpieczeństwa	PL	52

Bezpečnostní pokyny	CS	53
Biztonsági utasítások	HU	54
Varnostna navodila	SL	55
Bezpečnostné pokyny	SK	56
Ohutusjuhised	ET	57
Drošības norādījumi	LV	58
Saugos nurodymai	LT	59
Указания за безопасност	BG	60
Instrucțiuni de siguranță	RO	61

**DRY BLOCK HEATER 1
DRY BLOCK HEATER 2
DRY BLOCK HEATER 3
DRY BLOCK HEATER 4**

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A	Netzbuchse
B	KT-Buchse
C	Stativgewindebohrung
D	Kontaktstecker
E	Geräteschalter
F	Taste "Start/Stop"
G	Taste Timer "+"
H	Taste Timer "-"
I	Taste Timer/ Counter
J	Taste Heizung "+"
K	Taste Heizung "-"
L	LED Display Heizung
M	LED (Set=Sollwert)
N	LED externer Temperatursensor
O	LED Heizplatte
P	LED Display Timer

Item Designation

A	Power socket
B	KT socket
C	Threaded hole for stand
D	Contact plug
E	Mains switch
F	"Start/Stop" button
G	Button timer "+"
H	Button timer "-"
I	Button Timer/Counter
J	Button heating "+"
K	Button heating "-"
L	LED Display heating
M	LED (set=set point value)
N	LED external temperature sensor
O	LED heating plate
P	LED display timer

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EN

Contents

	Page
EU Declaration of conformity	11
Explication of warning symbols	11
Safety instructions	12
Intended use	13
Unpacking	13
Assembling the support rod	13
Operation	14
Setting the operating mode	14
Timer / Counter	14
Temperature control	15
Maintenance and cleaning	16
Accessories	16
Error codes	17
Technical data	18
Warranty	18

EU Declaration of conformity

We declare under our sole responsibility that this product corresponds to the directives 2014/35/EU, 2014/30/EU and 2011/65/EU and conforms with the following standards or normative documents: EN 61010-1, EN 61010-2-010, EN 61326-1, EN 60529 and EN ISO 12100.

A copy of the complete EU Declaration of Conformity or further declarations of conformity can be requested at sales@ika.com.

Explication of warning symbols

Indicates an (extremely) hazardous situation, which, if not avoided, will result in death, serious injury.

Indicates a hazardous situation, which, if not avoided, can result in death, serious injury.

Indicates a potentially hazardous situation, which, if not avoided, can result in injury.

Indicates the exposure to a hot surface.

Safety instructions

/// General information

- › **Read the operating instructions completely before starting up and follow the safety instructions.**
- › Keep the operating instructions in a place where it can be accessed by everyone.
- › Ensure that only trained staff work with the device.
- › Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- › The device must only be used in a technically perfect condition.

Risk of burns! Exercise caution when touching parts of the housing and the heating plate.

The heating plate can reach dangerous temperatures. Pay attention to the residual heat on the heating plate after switching off the device.


The device may only be transported when the heating plate has cooled down.

/// Device setup

The power switch of the device must be accessed immediately, directly and without risk at any time. If access to the power switch cannot be ensured, an additional emergency stop switch that can be easily accessed must be installed in the work area.

- › Set up the device in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.
- › The feet of the device must be clean and undamaged.
- › Ensure that the power cord set / temperature sensor cable does not touch the heating plate.
- › Observe the minimum distances:
 - between devices min. 100 mm,
 - between device and wall min. 100 mm.

/// Working with the device

Do not use the device in explosive atmospheres, it is not EX-protected.

With substances capable of forming an explosive mixture, appropriate safety measures must be applied, e.g. working under a fume hood.

To avoid body injury and property damage, observe the relevant safety and accident prevention measures when processing hazardous materials.


- › Only process media that will not react dangerously to the extra energy produced through processing. This also applies to any extra energy produced in other ways, e.g. through light irradiation.
 - › The safety temperature must be set in accordance with EN 61010-2-010 Chapter "Requirements for devices containing or using flammable liquids".
 - The surface temperature of the flammable medium that is exposed to air may not exceed its flash point.
- A danger usually arises if a medium is heated in open vessels.

- The surface temperature of the heating device (e.g. the mounting plate) may not exceed the value of $(t - 25) \text{ }^\circ\text{C}$ (= set value of the safety circuit) on the surface of the flammable medium and in contact with air, whereby t is the fire point of the liquid.

A danger usually arises if a medium is heated in glass vessels (glass breakage).

If a setting made by the user (medium temperature or safety temperature) could bring a flammable medium into a state in which the conditions mentioned above could be exceeded, additional measures must be introduced that will protect the user from this danger.

- › When handling hazardous substances or mixtures of substances which are toxic or are contaminated with pathogenic micro-organisms, the user should take appropriate precautions. Observe the national regulations, the biosafety level of your laboratory, the material safety data sheets.
- › For germs or biological material or risk group II or higher, please observe the "Laboratory Biosafety Manual" issued by the World Health Organization).

Wear your personal protective equipment in accordance with the hazard category of the media to be processed. There may

be a risk from:

- splashing and evaporation of liquids,
- release of toxic or combustible gases.

- › Please pay attention to avoiding hand injury when operate the device.

The heating plate and the blocks can be very hot and cause injuries. Be careful when touching the blocks.

Allow the blocks attachments to cool down completely before removing them.

- › Always use the tool supplied with the device to remove the blocks from the device.
- › Do not heat liquid media except with the block. Never pour liquid media directly on to the heating plate!

/// Accessories

- › Protect the device and accessories from bumps and impacts.
- › Check the device and accessories for damage before each use. Do not use damaged components.
- › Safe operation is guaranteed only with the use of original IKA accessories.
- › Ensure that the external temperature sensor is inserted into the medium to a depth of at least 20 mm when connected.
- › Always disconnect the plug before attaching accessories.
- › Accessories must be securely attached to the device and cannot come off by themselves. The centre of gravity of the assembly must lie within the surface on which it is set up.
- › Observe the operating instructions of the accessories.

/// Power supply / Switching off the device

The device will automatically restart in mode B following any interruption to the power supply.

- › The voltage stated on the type plate must correspond to the power voltage.
- › The device can only be disconnected from the power supply by pulling out the power plug or the connector plug.
- › The device must only be operated with the original power cord set.
- › The socket for the power cord must be easily accessible.
- › Socket must be earthed (protective ground contact).

/// Maintenance

- › The device must only be opened by trained specialists, even during repair. The device must be unplugged from the power supply before opening. Live parts inside the device may still be live for some time after unplugging from the power supply.
- › Only use original IKA spare parts!

/// Disposal instructions

- › The device, accessories and packaging must be disposed of in accordance with local and national regulations.

Intended use

/// Use

- For heating media in block heaters.

/// Area of use

Indoor environments similar to that a laboratory of research, teaching, trade or industry area.

The safety of the user cannot be guaranteed:

- if the device is operated with accessories that are not supplied or recommended by the manufacturer,
- if the device is operated improperly or contrary to the manufacture's specifications,
- if the device or the printed circuit board are modified by third parties.

Unpacking

/// Unpacking



- Please unpack the device carefully
- In the case of any damage a fact report must be sent immediately (post, rail or forwarder)

/// Scope of delivery


- IKA DRY BLOCK HEATER
- Depending on the IKA DRY BLOCK HEATER model 1, 2, 3, or 4
- Mains cable
- Temperature sensor PT 1000.60
- Operating instructions
- Tool for block removal

Assembling the support rod

- › Review the mounting and safety instructions of the IKA boss head clamp prior to using it.
- › The device must not be suspended from the support rod!

Risk of tipping! Make sure that the center of mass of the attached device does not protrude over the safe area indicated below with a dotted rectangle.

Operation

Commissioning	☞ Move device switch (E) to the OFF position
	☞ Insert the mains power cable into the power socket (A)
	☞ Move device switch (E) to the ON position
	⇒ The unit will be set to factory setting mode A (see "Setting operating mode")
Heating	☞ Set the target temperature using the heater buttons +/- (J) and (K)
	⇒ The value selected will be shown on the display (L)
	☞ Start the heating function by pressing the "Start / Stop" button (F)
	The set and actual temperatures are shown alternately on the display (L) via blinking LEDs (M) :
	(i) ⇒ When the heating is on, the LED (O) lights up.
(i) ⇒ During standby operation, the display (L) will show [HOT] if the base plate temperature is above 50 °C.	
Connecting external temperature sensor	☞ Move device switch (E) to the OFF position
	☞ Detach contact plug (D)
	☞ Attach PT 1000.60 temperature sensor to the socket (B)
	☞ Move device switch (E) to the ON position
	(i) ⇒ The actual temperature for the temperature sensor shown on display (L) will correspond to the temperature of the media. The LED "external temperature sensor" PT 1000 (N) will be lit.

Setting operating mode

The mode selected will be shown on the display when the device is started up.

/// Mode A

All settings will be stored if the device is switched off or disconnected from the power supply. The heating functions will be set to OFF when the device is powered on.

⇒ Please be noted that, during power on / Power failure, the heating function will be OFF.

⇒ Factory setting: Mode A

/// Mode B

All settings will be stored if the device is switched off or disconnected from the power supply. The heating function will be set to ON or OFF when the device is powered on, depending on the previous status of the device.

/// Changing Operating Mode

- ☞ Move device switch **(E)** to the OFF position
- ☞ Keep the selection button **(F)** depressed
- ☞ Move device switch **(E)** to the ON position
- ☞ Release the selection button **(F)**

(i) ⇒ The set mode is indicated on the Display **(P)** **[b]**

Timer / Counter

/// Selection of Timer / Counter

After the device has been switched on, the "Timer/ Counter" button **(I)** can be used to toggle between timer mode and counter mode.

⇒ The mode can be activated only when the heating function is off.

/// Timer

The Timer is selected by pressing Timer/ Counter button **(I)**. Pressing the buttons "+" **(G)** and "-" **(H)** allows selection of a time period of up to 99 hours and 59 minutes. Once the heating function **(F)** is activated, the pre-set time starts to run. When the total time has run down to 00 hours and 00 minutes, a warning buzzer sounds for 30 seconds and in addition the display flashes. The heating function is then deactivated.

/// Counter

When the heating function is started by pressing the button **(F)** the counter is activated. In counter mode the time is counted upwards from 00 hours 00 minutes until the heating function is stopped. The counter always starts from 00:00 hh:mm.

⇒ Once the counter reaches to 99:59, the counter restarts from 00:00. The heating function will not be stopped.

Temperature control

/// Recommended procedure

- › Always use the closed vessel to avoid evaporation of the samples. This procedure helps to achieve the highest possible temperature accuracy within the sample.
- › It is recommended to select proper heating blocks to ensure best results. The fitment of vials / tubes plays an important role for best results. There should not be any air gap between the tube cavity in the block and the tubes.
- › Make sure that external temperature probes if used, are immersed properly inside the tubes.
- › Please note the characteristics of the sample liquid used while making temperature adjustments.
- › Note that the liquid level inside the vials / tubes is lower than the height of the heating block.
- › Ensure proper contact between the heating plate and modular block for good thermal conductivity.
- › Blocks and heating plate surface must be free from any scratches / deformation.

/// Temperature Control Without External Temperature Sensor

When the KT socket **(B)** is connected to the contact plug **(D)**, the device switches its temperature control to the internal sensor control mode, in which the surface temperature of the heating plate is the set temperature, irrespective of the temperature of the medium or the block.

- › The final temperature inside the medium may vary from the set temperature.
- › The device must not be put into operation without the contact plug **(D)** /external temperature probe in the KT socket **(B)**.

/// Temperature Control With External Temperature Sensor

The preferable method for controlling the block temperature is with the external temperature sensor.

Insert external temperature sensor PT 1000.60 supplied along with the device to the socket **(B)**. Ensure external temperature sensor LED **(N)** is lit up. With this set up, the device will shift its control to external sensor control where the temperature on the sensor is considered for target temperature. The sensor probe can be inserted into the modular block (Same position on the block for block removal tool and location for external temperature) or into the medium directly. Below points need to be checked while using external sensor.

- › Probe must be inserted into the medium / block properly. Failure which can result in maximum heating of block.
- › Always use closed tubes for temperature accuracy.
- › If the sensor is immersed in the sample always note the characteristics of the sample example – evaporation rate, viscosity etc.
- › Ensure the tip of the sensor reaches the bottom of the tube or vial.

/// Temperature offset adjustment

In some test cases, due to the sample characteristics or the heat losses due to the vessel / block temperature variations may be observed in the sample with respect to the set temperature when the external probe is placed in the sample / Block. To measure this, insert an external probe in the sample and observe the difference in the achieved temperature in the sample. In this case, IKA Dry Block Heater has an option to introduce an offset to the set temperature. The procedure for offset adjustment is as follows.

- › Move device switch **(E)** to the OFF position
 - › Keep the Timer "-" button **(H)** depressed
 - › Move device switch **(E)** to the ON position
 - › Release the Timer "-" button **(H)**
 - › Now the display **(L)** will show "oF" and Display **(P)** will show the offset value
 - › The offset value can be changed from -5 °C (Set Temperature - 5 °C) to +5 °C (Set temperature +5 °C) using "Temp+" key **(J)** and "Temp-" **(K)** key. (Set Temperature ≤ 115 °C)
 - › Once the offset is registered with the device the value is stored in the memory until it is changed by the user
- Since the heat transfer coefficient between block and Vessels can differ (depending on size tolerances of the used vessels) we recommend to check the Temperature in more than one vessel.

/// Temperature uniformity test

Uniformity is tested in the factory by using IKA AB3 Test Block to ensure even heating with the following procedure.

- › Set up the unit in a stable ambient.
 - › The device external temperature sensor (PT 1000) is placed in the middle slot of the test Block as shown.
 - › Select multiple calibrated temperature measuring probes in other slots.
 - › Set the temperature of the unit.
 - › Allow the unit to reach the set temperature and stabilize the set temperature for 30 minutes.
 - › Take the temperature readings to calculate the uniformity.
- By using the same setup it is possible to measure stability and homogeneity.

Maintenance and cleaning

The device is maintenance-free. It is only subject to the natural wear and tear of components and their statistical failure rate.

/// Cleaning

- For cleaning disconnect the mains plug!
- Use only cleaning agents which have been approved by IKA to clean the devices:
These are: water (containing surfactant) and isopropyl alcohol.
- Wear protective gloves during cleaning the devices.
- Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.
- Do not allow moisture to get into the device when cleaning.
- Before using another than the recommended method for cleaning or decontamination, the user must ascertain with that this method does not destroy the device.

/// Ordering spare parts

When ordering spare parts, please indicate:

- device type.
- serial number, see type plate.
- Item and designation of the spare part, see: www.ika.com, spare parts diagram and spareparts list.
- software version.

/// Repair

Please only send devices in for repair that have been cleaned and are free of materials which might present health hazards.

For repair, please request the “**Safety Declaration (Decontamination Certificate)**” from IKA or use the download printout of it from IKA website at www.ika.com.

If your appliance requires repair, return it in its original packaging. Storage packaging is not sufficient when sending the device - also use appropriate transport packaging.

Accessories

Block Name	Suitable For	Number of Wells
DB 1.1	PCR tubes (0,5 ml)	30
DB 1.2	Eppendorf microtubes (1.5 ml)	20
DB 1.3	Microtubes (1,5 ml)	20
DB 1.4	Eppendorf microtubes (2,0 ml)	20
DB 1.5	Corning microtubes (2,0 ml)	20
DB 2.1	Conical tubes (15 ml)	12
DB 2.2	Conical tubes (50 ml)	5
DB 3.1	Test tube combination, 6 / 12 / 13 / 25 mm	6 / 5 / 3
DB 3.2	Centrifuge tube combination, 1,5 / 15 / 50 ml	4 / 3 / 2
DB 3.3	Microtube combination, 0,5 / 1,5 / 2,0 ml	6 / 10 / 5
DB 4.1	Round bottom tubes (6 mm)	30
DB 4.2	Round bottom tubes (10 mm)	24
DB 4.3	Round bottom tubes (12/13 mm)	16
DB 4.4	Round bottom tubes (12/13 mm)	20
DB 4.5	Round bottom tubes (15/16 mm)	12
DB 4.6	Round bottom tubes (20 mm)	8
DB 4.7	Round bottom tubes (25 mm)	6

Block Name	Suitable For	Number of Wells
DB 4.8	Round bottom tubes (35 mm)	4
DB 4.9	Round bottom tubes (17/18 mm)	12
DB 5.1	12 mm vials	30
DB 5.2	15 mm vials	20
DB 5.3	17 mm vials	12
DB 5.4	19 mm vials	12
DB 5.5	21 mm vials	9
DB 5.6	23 mm vials	8
DB 5.7	25 mm vials	8
DB 5.8	28 mm vials	6
DB 5.9	16 mm vials	15
DB 6.1	10 x 8 individual PCR tube strips (0,2 ml tubes)	80
DB 6.2	64 x individual PCR tubes (0,2 ml tubes)	64
DB 6.3	1 x 96 – well PCR plate (0,2 ml tubes)	96
DB 7.1	96 – or 384 – well plate	-
DB 8.1	2 x 6 cuvettes (12,5 mm)	-

› For accessories see www.ika.com.

Error codes

Any malfunctions during operation will be identified by an error message on the display (LED **(P)**).

Proceed as follows in such cases:

- Switch off device using the main switch **(E)**
- Carry out corrective measures
- Restart device

Error code	Cause	Effect	Solution
Error 03	Device internal temperature is too high	Heating function off	- Switch off the device and allow it to cool down. - Restart the device.
Error 25	Defective Heating / Switching element	Heating function off	- Switch off the device and allow it to cool down.
Error 26	Difference between temperature reading from safety sensor and control sensor is greater than 40K	Heating function off	- Switch off the device and allow it to cool down.

If the actions described fail to resolve the fault or another error code is displayed then take one of the following steps:

- contact the **IKA** service department,
- send the device for repair, including a short description of the fault.

Technical Data

		IKA DRY BLOCK HEATER 1	IKA DRY BLOCK HEATER 2	IKA DRY BLOCK HEATER 3	IKA DRY BLOCK HEATER 4
Operating voltage range	Vac	(220 ... 230) ± 10%			
	Vac	115 ± 10%			
Rated voltage	Vac	100 ± 10%			
	Vac	230 / 50 Hz			
	Vac	115 / 60 Hz			
	Vac	100 / 60 Hz			
Frequency	Hz	50 / 60			
Display		digital			
Temperature range	°C	ambient temperature +5 ... 120			
Timer	min/h	1 min ... 99 h 59 min			
Perm. duration of operation	%	100			
Perm. ambient temperature	°C	+5 ... +40			
Perm. humidity	%	80			
Protection acc. to DIN EN 60529		IP 21			
Protection class		I			
Overvoltage category		II			
Contamination level		2			
Operation at a terrestrial altitude	m	max. 2000			
Dimension (W x D x H)	mm	151 x 228 x 73	151 x 304 x 73	151 x 380 x 73	151 x 456 x 73
Weight	kg	1.3	1.7	2.0	2.0
Heating					
Power consumption (+10%) max. at 230 Vac	W	165	250	330	412
115 Vac		165	250	330	412
100 Vac		165	250	330	412
Temperature sensor PT 1000 variation DIN EN 60751 Kl. A	K	≤ ± 0.1			
Temperature stability within the blocks at 37 °C	°C	± 0.2	± 0.2	± 0.2	± 0.2
Temperature stability within the blocks at 60 °C	°C	± 0.4	± 0.4	± 0.4	± 0.4
Temperature homogeneity within the blocks at 37 °C	K	± 0.2	± 0.2	± 0.3	± 0.3
Temperature homogeneity within the blocks at 60 °C	K	± 0.4	± 0.4	± 0.6	± 0.6
Heating rate (PT 1000 in solid block)	K/min	5.0	4.5	4.0	4.0
Fixed safety circuit					
Safety temperature limit	°C	150			

Test parameters: with test blocks

Subject to technical changes!

Warranty

In accordance with **IKA** warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual.

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IKA India Private Limited
IKA England LTD.
IKA Vietnam Company Limited
IKA Works (Thailand) Co. Ltd.