STICKMAX II

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DE

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
Safety instructions

EN

- › Always wear personal protective equipment in accordance with the hazard class of the media being worked with.
- › The surface of the adhesive pad must be free of grease, dirt, moisture, and oil.
- › Check all vessels placed on the pad for damage. Never use damaged vessels.


Intended use

- › The STICKMAX II adhesive pad, with a diameter of 82 mm, makes the use of fixing clips unnecessary in many situations. It is designed for holding vessels (e.g. Erlenmeyer flasks, Petri dishes). The vessels used must have a flat, even base.

Assembly

- › Degrease the MS 3.3 Universal attachment (e.g. with isopropanol).
- › Remove the protective foil on the bottom side.
- › Place the adhesive pad on the universal attachment with the "STICKMAX II" lettering facing downwards. Press them lightly to increase the adhesive effect.
- › Then remove the protective foil on the top side.
- › Press the vessel lightly to increase the adhesive effect.

Use, handling and cleaning

Permissible ambient temperature	5 ... 80 °C
Maximum shaker speed	500 rpm
Maximum load	0.5 kg

- › Avoid strong or supplemental UV light.
- › Clean the adhesive pad of the attachment with a surfactant-containing water moistened, lint-free cloth.
- › Use, for example, 70 percent ethanol to disinfect.
- › To detach the set-up vessel from the adhesive base, press the upper edge of the vessel sideways in any direction (**Fig. 1**). While doing so, fix the universal attachment with your other hand.
- › The adhesive pad can be removed from the shaker table or the universal attachment by pulling it off sideways at an angle.
- › The adhesive effect of the adhesive base can be restored by cleaning with water.