# IKA

designed for scientists

**RC 2 lite**

ESPAÑOL


Fig. 1

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IKA Japan K.K.  
Phone: +81 6 6730 6781  
eMail: info\_japan@ika.ne.jp

#### **INDIA**

IKA India Private Limited  
Phone: +91 80 26253 900  
eMail: info@ika.in

#### **UNITED KINGDOM**

IKA England LTD.  
Phone: +44 1865 986 162  
eMail: sales.England@ika.com

#### **VIETNAM**

IKA Vietnam Company Limited  
Phone: +84 28 38202142  
eMail: sales.lab-vietnam@ika.com

#### **THAILAND**

IKA Works (Thailand) Co. Ltd.  
Phone: +66 2059 4690  
eMail: sales.lab-thailand@ika.com

#### **TURKEY**

IKA Turkey A.Ş.  
Phone: +90 216 394 43 43  
eMail: sales.turkey@ika.com

---

Discover and order the fascinating products of IKA online:  
[www.ika.com](http://www.ika.com)

# IKA

designed for scientists

## IKA Vacstar digital

Operating instructions  
Source language: German

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Fig. 1

### Pos. Description (EN)

A	Handle
B	Display
C	Rotating/pressing knob
D	Main switch (on/off)

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Fig. 1

### Pos. Description (EN)

E	Handle securing screw
F	Vacuum controller VC 10 stand hole
G	Stand securing screw
H	Fan / ventilation slit
I	Hose connection for suction line d= 8 mm (INLET)
J	Accessory fastening (emission condenser)
K	Connection for vacuum controller VC 10 (Mini DIN)
L	USB interface
M	RS 232 interface
N	Water valve connection
O	Hose connection for pressure line d= 8 mm (OUTLET)
P	Power supply cable connection

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	Page
Device setup .....	2/3
Declaration of conformity .....	5
Explication of warning symbols .....	5
Warranty .....	5
Safety instructions.....	6
Unpacking .....	7
Correct use .....	8
Useful information .....	8
Setting up.....	9
Commissioning .....	12
Interfaces and outputs .....	13
Maintenance and cleaning.....	14
Error codes .....	15
Accessories .....	15
Product contact parts.....	15
Technical data.....	16

## Declaration of conformity

We declare under our sole responsibility that the product to which this declaration relates is in conformity with directives 2014/35/EU, 2006/42/EC, 2014/30/EU and 2011/65/EU and conforms with the following standards or normative documents: EN 61010-1, EN 60529, EN 61326-1 and EN ISO 12100.

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## Explication of warning symbols

 <b>DANGER</b>	Indicates an (extremely) hazardous situation, which, if not avoided, will result in death, serious injury.
 <b>WARNING</b>	Indicates a potentially hazardous situation, which, if not avoided, can result in death, serious injury.
 <b>CAUTION</b>	Indicates a potentially hazardous situation, which, if not avoided, can result in injury.
 <b>NOTE</b>	Indicates practices which, if not avoided, can result in equipment damage.


## Warranty

According to **IKA**'s Terms and Conditions of sale and delivery, this product is covered by a warranty for a period of 24 months. Please contact your dealer for any warranty claims. If you wish, you can send the device directly to our factory. Please provide the sales invoice and state the reasons for your guarantee claim. In this case, you are responsible for shipping and handling costs.

The warranty does not cover wearing parts, nor defects that are the result of improper use, insufficient care and maintenance or failure to follow the instructions in this operating manual.

## Safety instructions

### For your protection

**Please read the instruction manual in full before use and follow the safety instructions.**

- Keep the instruction manual in a place where it can be accessed easily.
- Ensure that only trained staff use the device.
- Be sure to comply with all safety instructions, directives and all matters of health, safety and accident prevention in the workplace.
- The device and all parts of the device must not be used on people or animals.

**Always wear personal protective equipment in accordance with the hazard class of the media being worked with - otherwise there are dangers due to spraying liquids, and release of toxic or flammable gases.**

- Do not expose human or animal body parts to vacuum.
- Do not work with the device underwater or underground.

### Device configuration

**The vacuum pump IKA Vacstar digital is not designed to be set up in potentially explosive atmospheres.**

- Set up the device in accordance with the chapter "Setting up" and connect the connection lines and interfaces as described.
- Set up the device on a stable, even, non-flammable surface.

- Never work with a faulty or incorrectly connected device.
- Set up the vacuum pumps in appropriate and functioning extractor hoods, or set up an exhaust gas line to them. Ensure that the exhaust gas line cannot be kinked! The maximum permitted exhaust gas line length is 2 metres.
- Occurrence of explosive mixtures; if appropriate connect inert gas for ventilation and/or thinning.

### Accessories

- Safe operation can only be ensured when working with accessories as described in the "Accessories" section.
- Carefully observe the operating instructions for additional devices (e.g. rotary evaporators, vacuum controllers), with which the vacuum pump IKA Vacstar digital is operated.
- The pressure at the gas inlet and outlet must not exceed 1100 mbar.

- Elastic elements may be pressed together under vacuum.
- Only use flexible hose lines.
- Observe your emergency measures for power failure and ensure that the device is put in a safe state (see the chapter: Commissioning, Operating mode).

### Work with the device

**You can prevent the release of solvent vapours into the atmosphere using a downstream emission condenser.**

**The vacuum pump is not suitable for use with self-igniting materials, materials that are flammable without an air supply, or explosive materials.**

**Do not operate the pump while it is open. Otherwise there is a danger of severe injuries as a hand may be pulled into the moving parts.**

**Inhaling or coming into contact with media such as poisonous liquids, gases, spray mists, vapours, dusts or biological materials can endanger the health of the user. Ensure that all connections are well sealed and free of leaks if you are working with such media.**

- The vacuum pump IKA Vacstar digital must only be operated under the conditions described in the chapter "Technical data".
- Prevent release of the materials listed above. Take measures to protect staff and the environment.
- Pay attention to possible interactions or chemical or physical

reactions when working with media at reduced pressure and increased temperature.

- There can be electrostatic processes between the medium and the device which can lead to direct danger.
- Some medium may be released due to the residual leakage rate of the device.
- Before commissioning check that all the housing parts are present and fastened to the device.
- Do not lift the pump if the handle is loose and the handle securing screw (E) is loose or missing!
- Only switch the pump on if the pump is standing vertically.
- Connect the hose connections (INLET-OUTLET) and interfaces in accordance with the labelling on the device and the operating instructions.
- Please note that the temperature of the evacuated medium is below its ignition temperature. The pumping process (compression) increases the temperature of the medium additionally.
- Ensure that vapours containing solvent can be sucked into the pump.
- Do not use the pump to create pressure.
- Observe the permitted pressure at the inlet and outlet side; see the chapter "Technical data".
- The gas flow must only be regulated/throttled in the suction-side line.
- Use a solenoid valve or a check valve in the suction line if there are several load units.
- When using an emission condenser ensure that the coolant is flowing freely.

### Protection of the device

**Ensure that the substances used with the device are compatible with the materials of the device parts that come into contact with the product; see the chapter "Product contact parts".**

**The specified settings on the rating plate must coincide with the actual power supply.**

- The device is only disconnected from the power supply when the power or device plug is removed.
- The socket for the power supply cable must be easy to access.
- Covers or parts that can be removed from the device without tools must later be refitted to ensure safe operation, unless anything else is connected at this point. This will prevent the infiltration of foreign objects, liquids and other contaminants.
- Unwanted liquids that have entered can be removed by extracting air during no-load operation.
- Avoid knocks or impact to the device.

- The device must only be opened by trained, skilled personnel (who have been authorised to do so).
- The device must not be modified without authorisation from IKA.
- To ensure sufficient cooling for the vacuum pump IKA Vacstar digital, the ventilation slits on the housing must not be covered.
- Only use original replacement parts for maintenance to ensure that the device works reliably.
- Watch out for water condensation inside and outside the device. Warm up the device first, if it has been brought in from a cold environment.
- Never fasten the vacuum pump over a heating bath.
- Ensure that no solids and/or liquids can enter the diaphragm of the vacuum pump IKA Vacstar digital through the suction line of the pump. This damages the diaphragm and other internal parts of the pump. This reduces delivery capacity and the final pressure can no longer be achieved. Deposits may be produced in the interior, reducing service life and causing leaks.

## Unpacking

### Unpacking

- Unpack the device carefully.
- Any damage should be notified immediately to the shipping agent (post office, railway network or transport company).

### Scope of delivery

Diaphragm vacuum pump IKA Vacstar digital	
Power supply cable Laptop Euro	
Power supply cable Laptop GB	
Power supply cable Laptop CH	
USB A - USB B cable	
Vacuum hose 1m	

Silencer	
Connecting cable	
Operating instructions	
Warranty card	

## Correct use

### Use

#### WARNING

**The device must not be used for:**  
 - **evacuating biological habitats,**  
 - **evacuating explosive, corrosive or similar gases,**  
 - **throughput/use of liquids**

**The device is not suitable for:**  
 - **pumping up objects**  
 - **pumping liquids**  
 - **use in potentially explosive atmospheres**  
 - **underwater**  
 - **underground**

#### Standalone operation

The vacuum pump can be operated without a controller. The delivery capacity can be set using the adjustable speed. Generally speaking, the smaller the volume to be evacuated, the lower the speed should be. The final pressure is not controlled.

#### Operation with accessories

Along with the accessories recommended by **IKA**, the vacuum pump **IKA** Vacstar digital is suitable for controlled evacuation to a specified final pressure e.g. for operating an **IKA** rotary evaporator. It can also be used for classical separation, filtration or drying in the laboratory.

With an appropriate vacuum controller (e.g. VC 10) or with a rotary evaporator (e.g. RV 10 auto), the Vacstar digital can be operated automatically, for automatic boiling point detection, scheduled pressure-time curves or programs from the solvent library, for example.

Intended use: Tabletop device

The laboratory diaphragm vacuum pump **IKA** Vacstar digital is used together with other laboratory devices to create a vacuum.

For specific vacuum control, an additional vacuum controller **IKA** VC 10 or **IKA** rotary evaporator RV 10 auto is required.

### Area of use (only indoors)

The device is suitable for use in residential areas and all other areas.

The safety of the user cannot be ensured:

- if the device is used in conjunctions with accessories not made or recommended by the manufacturer!

- if modifications are made to the device or parts of the device by third parties!
- if the device is operated in contravention of the safety instructions!

## Useful information

#### DANGER

**You can prevent the release of solvent vapours into the atmosphere using a downstream emission condenser.**

#### CAUTION

**Special care must be taken when working with gas mixtures that include condensable gases (e.g. solvents). If diaphragm pumps are used with such gases they condense in the pump's air flow. This leads to a build-up of pressure in the vacuum chamber, which in turn influences the effectiveness and service life of the diaphragms and valves.**

**If liquid or condensate enters the vacuum chambers, the specified minimum working pressure is no longer achieved.**


**To protect the inner valves and diaphragms against liquid, use an upstream condensate separator, such as a Woulff bottle.**

#### CAUTION

**It is generally recommended that you run the pump with no load at the end of the process (no-load operation, approx. 3 up to max. 5 minutes) to dispel any residues in the pump.**

**If applicable, vent the recipient. To do this, remove the suction line during manual operation.**

#### NOTE

**Observe the operating instructions of the relevant devices.**

Diaphragm vacuum pumps are oscillating displacement pumps. The diaphragms expand the volume of the vacuum chamber to produce an underpressure, which sucks air into the chamber. The air is then pressed out of the vacuum chamber by subsequent compression. Changes between the inlet and outlet are controlled by diaphragm valves.

The physical unit for pressure is 1 Pascal [Pa].

However people often still use the unit millibar [mbar] colloquially.

1 mbar = 100 Pa  
 1 bar = 10<sup>5</sup> Pa  
 1 Pa = 10<sup>-5</sup> bar  
 1 Pa = 0.01 mbar

A distinction is made between different ranges in vacuum technology.

Rough vacuum:	10 <sup>5</sup> to 10 <sup>2</sup> Pa	(1000 to 1 mbar)
Fine vacuum:	10 <sup>2</sup> to 10 <sup>-1</sup> Pa	(1 to 10 <sup>-3</sup> mbar)
High vacuum:	10 <sup>-1</sup> to 10 <sup>-5</sup> Pa	(10 <sup>-3</sup> to 10 <sup>-7</sup> mbar)
Ultra high vacuum	< 10 <sup>-5</sup> Pa	(< 10 <sup>-7</sup> mbar)

The periodic movements of the diaphragms create a vacuum at the pump inlet. Diaphragm vacuum pumps work with a rough vacuum of up to 10<sup>2</sup> Pa. Lower pressures cannot be achieved with these pumps.

To create fine, high and ultra high vacuums, oil-sealed vacuum pumps, screw pumps, diffusion pumps or turbomolecular pumps must be used, for example. Diaphragm vacuum pumps are often used as pumps for the preliminary stage in these contexts.

Compressing and releasing the pressure on air creates a noise. The pump is quieter the less air is transported. This means that there is less operating noise during normal operation if the pressure is reduced.

A fitted exhaust hose with a sound absorber can also be used to reduce noise.


In automatic operation with controlled speed the pump reduces its speed when it reaches the target value. If the leakage rate in the entire system is low enough, the speed may even be set to "zero".

When combined with the vacuum controller VC 10 or the rotary evaporator with integrated vacuum controller RV 10 auto, the pump can also be used for "two-position control" (valve-controlled) or "analogue speed control" of the vacuum.

### Two-position control

When the target value is reached, the suction line is interrupted. Due to the small delay between the times when the pressure value is detected, when comparison is made with the target value, when a vacuum valve is switched off and when the suction line of the vacuum-producing system are switched off, it is easy to go below the set target value. If the pressure in the system goes up again due to a natural leakage rate, the suction line valve opens.


Schematic view of two-position vacuum control

### Analogue speed control

With this control, the speed of the pump, and therefore its suction power, is reduced the closer the measured pressure reaches the target value. Once the target value is reached, the pump only operates according to the leakage rate of the system. This enables quiet functioning, precise vacuum control and a long service life for the moving pump parts. Exact automatic boiling point recognition is possible with this type of control, meaning that the system reaches and keeps the solvent at its boiling point in automatic mode.

Schematic view of speed vacuum control

## Setting up

#### NOTE

**Ensure that you follow the general rule, that the recipient (load/vacuum container/glass condenser) should always be connected to the suction line at its highest point. This reduces the risk of liquids entering the vacuum pump.**

**Liquid in the pump chambers makes the pump properties worse.**

**Connect the hose connections (INLET/OUTLET) and interfaces in accordance with the labelling on the pump and the operating instructions.**

The connections are barbed hose connectors. Ensure that the hoses have hose clips if applicable. Lay the hoses in such a way that no condensate can flow into the pump.

Ensure that the air can enter the ventilation slit (H) unhindered.

Install a condensate separator (e.g. a Woulff bottle) in the suction line before the intake spigot of the pump to protect against water coming in.

If solvent vapours are sucked in, an emission condenser (accessory) after the pump helps to condense it, and prevents it being released into the atmosphere.

## Connection of interfaces

### NOTE

Observe the relevant connections (see Fig. 1).

O: Hose connection d= 8 mm OUTLET

Link this connection to the emission condenser with a vacuum hose, or fit a sound absorber to the end.

### NOTE

Put the end of the hose in the extractor hood!

Check there is a free outlet on the pressurised side!

Do not use a throttle on the pressurised side and do not close the outlet! Connect the exhaust line to this connection.

K: Connection for vacuum controller VC 10 (Mini DIN)

You can connect the vacuum pump and the vacuum controller VC 10 or the rotary evaporator RV10 auto with the analogue connection cable for precise speed-controlled vacuum control. The vacuum controller detects the pump and switches to speed-vacuum control mode. Two-position control is deactivated.

L: USB interface

Connect the vacuum pump IKA Vacstar digital to a PC with a USB A - USB B cable. Any device software updates can be loaded using the IKA FUT software tool.

M: RS 232 interface

You can connect the IKA Vacstar digital vacuum pump to a PC using an RS 232 interface cable. The pump can be operated in conjunction with other devices with labworldsoft® laboratory device software. For further information see the chapter "Interfaces and outputs".

N: Water valve connection

Connect the optional water throttle valve RV 10.5001 to the diaphragm vacuum pump. The water flow to the emission condenser is controlled using the water valve. The valve is opened as soon as the pump is switched on.

I: Hose connection for suction line d= 8 mm INLET

Connect the suction line to this connection.

Link this connection to the recipient (rotary evaporator condenser, laboratory reactor etc.) with a vacuum hose.

P: Power supply cable connection

Check that the voltage information on the rating plate matches your mains supply. Connect the power supply cable to supply electricity.

## IKA Vacstar digital with vacuum controller VC 10

### NOTE

Cannot be operated with optional water valve RV 10.5001 in conjunction with VC 2.4.


Refer to the operating instructions for the vacuum controller VC 10!

When using the optional pump cut-out unit VC 2.4 (ref. 6) the pump is switched off when the target value is reached (operating mode B).

This device configuration is also possible without the pump cut-out unit VC 2.4. The pump runs on. For good control, the pump should be operated at a low speed.


When the target value is reached, the suction line is closed by the vacuum valve integrated in the VC 10. The pump speed can be changed by manually setting the speed.

- 1 Recipient (load, e.g. rotary evaporator, reactor)
- 2 Vacuum separator (e.g. Woulff bottle)
- 3 Vent valve
- 4 Pressure sensor
- 5 Vacuum valve/ball valve
- 6 Pump cut-out unit VC 2.4
- 7 Emission condenser

## Speed regulation

IKA Vacstar digital with vacuum controller VC 10 or rotary evaporator RV 10 auto


### NOTE

Observe the operating instructions for the devices.

Automatic setting of the operating mode "Speed regulation" as soon as the analogue connecting cable (ref. 8) connects the pump with the vacuum controller or rotary evaporator.

When the target value is reached, the pump speed is automatically reduced to a minimum of "0 rpm". Once the target value is reached, the pump applies suction according to the leakage rate of the vacuum system. Automatic boiling point recognition is possible.

- 1 Recipient (load, e.g. rotary evaporator, reactor)
- 2 Vacuum separator (e.g. Woulff bottle)
- 7 Emission condenser
- 8 Analogue connecting cable

## Commissioning

**NOTE** Observe the operating instructions for the devices.

Check if the voltage specified on the type plate matches the available mains voltage. If these conditions are met, the device is ready for operation when it is plugged in to the mains.

 The socket used must be earthed (grounded contact).

Otherwise safe operation is not guaranteed or the device may be damaged. Observe the ambient conditions indicated under "Technical data".

### Switching on

- Switch on the device using the main switch (D). All display segments are shown to enable a visual check.

8.8.8.8.

- Software version number display (e.g. vers. 1.00).

0.03

- Operating mode display (A; B). (Factory setting: Operating mode A)

A

The operating mode can be changed when you switch on the device.

#### Operating mode A:

In this operating mode, the set target value is not saved when the current run comes to an end or the device is switched off. No automatic restart after mains power is cut off.

#### Operating mode B:

In this operating mode, the set target value is saved when the current run comes to an end or the device is switched off, and the value can be changed.

Operating mode B is required particularly in combination with the pump cut-out unit VC 2.4.

Automatic restart after mains power is cut off.

#### Changing operating mode:

- Press the rotating/pressing knob (C) and switch on the device at the main switch (D) at the same time to change between the two operating modes.

A

b

In standalone operation or two-position operation in combination with the vacuum controller VC 10 the target speed can be set by turning the rotating/pressing knob (C). Push to start the pump and push again to stop it. A flashing light on the display (B) indicates when the device is running.

### Switching on with several interfaces

If the pump is operated via the analogue connecting cable with the vacuum controller VC 10 or the rotary evaporator RV 10 auto, "VC 10" is shown on the display.

VC 10

If the vacuum controller is connected to a PC, "PC" is shown on the display (B).

PC

### Switching off

Switch off the device using the main switch (D).

 Unplug the appliance from the mains socket to disconnect it from the mains.

## Interfaces and outputs

**NOTE** Please comply with the system requirements together with the operating instructions and help section included with the software.

The device can be operated in "Remote" mode via an RS 232 or USB interface using the laboratory software labworldsoft®. The RS 232 interface (M) at the back of the device is fitted with a 9-pole SUB-D jack which can be connected to a PC. The pins are assigned serial signals.

### USB Interface

The Universal Serial Bus (USB) is a serial bus system which allows the device to be connected to the PC. Devices that support USB can be connected to each other whilst they are running (hot plugging) and provide automatic recognition of the connected devices and their properties.

Use the USB interface in conjunction with labworldsoft® for operation in "Remote" mode and for updating the firmware using the "Firmware update tool".

### Installation

Before the device is connected with the PC using the USB data cable, the USB driver must be installed.

The USB driver can be downloaded from the website:

### Serial interface RS 232 (V 24)

Configuration:

- The functions of the interface circuit between the device and the automation system are a selection from the signals specified in the EIA standard RS232 as per DIN 66020 Part 1.
- Standard RS 232, corresponding to DIN 66259 Part 1 is valid for the electric characteristics of the interface circuits and assignment of signal states.
- Transmission process: Asynchronous character transmission in start-stop operation.
- Transmission type: Full duplex.
- Character format: Character composition according to data format in DIN 66022 for start-stop operation. 1 start bit, 7 character bits, 1 parity bit (even), 1 stop bit.
- Transmission speed: 9600 Bits/s.
- Data flow control: none
- Access method: Data transmission from the device to the computer only occurs after a request from the computer.

### Command syntax and format

The following points should be noted for the instruction set:

- Commands are generally sent from the computer (master) to the device (slave).
- The device only responds to requests from the computer. Even error messages are not sent spontaneously from the device to the computer (automation system).
- The commands are transmitted in capital letters.
- Commands and parameters, as well as consecutive parameters, must be separated by at least one space (code: hex 0x20).
- Each individual command (including parameters and data) and all responses are completed with CRLF (code: hex 0x20 hex 0x0d hex 0x0a) and can have a maximum length of 50 characters.
- The dot is used for decimal separators in a floating-point value (code: hex 0x2E).

The details given above generally comply with the recommendations of NAMUR (NAMUR recommendations for the design of electrical plug-in connectors for analogue and digital signals in laboratory MSR devices. Rev. 1.1).

The NAMUR commands and the additional IKA-specific commands are only used as low-level commands for communication between the device and the PC. With an appropriate terminal or communication program, these commands can be transmitted directly to the device.

### NAMUR commands Function

NAMUR commands	Function
IN_NAME	Read device name
IN_PV_4	Read current speed value
IN_SP_4	Read rated speed value
OUT_SP_4	Adjust the rated speed value
START_4	Start evacuation
STOP_4	Stop evacuation
RESET	Switch to normal operating mode
OUT_NAME	Set device name
OUT_SP_42@n	Setting WD safety limit speed with set value echo
OUT_WD1@m	Watchdog mode 1
OUT_WD2@m	Watchdog mode 2

### "Watchdog" functions; monitoring of the serial data flow

If, once this function has been activated (see NAMUR commands), there is no retransmission of the command from the computer within the set time ("watchdog time"), the Evacuation function is switched off in accordance with the set "watchdog" mode or is returned to previously set target values. The data transmission may be interrupted by, for example, a crash in the operating system, a power failure in the PC or an issue with the connection cable between the computer and the device.

#### "Watchdog" mode 1

If event WD1 should occur, the evacuation function is switched off and ER 2 is displayed. Set watchdog time to m (20 - 1,500) seconds, with watchdog time echo. This command launches the watchdog function and must be transmitted within the set watchdog time.

#### "Watchdog" mode 2

If there is an interruption in data communications (longer than the set watchdog time), the speed target value is changed to the set WD safety speed limit. The warning PC 2 is displayed. The WD2 event can be reset with OUT\_WD2@0 - this also stops the watchdog function.

Set watchdog time to m (20 - 1,500) seconds, with watchdog time echo. This command launches the watchdog function and must be transmitted within the set watchdog time.

### Accessories: PC 1.1 cable (device to PC)

Required for connecting the 9-pin socket to a PC.


Fig. 10

### USB A - USB B cable

Required for connecting the USB interface (L) to a PC.


Fig. 11

### Connecting cable

Required for connecting VC 10 and RV 10 auto with the diaphragm vacuum pump IKA Vacstar digital (speed operation).


Fig. 12

## Maintenance and cleaning

The device is maintenance-free. It is only subject to the natural wear and tear of components and their statistical failure rate.

### Cleaning

For cleaning disconnect the mains plug!

Use only cleaning agents which have been approved by IKA to clean IKA devices.

These are water (with tenside) and isopropanol.

Wear protective gloves during cleaning the devices.

Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.

Do not allow moisture to get into the device when cleaning.

Before using another than the recommended method for cleaning or decontamination, the user must ascertain with IKA that this method does not destroy the device.

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- device type.
- serial number, see type plate.
- position number and description of spare part, see [www.ika.com](http://www.ika.com).

### Repair

**Please only send devices in for repair that have been cleaned and are free of materials which might present health hazards.**

For repair, please request the "Safety Declaration (Decontamination Certificate)" from IKA or use the download printout of it from IKA website at [www.ika.com](http://www.ika.com).

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If an error occurs, this is shown in an error on the display (B), e.g. Error 4.

Proceed as follows:

- Switch the device off on the device switch,
- Take corrective measure,
- Switch device back on.

Error	Cause	Reason	Solution
Er 3	Overheat	The pump does not run	Switch off the pump, disconnect it from the mains, let it cool down.
Er 4	Motor blocked/overload	The pump does not run	Switch off the pump, disconnect it from the mains Let the pump cool down. Ventilate the vacuum system to normal pressure.
Er 8	Hall sensor defective	The pump does not run	Service
Er 9	BLP self-test/EEPROM error	The pump does not run	Service
Er 10	Outlet Choked	The pump does not run	User needs to check the outlet
Er 11	Membrane stuck	The pump does not run	Air the pump and run at low speed for self-cleaning
Er 12	Connection problem	The pump does not run	Service: check the motor cable connection

In the event that the fault cannot be eliminated using the measures described or another error code is displayed:

- Please contact the Service Department,
- Send the device and a brief description of the fault.

## Accessories

- Vacuum controller **IKA** VC 10
- Vacuum control valve **IKA** VCV 1 and VCV 2
- Check valve **IKA** VC 10.300
- Vacuum hose **IKA** VH.SI.8
- PC 1.1 cable (RS 232)
- Maintenance kit
- Vacuum safety emission condenser VSE 1
- Woulff bottle/water valve set VSW1

For further accessories see [www.ika.de](http://www.ika.de)

## Product contact parts

Name	Material
Connection spigot	PPS
Connection piece	PPS
Distributor	PPS
Head piece	PPS
Diaphragm valve	FFPM
Diaphragm	NBR/PTFE
Clamp	PPS

## Technical data

	Base unit	Value
Delivery capacity max.	<b>m<sup>3</sup>/h</b>	1.32
Delivery capacity max.	<b>L/min</b>	22
Final pressure	<b>mbar</b>	2
Suction steps		4
Cylinder		4
Connection diameter suction side	<b>mm</b>	8
Connection diameter pressure side	<b>mm</b>	8
Gas ballast valve		no
Inlet pressure min.	<b>mbar</b>	2
Inlet pressure max.	<b>mbar</b>	1030
Boiling point recognition		no
Solvent library		no
Two-position control		yes
Analogue speed-vacuum control		yes
Speed setting option		rotating/pressing knob
Speed min.	<b>rpm</b>	285
Speed max.	<b>rpm</b>	1200
Display		LED
Vacuum sensor		no
Vent valve		no
Condensate separator		no
Emission condenser		no
Volume at minimum pressure	<b>dB(A)</b>	54
Product contact material		PTFE, FFPM, PPS, NBR
Housing material		Cast alu, coated
Appropriate for clean room		no
Explosion-proof		no
Width	<b>mm</b>	150
Depth	<b>mm</b>	370
Height	<b>mm</b>	375
Weight	<b>kg</b>	11.5
Permissible ambient conditions		5 °C to 31 °C at 80 % relative humidity. 32 °C to 40 °C decreasing linearly to a maximum 50 % relative humidity.
Protection class according to DIN EN 60529		IP 20
RS 232 interface		yes
USB interface		yes
Analogue output		no
Voltage	<b>V</b>	100-240
Frequency	<b>Hz</b>	50/60
Power consumption	<b>W</b>	130
Standby power consumption	<b>W</b>	1.5

	Base unit	Value
Protection class		I
Overvoltage category		II
Degree of contamination		2
Use above max. sea level	<b>m</b>	2000
Firmware update		yes
<i>Note</i>		
Protective earth conductor	The protective earth conductor (PE) is not connected to the metal housing. (Internally encapsulated DC power supply)	

*Average values.  
Subject to technical changes!*

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55  
80

**EN**

- A RV 10 auto drive
- B Heating bath
- C Bracket
- D Woulff bottle
- E Glassware

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
#### EN

- A** "Power" switch
- B** "Power" key
- C** Lift position "▲" key
- D** Lift position "▼" key
- E** "Rotating speed" key
- F** "Vent" key
- G** "Vacuum" key
- H** "Menu" key
- J** "Back" key
- K** Start/Stop knob

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## Contents

	Page		
Device setup	3	Interfaces and outputs	25
EU Declaration of conformity	5	Maintenance and cleaning	26
Explication of warning symbols	5	Accessories	26
Safety instructions	6	Error messages	26
Intended use	7	Warranty	28
Unpacking	8	Technical data	29
Useful information	8	Drawings for mounting the glassware	105
Vacuum control	9	Drawings for connecting the hose system	107
Setting up	10	Solvent table (excerpt)	110
Commissioning	16		

## EU Declaration of conformity

We declare under our sole responsibility that this product corresponds to the directives 2014/35/EU, 2006/42/EC, 2014/30/EU and 2011/65/EU and conforms with the following standards or normative documents: EN 61010-1, EN 61010-2-051, EN 61326-1, EN 60529 and EN ISO 12100.

A copy of the complete EU Declaration of Conformity can be requested at [sales@ika.com](mailto:sales@ika.com).

## Explication of warning symbols

Indicates an imminently hazardous situation, which, if not avoided, will result in death, serious injury.

Indicates a potentially hazardous situation, which, if not avoided, can result in death, serious injury.

Indicates a potentially hazardous situation, which, if not avoided, can result in injury.

Indicates practices which, if not avoided, can result in equipment damage.

Danger - Reference to exposure to a hot surface!

## Safety instructions

### For your protection

- Read the operating instructions completely before starting up and follow the safety instructions.
- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations. When working under a vacuum in particular!
- Wear your personal protective equipment in accordance with the hazard category of the medium to be processed. There may be a risk of:
  - splashing liquids,
  - body parts, hair, clothing and jewellery getting caught,
  - injury as a result of glass breakage.

### DANGER

Inhalation of or contact with media such as poisonous liquids, gases, spray mist, vapors, dusts or biological and microbiological materials can be hazardous to user.

- Set up the device in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.
- Ensure that there is sufficient space above the device as the glass assembly may exceed the height of the device.
- Prior to each use, always check the device, accessories and especially the glass parts for damage. Do not use damaged components.
- Ensure that the glass assembly is tension-free! Danger of cracking as a result of:
  - stress due to incorrect assembly,
  - external mechanical hazards,
  - local temperature peaks.
- Ensure that the stand does not start to move due to vibrations respectively unbalance.
- Beware of hazards due to:
  - flammable materials,
  - combustible media with a low boiling temperature.

### CAUTION

Only process and heat media that has a flash point higher than the adjusted safe temperature limit of the heating bath that has been set. The safe temperature limit of the heating bath must always be set to at least 25 °C lower than the fire point of the media used. (acc. to EN 61010-2-010)

- Do not operate the appliance in explosive atmospheres, with hazardous substances or under water.
- Only process media that will not react dangerously to the extra energy produced through processing. This also applies to any extra energy produced in other ways, e.g. through light irradiation.
- Tasks with the device must only be performed when operation is monitored.
- Operation with excess pressure is not permitted (for cooling water pressure see "Technical Data").
- Do not cover the ventilation slots of the device in order to ensure adequate cooling of the drive.
- There may be electrostatic discharges between the medium and the drive which could pose a direct danger.
- The appliance is not suitable for manual operation (except lift movement).
- Safe operation is only guaranteed with the accessories described in the "Accessories" chapter.

- Refer to the operating instructions for the accessories, e.g. vacuum pump, heating bath.
- Position the positive pressure outlet of the vacuum pump under a fume hood.
- Only use the device under an all side-closed exhaust, or a comparable protective device.
- Adapt the quantity and the type of distill and to the size of the distillation equipment. The condenser must work properly. Monitor the cooling water flow rate at the condenser outlet.
- The glass equipment must always be ventilated when working under normal pressure (e.g. open outlet at condenser) in order to prevent a pressure build-up.
- Please note that dangerous concentrations of gases, vapors or particulate matter can escape through the outlet at the condenser. Take appropriate action to avoid this risk, for example, downstream cold traps, gas wash bottles or an effective extraction system.
- Evacuated glass vessels must not be heated only on one side; the evaporating flask must rotate during the heating phase.
- The glassware is designed for operation under a vacuum of up to 1 mbar. The equipment must be evacuated prior to heating (see chapter "Commissioning"). The equipment must only be aired again after cooling. When carrying out vacuum distillation, uncondensed vapors must be condensed out or safely dissipated. If there is a risk that the distillation residue could disintegrate in the presence of oxygen, only inert gas must be admitted for stress relief.

### CAUTION

Avoid peroxide formation. Organic peroxides can accumulate in distillation and exhaust residues and explode while decomposing! Keep liquids that tend to form organic peroxides away from light, in particular from UV rays and check them prior to distillation and exhaust for the presence of peroxides. Any existing peroxides must be eliminated. Many organic compounds are prone to the formation of peroxides, e.g. dekalin, diethyl ether, dioxane, tetrahydrofuran, as well as unsaturated hydrocarbons, such as tetralin, diene, cumene and aldehydes, ketones and solutions of these substances.

### DANGER

The heating bath, tempering medium, evaporation flask and glass assembly can become hot during operation and remain so for a long time afterwards! Let the components cool off before continuing work with the device.

### WARNING

Avoid delayed boiling! Never heat the evaporating flask in the heating bath without switching on the rotary drive! Sudden foaming or exhaust gases indicate that the flask content is beginning to decompose. Switch off heating immediately. Use the lifting mechanism to lift the evaporation flask out of the heating bath. Evacuate the danger zone and warn those in the surrounding area!

### Safety lift

When the device is switched off or the power supply disconnected, the internal safety lift removes the evaporating flask from the heating bath.

The safety lift at loss of power is designed for a maximum total weight (glassware and solvent) of 3.1 kg.

Sample calculation of the maximum load with vertical glassware and a 1 litre flask:

Condenser + receiving flask + evaporating flask + fittings =  
1200 gr + 400 gr + 280 gr + 100 gr = 1980 gr

Maximum loading of solvent = 3100 gr – 1980 gr = 1120 gr

Due to the design, operation of the safety lift cannot be guaranteed for higher loads!

When using other types of condensers such as dry ice or intensive condensers, and also when using return distillation distributors with slip-on condensers, it may be necessary to reduce the load by the amount of the added weight of the glass apparatus. Thus, prior to distillation, check whether the lift goes up without power when laden with the glassware and distillation material.

The safety lift must be checked daily prior to operation. Use the motor to move the lift to the lowest position and press the "Power" key on the front plate or the main switch on the back right side of the device using the maximum total weight of 3.1 kg. → The evaporating flask is lifted out of the heating bath.

When using for the first time or after a long period of inactivity, proceed as follows: use the motor to move the lift several times to the top and bottom positions by pressing the lift function keys prior to interrupting the power supply.

If the safety lift does not work, please contact the **IKA** Service department.

For the evaporating equipment (evaporating flask plus contents), the maximum permissible weight is 3.0 kg! Loads greater than this risk breakage of glass at the vapor tube!

Ensure that the safety lift has been powered down.

When working with large loads, always use low speeds. Unbalanced loads can result in breakage of the vapor tube!

- A vacuum may be formed inside the glassware in the case of power outage. The glassware must be vented manually.

### CAUTION

Never operate the device when the evaporation flask is rotating and the lift is raised. Always lower the evaporation flask into the heating bath first before starting the rotation drive. Otherwise hot tempering medium may be sprayed out!

- Set the speed of the drive so no tempering medium is sprayed out as a result of the evaporation flask rotating in the heating bath. If necessary reduce the speed.
- Do not touch rotating parts during operation.
- Imbalance may result in uncontrolled resonance behavior of the device or assembly. Glass apparatus may be damaged or destroyed. In the event of unbalance or unusual noises, switch off the appliance immediately or reduce the speed.
- The appliance does not start up again automatically following a cut in the power supply.
- The device can only be disconnected from the mains supply by pulling out the mains plug or the connector plug.
- The socket for the mains cord must be easily accessible.
- Always check evaporation processes before starting automatic operation. Automatic operation must not be used for unfamiliar evaporation processes. The following automatic modes are available in addition to manual mode:
  - 100% distillation,
  - Volume-/quantity-dependent distillation.

### For protection of the equipment

- The voltage stated on the type plate must correspond to the mains voltage.
- Socket must be earthed (protective ground contact).
- Removable parts must be refitted to the appliance to prevent the infiltration of foreign objects, liquids, etc.
- Protect the appliance and accessories from bumps and impacts.
- The appliance may only be opened by experts.

## Intended use

### • Use

Together with the accessories recommended by **IKA**, the device is suitable for:

- quick and gentle distillation of liquids,
- evaporation of solutions and suspensions,
- crystallization, synthesis or cleaning of fine chemicals,
- drying of powder and granulate material,
- recycling of solvents.

Mode of operation: Tabletop device

### • Area of use

Indoor environments similar to that a laboratory of research, teaching, trade or industry area.

The safety of the user cannot be guaranteed:

- If the device is operated with accessories that are not supplied or recommended by the manufacturer;
- If the device is operated improperly or contrary to the manufacturer's specifications;
- If the device or the printed circuit board are modified by the third parties.

### • Further remarks

Legal bases and country-specific characteristics, especially for medicine and food:

The operator is responsible for compliance with the laws and regulations applicable to him.


## Unpacking

### • Unpacking

- Unpack the device carefully.
- Any damage should be notified immediately to the shipping agent (post, rail or logistics company).

### • Scope of delivery

	Drive RV 10 auto	Heating bath HB digital	Vertical glassware RV 10.1	Vertical glassware RV 10.10 coated	Condensate flask 1	Bracket 2	Cooler locking device, compl. 3	Vacuum hose (2 x 0.55 m) 4	Water discharge hose (1 x 1 m) 5	Vapour tube 6	Ring spanner 7	Handhold 8	RV 10.4002 vacuum valve for laboratory vacuum system 9	IKA Vacstar digital vacuum pump 10	Operating instructions
RV 10 control V	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x
RV 10 control V-C	x	x		x	x	x	x	x	x	x	x	x	x		x
RV 10 control FLEX	x	x			x	x	x	x	x	x	x	x	x		x
RV 10 auto V	x	x	x		x	x	x	x	x	x	x	x			x
RV 10 auto V-C	x	x		x	x	x	x	x	x	x	x	x			x
RV 10 auto FLEX	x	x			x	x	x	x	x	x	x	x			x
RV 10 auto pro V	x	x	x		x	x	x	x	x	x	x	x			x
RV 10 auto pro V-C	x	x		x	x	x	x	x	x	x	x	x			x
RV 10 auto pro FLEX	x	x			x	x	x	x	x	x	x	x			x

## Useful information

Distillation is a thermal separating process for liquid compounds based on substance-specific, pressure-dependent boiling points through evaporation and subsequent condensation.

The boiling point temperature decreases with decreasing external pressure which means that work is usually done under reduced pressure. In this way the heating bath can be maintained at a constant temperature (e.g. 60 °C). Using the vacuum, the boiling point is set with a steam temperature of approx. 40 °C. The cooling water for the condensation cooler should not be warmer than 20 °C (60-40-20 rule).

A chemical resistant membrane pump with a vacuum controller should be used to create the vacuum. The pump is protected from solvent residue by the addition of a condensate flask and/or a vacuum separator. Working with a jet pump to create a vacuum can only be recommended to a limited extent as the solvents may contaminate the environment when using these systems.

Speed, temperature, piston size and system pressure all affect the evaporator air capacity.

The optimum capacity of the flow-through condenser is approx. 60%.

This corresponds to condensation on approx. 2/3 of the cooling coil.

With larger capacities there is the risk that the uncondensed solvent vapor will be extracted.

In order to prevent this, the condenser loading safety monitoring system can be activated. Before activating the safety monitoring, first carry out several distillations, noting the value given under menu option "EVAPORATOR → Advanced → Max power of cooler". Now calculate the heat rejection rate using the following formula:

$$\text{Cooler power limit} = \text{Max power of cooler} * 1.3$$

Enter this value under menu option "EVAPORATOR → Advanced → Cooler power limit". The device will now stop the measurement and issue an error message if the threshold for the condenser is exceeded. If the configuration of the condenser equipment or the distillation process is modified, then it may be necessary to re-evaluate the maximum condenser loading.

The device is equipped with a piston-operated safety mechanism. The glass apparatus may contain a vacuum following interruption to the power supply; vent the system before switching back on. If the power cuts out, the evaporation flask is automatically lifted out of the heating bath by an integrated gas spring.

**CAUTION**

The safety lift must be checked daily prior to operation. For more information see the section "Safety Instructions – Safety Lift"!

In addition to offering a full range of manual and semi-automatic evaporation operating modes, the RV 10 control rotary evaporator can also be used for fully automatic and quantity-controlled evaporation processes. Standard fittings for fully-automatic operation include a vacuum controller, a cooling water differential temperature measurement sensor and a cooling water flow rate meter. The vacuum can be controlled in both two-position and speed-controlled pumping operation with the help of the integrated vacuum controller. The device is designed for operation with a cooling water supply system (e.g. laboratory thermostat), but can also be run off a water supply line. Please refer to the Technical Data for information on cooling water pressure, temperature constancy, flow rate and local regulations.

We recommend using the RV 10.5001 water regulator valve when operating directly off a water supply line.

This valve allows the cooling water flow rate to be adjusted and automatically cuts off the cooling water supply on completion of the distillation process.

The RV 10.5002 water filter can be used to prevent particles of dirt passing from the water supply line into the valve mechanism. If pressure reduction is required, we recommend fitting the RV 10.5003 pressure regulator valve in the pipe immediately downstream of the pick-off point.

**Automatic volume-dependent distillation:** The apparatus must be brought up to operating temperature. This is achieved by conducting a test distillation.

**RV 10 auto:** The rotary evaporator is fitted as standard with speed control pump, internal venting valve and a pressure measuring chamber.

**RV 10 control (for vacuum):** The external two-position vacuum control is also available for the application which in-house vacuum system is available. Vacuum valve RV 10.4002 must be used.

## Vacuum control

### RV 10 control

A vacuum is made in the glassware with the help of a vacuum pump. The vacuum pump operates at constant speed which generally cannot be adjusted.

Once the target value is reached the suction line is interrupted. The device will fall slightly short of the target value setting due to the minimum time difference from detection of the pressure value, the target value comparison for switching the vacuum valve and the suction line of the vacuum generating system. If the pressure in the system increases again due to the natural leakage rate, the valve opens the suction line.


The switching frequency and accuracy can be set via the hysteresis setting.

With the VC 2.4 pump control connected, the power supply to the pump is also switched on and off accordingly. This reduces the noise level and increases the service life of the pump membranes and motor.

The actual value fluctuates in size of the set hysteresis around the target value.

Automatic boiling point recognition is not accurately possible when using two-position control.

### Schematic view of a two-position vacuum control

### RV 10 auto

Imprecision of the two-position control is avoided by using the speed-vacuum control.

However a vacuum pump with a controlled speed is needed for speed-vacuum control; the RV 10.4002 vacuum valve is not required and must not be connected! To do so, connect the vacuum pump (e.g. IKA Vacstar digital pump) to the port on the rear of the rotary evaporator.

With this kind of control, the speed of the pump, and therefore its suction power, is reduced the closer the measured pressure reaches the target value.

Once the target value is reached, the pump only operates according to the leakage rate.


It is possible to achieve quieter operation and more exact vacuum control.

The speed-vacuum control on the RV 10 is controlled automatically as soon as an appropriate vacuum pump is attached.

Automatic boiling point recognition is possible with this type of control, i.e. the system reaches and keeps the solvent at its boiling point in automatic mode.

As it is not necessary to know the boiling point of the solvent when in automatic boiling point recognition mode (unlike in the case of volume-dependent distillation), these two operating modes are mutually exclusive.

### Schematic view of speed-vacuum control


## Setting up

### Drive RV 10 auto

Loosen transportation lock!


**CAUTION**

- Hold the lift with your hand on the height position and remove the thumb screw (1) by turning counterclockwise on the back of the appliance.
- Once the transportation lock has been removed, the lift moves slowly to its upper end position. The distance is approx. 140 mm.
- Connect the supplied power cable to the connection socket (2).

Adjustable base

**Note:** Only one base foot is adjustable.

Fix the bracket


Insert the bottle and attach the supplied hose connectors to the bottle

Set the drive at an angle of approx. 30°

- Remove the clamping device for the angle setting of the rotation drive on the right side of the lift by rotating the knurled screw counterclockwise (gently press and rotate the knurled screw (1) at the same time to move it further in or out).
- Set the drive at an angle of approx. 30°.

- Then secure the rotation drive from being accidentally turned by tightening the knurled screw in a clockwise direction.

### Heating bath

**CAUTION**

Refer to the chapter "Operation" in the heating bath instruction manual!

- Place the heating bath on the stand of the rotation drive and push it into the left position.

**Note:** Data is exchanged between the drive unit and the heating bath by means of an infrared link (1). Please note that reliable communication is only guaranteed when the infrared beam has clear line-of-sight to the detector.




### Glassware

**NOTE:** Please read the operating instructions of the glassware for the safe handling of laboratory glassware!

- There is a locking knob for locking or unlocking the locking device.

If a red mark is shown, the locking device is unlocked, otherwise it is locked.

To lock or unlock the locking device, push the locking knob to the end position.

- Open the locking device on the drive head by turning it 60° counterclockwise. Then you can see a red mark.
- Feed the vapour tube in until it stops.
- Then, lock the locking device by turning it clockwise by 60°.
- Pushing the lock knob to the end position, the red mark should be covered and invisible.
- The vapour tube is not allowed to be pulled out!
- Check the correct axial locking device on the vapour tube.
- Keep the red mark invisible.

Fit the condenser seal



- Place the condenser seal in the condenser receptacle and fit the glassware to the device according to the assembly instructions.

#### First use


- Insert vapor tube (1).
- Insert the vacuum seal (2). Pay attention to the correct position of the vacuum seal!
- Slide the union nut (3b) over the flange on the condenser (3a).
- Also slide the annular spring (3c) over the flange on the condenser (3a).
- Position the condenser (3a) on top of the seal (2).
- Tighten the union nut (3b) by hand onto the threaded flange. Run the device for 20 minutes at 120 rpm. Finally, retighten the union nut (3b) by hand.

**Note:** Follow the mounting instructions for the glassware.

#### Assembling the vertical glassware cooler locking device

- Assemble the condenser locking device according to the diagram.
- Mount the plate (1) using the thumb screw (2).
- Put the support rod (4) on the plate (1) and attach it with the screw nut (3).
- Attach the rubber protector (5).
- Fasten the Velcro® band (6) to the support rod (4).
- Secure the vertical glassware with the Velcro® band (6).

No stand is required if the condenser has been properly fitted and secured by tightening the condenser union nut on the rotary drive.

The purpose of the stand is solely to prevent the condenser from twisting.

**CAUTION** If the stand is fitted incorrectly, the strong forces at the Velcro strap can cause stresses within the glass, which in turn can cause damage to the glass condenser.

After the condenser has been fitted, the stand is attached to the rotary drive. Make sure that the condenser is installed parallelly to the lift body.

#### Mounting the glassware

**Note:** Observe "Drawings for mounting the glassware" on page 105.

Item	Designation	Quantity					
		RV 10.1 non-coated RV 10.10 coated	RV 10.2 non-coated RV 10.20 coated	RV 10.3 non-coated RV 10.30 coated	RV 10.4 non-coated RV 10.40 coated	RV 10.5 non-coated RV 10.50 coated	RV 10.6 non-coated RV 10.60 coated
1	Receiving flask	1	1	1	1	1	1
2	Clamp (stainless steel)	1	1	1	1	1	1
4	Stopcock	1	1	1	1	1	1
5	Pipe tube	1	-	1	1	1	1
6	Condenser	1 Vertical condenser	1 Diagonal condenser	1 Vertical-intensive condenser	1 Dry ice condenser	1 Vertical condenser	1 Vertical-intensive condenser
7	Connection	1 Vacuum connection	1 Introduction sleeve	1 Vacuum connection	-	1 Vacuum connection	1 Vacuum connection
8	Clamp (plastic)	1	1	1	-	1	1
10	Evaporation flask 1000 ml	1	1	1	1	1	1
11	Ball joint clamp RV 05.10	1	1	1	1	1	1
12	Screw joint cap	4	4	4	2	4	4
13	Hose connection	4	4	4	2	4	4
14	Manifold	-	-	1	-	1	1
15	Cover cap	-	-	-	1	-	-
16	Cartridge	-	-	-	1	-	-
17	Ring, slotted	-	-	-	1	-	-
18	Washer	1	-	1	1	1	1
19	O-ring	-	-	-	1	-	-
20	Plug, high vacuum	-	-	-	-	1	1


#### Note:

The PTFE tube (5) with the washer (18) can be mounted optionally on the vertical glass sets.

It serves for filling the evaporator flask when there is a vacuum in the glass set.

Solvent can be sucked through the PTFE tube into the evaporator flask by opening the stopcock (4).

#### Mounting the washer


**Note:** Pay attention to the correct position of the washer.

#### Removing the condenser


- Use the ring spanner provided to loosen union nuts that are tightly fitted.
- Loosen the union nut by turning counterclockwise.
- Remove the Velcro®.

**Note:** Ring spanner is only used for removing the condenser, while fixing the condenser please use hands.

## Hose system

**Note:** Observe "Drawings for connecting the hose system" on page 107.


## RV 10 control operation

- Fit RV 10.4002 vacuum, in the bracket provided for this purpose and connect the vacuum hose to the valve.

## RV 10 auto operation

- No additional valve is required for vacuum normal operation with a speed-vacuum control pump. Connect the pump's vacuum inlet directly to the RV 10 pump hose connection.
- Connect the pump controller cable to RV 10 control.


**Note:** When connecting the pump controller cable, make sure that the snap-and-lock plug has been inserted into the stop position inside the pump controller cable connection.


## Connect RV 10 temperature sensor (dT)

- Connect the temperature sensor to the lower socket (dT). Always check for correct direction of the arrow.

**Note:** When connecting the temperature sensor, make sure that the snap-and-lock plug has been inserted into the stop position inside the temperature sensor connection.


## Connect RV 10.4002


## Water

- Connect the water inlet to the water supply (with optional RV 10.5001 water regulator valve). Please refer to the technical data relating to the water supply. The RV 10.5001 water regulator valve is not suitable for use on the condenser unit because it would reduce the flow rate too greatly.
- Connect the water hoses to the glass condenser (short hose (1) = lower outlet, long hose (2) = upper inlet) and fix hoses in place.


- Connect the water outlet hose (supplied) by pushing it onto the nipple until the stop inside the connection socket is reached.
- Note:** Ensure that the condenser inlet and condenser outlet are connected correctly. Automatic distillation cannot be carried out if the water inlet and outlet hoses are not fitted according to the instructions and illustrations, since this would cause the temperature readings for the inlet and outlet to be measured incorrectly.


- Release the hose connection using the supplied tool.

- Insert the valve connector cable (RV 10.5001, RV 10.4002) into the appropriate socket. Connect the speed-vacuum control pump. As a result, the RV 10 control automatically switches to normal speed-vacuum control operation.


- Connect "LOAD" to the Condensate flask and the Condensate flask to the vacuum connector on the condenser using the vacuum hoses supplied. Please note that the vacuum must always be connected at the highest connection point on the condenser.


## Commissioning

### Working screen at the time of delivery

The start screen appears for a few seconds after the device is switched on. The device name and the software versions are displayed.

After this, the working screen appears automatically on the display.

An information prompt then appears to download the IKA firmware update tool.


### Explanation of symbols on the working screen

The symbols displayed change depending on the status and settings of RV 10 auto.

Symbol	Designation	Description
	Bluetooth®	This symbol means the RV 10 is communicating via Bluetooth® with a PC or tablet. The symbol no longer appears if no Bluetooth® communication is being performed.
	USB	This symbol means RV 10 is communicating via a USB cable. The symbol no longer appears if no USB cable is being used for communicating with the station.
<b>LIMIT</b>	Limit	This symbol indicates if a minimum or maximum limit of a set value is reached.
	Temperature sensor	This symbol appears when the temperature display is active in the display.
<b>PC</b>	PC control	This symbol means that device is connected to a computer and that it is being controlled by the computer.
<b>PR</b>	Program controlled	This symbol means that device is being controlled by a program.
	Rotation	This symbol indicates that the rotation is in continuous mode and indicates the direction of rotation.
	Drying	This symbol indicates that the drying option was selected (only in 100 % mode).
	Flow rate	This symbol indicates that a fluid (water) is running through the system.
	Cooler connected	This symbol indicates that a cooler is connected and running.
	Venting	This symbol indicates that the venting valve is open. If this symbol is not shown the venting valve is closed.
<b>vac</b>	Vacuum	This symbol indicates that the pressure inside the RV 10 system is lower than the atmosphere pressure.
	Detecting	This symbol is shown during the automatic evacuation mode.
	Boiling	This symbol indicates that a boiling point was found (only in automatic mode).

## Menu structure

				Factory settings		
	Advanced	Flow control.....	<input checked="" type="checkbox"/>			
		Bath medium.....		Water		
		Cooler power limit.....		900 w		
		Max power of cooler.....		0 w		
	Rotation	Interval mode	Interval	Ccw Value.....	00:10	
				Cw Value.....	00:10	
			Activate.....		-	
		Speed limit	Minimum.....		5 rpm	
			Maximum.....		300 rpm	
			Display.....		-	
	Timer function	Time hh:mm:ss.....		00:00:00		
		Beep after timeout.....		-		
		Activate.....		-		
	Display.....		-			
	Display	Timer.....	<input checked="" type="checkbox"/>			
		Flow rate.....	<input checked="" type="checkbox"/>			
		Δ T.....	<input checked="" type="checkbox"/>			
		Cooler power.....	<input checked="" type="checkbox"/>			
		Display.....		-		
	Programs	Program 1-10	Select.....			
			Edit	Edit	Pressure.....	
					Rotation (Speed).....	
				hh:mm:ss.....		
			Insert.....			
		Delete.....				
		Last Measurement	Save as.....			
		Edit	Edit	Pressure.....		
				Rotation (Speed).....		
			hh:mm:ss.....			
	Insert.....					
	Delete.....					
	Password	Value.....		0		
		Value.....		0		
		Value.....		0		
	Continue after power failure	Activate.....		-		
		Time mm:ss.....			00:10	


				Factory settings	
	Service	Valves	Venting valve.....	-	
			Vacuum valve.....	-	
			Water valve.....	-	
		Pump.....		-	
	Settings	Language	English.....	<input checked="" type="checkbox"/>	
			Deutsch.....		
			Français.....		
			Español.....		
			Italiano.....		
			Português.....		
			Русский язык.....		
			中文.....		
			한국어.....		
		日本語.....			
		Units	Pressure	mbar.....	<input checked="" type="checkbox"/>
				hPa.....	
				mmHg.....	
				Torr.....	
		Amount	mL.....	<input checked="" type="checkbox"/>	
			g.....		
		Display	Background	Black.....	<input checked="" type="checkbox"/>
				White.....	
Brightness Value.....			100%		
Firmware update info.....		-			
Sound	Volume Value.....		100%		
	Key Tone.....		-		
Factory settings	Menu values.....		-		
	Programs.....		-		
	Solvent library.....		-		
	All.....		-		
Communication	Device name.....		RV 10 auto		
	Bluetooth.....		<input checked="" type="checkbox"/>		
	Labworldsoft 5 protocol.....		<input checked="" type="checkbox"/>		
Information	Firmware update info.....		-		
	Display version.....		-		
	Logic version.....		-		
	Pressure max.....		1100 mbar		
	Pressure min.....		1 mbar		
	Pump % max.....		100%		
	Pump % min.....		0%		
	Rot. speed max.....		300 rpm		
	Rot. speed min.....		5 rpm		
	Operating hours.....		0 hours		

Change the menu items, values or settings by turning the rotary knob and confirm by pushing it.

## Menu (details)

### 1. Evaporator

#### Modes

**Automatic:** In this mode, the system will detect the actual evaporation by checking the temperature difference between water input and water output. The default setting of the medium temperature is 60 °C, the user can set the medium temperature manually. If the expected evaporation temperature is higher than 90 °C, please change the medium from water to oil, and use high temperature medium, e.g. silicone oil. In this mode, measurement will start automatically when the heating bath reaches the set temperature, the temperature difference will be below a threshold value, the water flow rate must be within 30-100 L/h, and the lift moves down automatically. You can change the speed manually.

**Manual:** In this mode the user can set all values manually. When the measurement is started the system is evacuated until it reaches the set pressure value. The lift has to be moved manually.


**Pump %:** In this mode the pump can be run in a long term, setting a value between 100 % and 1% of the possible pump speed.

**Volume:** This mode is used to distill a specific amount of the used solvent. The measured values of the cooling water flow rate and the cooling water temperature difference from the basis for a heat balance calculation used to determine the quantity of distillate at every stage of the distillation. The distillation process is stopped when the specified quantity of distillate has been reached.

**Solvent:** Select the solvent you want to distill.  
**Target:** Enter the amount, which you want to distill off.  
**Efficiency:** Balancing requires a precise thermal efficiency to be entered.  
**Adjustment:** The efficiency is dependent on a combination of factors that affect the distillation process, so it is initially estimated. The first time a distillation process is carried out should therefore be used for adjustment. The actual efficiency will be recalculated after the measurement based on the target volume and the actually distilled volume.

**Program:** This entry is a link to the "Programs" menu.  
**100 %:**  
**Solvent:** The distillation process is stopped when the measured cooling water temperature difference falls below a threshold value, i.e. as soon as a solvent has been completely distilled off.  
**Drying:** The cooling water temperature difference is not monitored if this option is activated, e.g. use in processes for drying powdered media.

#### Cleaning

**Start after...:** Set the minimum time that a measurement has to run before the cleaning after the measurement is activated. The cleaning process starts when the measurement has been running for at least the set time.


**Duration:** Set the duration of the cleaning process.

**Pump speed:** The set pump speed in the cleaning menu is used to set the speed of the pump during the cleaning process.

**NOTICE** Higher speed causes partial vacuum in the glassware if the system is closed. High speed cleaning affects the pump membrane lifetime.

**Clean now:** This starts the cleaning process manually.


#### Limits

**Manual:** This option allows user to set the permitting operating vacuum range.

**Pump %:** This option allows user to set the permitting operating pump speed range.

#### Hysteresis


The (vacuum) hysteresis value describes the difference between the closing and opening pressures for the vacuum valve. Setting the hysteresis value too low can have a negative effect on the vacuum pump and the vacuum valve. The hysteresis function is only used for indoor vacuum or uncontrollable vacuum pumps.

**Automatic:** Hysteresis is always 10 % of the actual pressure.

**Manual:** Manual value specification of hysteresis.


#### Solvent Library

User can select the solvent which will be distilled, the parameters (such as, Name, Formula, Bath temperature, Rotation speed, Boiling point, Efficiency, Heat capacity, Enthalpy, Density and i-factor) will be shown in the library. User can set the evaporating bath temperature, efficiency and the rotation speed. It is possible to add up to five new customized solvents to the solvent list. The customized solvents appear in orange and can be found on the top of the solvent library.


**Name and formula:** The name and the formula help to identify the solvent.

### Calculator


It helps you to calculate the boiling point (vapor temperature) of the solvent, heating bath temperature and pressure. If one of the 3 values has been defined, the other 2 values will be given.

#### Actions after distillation


It is possible to define actions when the measurement is stopped (manually or automatically).

### 2. Rotation

### Advanced

**Flow control:** This setting can be deactivated in the manual mode if dry ice is used instead of a cooler.

**Bath medium:** Select either water or oil as the medium for the heating bath. If water is the medium, please select water, the temperature will be from 20 °C to 90 °C. If oil is the medium, please select oil, the temperature range will be from 20 °C to 180 °C.

**Cooler power limit:** This setting is used to limit the cooler power. If the actual cooler power is higher than the cooler power limit a warning appears and the measurement is stopped.

**Max power of cooler:** This value shows the maximum cooler power value of a measurement. This value serves as information and cannot be changed.

#### Interval mode

The rotation can be run in an interval mode of alternating the rotation direction.


**Interval:** The interval sets the time how long the rotation takes place in one direction.

**Activate:** This entry activates the interval mode for the rotation.

#### Speed limit

The "Speed limit" menu allows the user to set the desired maximum and minimum speed limit for the rotation.

### 3. Timer

**Note:** The user can stop the measurement procedure before expiry of the set time. In this case the countdown of the timer is interrupted.

**Beep after timeout:** This menu item allows the user to turn on/off the sound that occurs when the timer reaches 00:00:00. A tick shows that the option is activated.

**Activate:** In this menu the user can activate or deactivate the timer function. A tick shows that the option is activated.

#### Display

In the "Timer" menu the user can specify that the timer is displayed on the display/working screen. A tick shows that the option is activated.

#### Timer function


**Time:** This setting allows the user to specify the actual time for the measurement procedure in the manual mode. A default time can also be set for the timer, which allows the user to start the measurement task for a standard time. The device stops automatically after expiry of the set time, and the set time used for the measurement procedure appears on the display.

### 4. Display

In this menu the user can set which values (Timer, Flow rate,  $\Delta T$  and/or Cooler power) should be shown on the main screen. A tick shows that the option is activated.

## 5. Program

### Program 1 - 10

Under menu "Programs", 10 user-defined pressure-rotation (rpm)-time profiles can be created. In addition, the user can define whether the interval mode is activated or not in the programs.

If the interval mode is activated, the run time/stop time value from "Interval Mode" setting will be taken.

**Note:** If the user needs to activate the interval mode in one segment of a program, he should set the "Ccw/cw time" in menu option "Interval", meanwhile activate the "Ccw/cw" function in menu option "Interval Mode" (see section "Rotation").

**Select:** Select the program.

**Edit:** Edit the selected program parameters. Start to edit the selected program parameters by pressing menu option "Edit" with start/stop knob. The user can edit, insert or delete one selected program segment in the program. When user edits the program time for at least one segment, a tick (✓) for respective program will appear.

**Delete:** Delete the selected program. If a selected program is deleted by pressing on menu option "Delete" with start/stop knob, all the program parameters will be emptied. The tick (✓) disappears.

### Last measurement

**Save as:** Save the distillation sequence as a program.

**Edit:** Edit the selected program parameters.

Start to edit the selected program parameters by pressing on menu option "Edit" with start/stop knob. The user can edit, insert or delete one selected program segment in the program. When the user edits the program time for at least one segment, a tick (✓) for the respective program will appear.

**Delete:** Delete the selected program. If a selected program is deleted by pressing on menu option "Delete" with start/stop knob, all the program parameters will be emptied. The tick (✓) disappears.

### Details for editing the program

When a program is being edited, the following screen appears.

No.	Pressure	Rotation	Time
01	500 mbar	40 rpm	00:00:20
02	400 mbar	50 rpm	00:00:20
03	300 mbar	60 rpm	00:00:20
04	200 mbar	70 rpm	00:00:20
05	100 mbar	80 rpm	00:00:20


In this program, user can define up to 10 segments. The selected segment is highlighted. Then, the user can edit, insert or delete a segment in this program. When the "Back" button is pressed after editing, the program is saved automatically.

**Edit:** When the background of the selected value turns yellow the user can change the pressure, the speed or the time.


**Insert:** With the insert option, a new segment will be insert below the selected segment.

**Delete:** When deleting a highlighted segment, the settings will be emptied. The section jumps to next segment.


### Sample for editing the program

### Example of saving the last measurement

## 6. Safety

### Continue after power failure

If this option is activated the measurement will continue after a break caused by a power failure. This option is only available in Automatic, 100 % and Volume mode.


**Activate:** If tick is shown the option is activated.

**Time:** The measurement will be continued if the the power comes back within the set time.

### Password

It is possible to set a password to secure the main menu from setting changes. To deactivate the password set the password 0 0 0.

## 7. Service

### Valves

In this menu the connected valves can be opened or closed individually.

### Pump


User can start or stop the pump manually without actual distilling.

To check functionality for servicing purposes, use the Service menu to directly active or deactivate included and optional valves and vacuum pump.

### Temperature adjustment

User can set the different temperature to zero when there is no vapor liquefaction. The correct temperature adjustment helps the auto distillation mode to start and stop the distilling.

## 8. Settings

### Languages

The "Languages" option allows the user to select the desired language.

### Units

The "Units" option allows the user to select the desired unit for displaying the pressure and the amount.

### Display

The "Display" option allows the user to change the background color and brightness of the working screen, and display "Firmware Update Info".

### Sound

The "Sound" option allows the user to set the volume and activate/deactivate the key tone.

### Factory Settings

This option allows the user to reset Menu values, Programs, Solvent library, or all the system settings to the original values at time of delivery.

### Communication

**Device name:** The device name is used to identify a device.

**Bluetooth:** The "Bluetooth" option allows the user to activate/deactivate the "Bluetooth" function.

**Labworldsoft 5 protocol:** To use the RV 10 auto in Labworldsoft 5 this option has to be active.

### Information

This option offers the user an overview of the most important system settings of the device.

## Setting the lower end stop

### CAUTION

Depending on the size of the piston, the setting angle of the rotation drive and the position of the heating bath and lift, the evaporation flask can be in contact with the bottom of the heating bath.

### CAUTION

**Glass may break!**

Limit the lower lift position with the variable end stop.

The lift can be moved to any desired position in manual mode using the "▲" and "▼" keys. There is no automatic cut-off if a collision occurs.

- Press the "▼" key until the lift reaches the desired position.


**Note:** The evaporation flask should be 2/3 immersed in the heating bath.

- To adjust the stop element (1), press the center button (2) on the front of the lift.

- Move the stop (1) to the desired position

- Press the "▲" key until the lift has reached the upper end stop.

**Note:** The path is limited from 0 – 6 cm.

Check that the limit stop is functioning properly:

- Lower the lift by holding down the "▼" key.
- The drive stops automatically when the preset desired lower end position has been reached.
- Move the drive back to the upper position.

To activate the safety lift following a long period of inactivity, use the motor to move the lift to the top or bottom position several times prior to starting the distillation. (See the section entitled Safety Instructions - Safety Lift!)

## Fill the evaporation flask

**Manual filling:** You can fill the evaporation flask manually prior to creating the vacuum. The evaporation flask should not be filled more than half its volume.

**Automatic filling:** Prior to filling the evaporation flask, a vacuum controller is used to regulate the glass apparatus to the target pressure.

- Now fill the evaporation flask using the back feed line.
- Due to the vacuum present, the solvent is suctioned into the evaporation flask. This enables you to keep solvent loss due to suctioning to a minimum.

### CAUTION

The maximum allowed load (evaporation flask and contents) is 3 kg.

## Setting up the heating bath

### CAUTION

**Refer to the operating instructions for the heating bath IKA HB!**

- Move the lift to the bottom position and check the position of the heating bath in relation to the evaporation flask. When using larger evaporation flasks (2 or 3 litres) or depending on the angle of the rotation drive, you can move the heating bath 50 mm to the right.

- Fill the heating bath with the tempering medium until the evaporation flask is surrounded by tempering medium to 2/3 of its volume.

- Switch on the rotation drive and slowly increase the speed

**Note:** Avoid creating waves.

- Switch on the heating bath using the main on/off switch.

**Note:** Avoid stress on the glass due to different evaporation flask and heating bath temperatures when lowering the evaporation flask into the heating bath!

**Note:** If non-original accessories are used that are not supplied by IKA, then it is possible that the 50 mm travel range provided by the heating bath will not be sufficient. This applies in particular when using 3 L evaporator flasks with foam brake.

Use the IKA RV 10.3000 extension plate to extend the heating bath travel range by 150 mm.

## Interfaces and outputs

The device can be operated in "Remote" mode via an RS 232 or USB interface using *labworldsoft*® laboratory software.

The RS 232 interface at the back of the device, is fitted with a 9-pole SUB-D jack and can be connected to a PC. The pins have serial signals. The USB interface is located on the left-hand side of the display on the drive unit and can be connected to a PC using the USB cable supplied.

**Note:** Please comply with the system requirements together with the operating instructions and help section included with the software.

### USB Interface

The Universal Serial Bus (USB) is a serial bus for connecting the device to the PC. Equipped with USB devices can be connected to a PC during operation (hot plugging). Connected devices and their properties are automatically recognized. Use the USB interface in conjunction with *labworldsoft*® for operation in "Remote" mode and also to update the firmware.

Installation:

First, download the latest driver for IKA devices with USB interface from: <http://www.ika.com/ika/lws/download/usb-driver.zip>. Install the driver by running the setup file. Then connect the IKA device through the USB data cable to the PC. The data communication is via a virtual COM port. Configuration, command syntax and commands of the virtual COM ports are as described in RS 232 interface.

### Serial interface RS 232

Configuration:

- The functions of the interface connections between the stirrer machine and the automation system are chosen from the signals specified in EIA standard RS 232 in accordance with DIN 66 020 Part 1.
- For the electrical characteristics of the interface and the allocation of signal status, standard RS 232 applies in accordance with DIN 66 259 Part 1.
- Transmission procedure: asynchronous character transmission in start-stop mode.
- Type of transmission: full duplex.
- Character format: character representation in accordance with data format in DIN 66022 for start-stop mode. 1 start bit; 7 character bits; 1 parity bit (even); 1 stop bit.
- Transmission speed: 9600 bit/s.
- Data flow control: none
- Access procedure: data transfer from the device to the computer takes place only at the computer's request.

Command syntax and format:

The following applies to the command set:

- Commands are generally sent from the computer (Master) to the device (Slave).
- The lab device sends only at the computer's request. Even fault indications cannot be sent spontaneously from the machine to the computer (automation system).
- Commands are transmitted in capital letters.
- Commands and parameters including successive parameters are separated by one space (Code: hex 0x20).
- Each individual command (incl. parameters and data) and each response are terminated with CR LF (Code: hex 0x0d hex 0x0a) and have a maximum length of 128 characters.
- The decimal separator in a number is a dot (Code: hex 0x2E).

The above details correspond as far as possible to the recommendations of the NAMUR working party (NAMUR recommendations for the design of electrical plug connections for analogue and digital signal transmission on individual items of laboratory control equipment, rev. 1.1).

The NAMUR commands and the additional specific IKA commands serve only as low level commands for communication between the rotary evaporator and the PC. With a suitable terminal or communication program these commands can be transmitted directly to the rotary evaporator. The IKA software package, *labworldsoft*®, provides a convenient tool for controlling the rotary evaporator and collecting data under MS Windows, and includes graphical entry features, for motor speed ramps for example.

The following table summarizes the (NAMUR) commands understood by the IKA control equipment.


Abbreviations used:

m =	Numbering parameter (integer)
X = 2	Temperature heating bath
X = 3	Safety temperature heating bath
X = 4	Rotation speed
X = 60	Interval time (1 - 5999 seconds, 1 <= m >=5999)
X = 61	Timer (1 - 1440 minutes, 1 <= m >=1440)
X = 62	Upper lift position (OUT_SP_62 1-> drive lift up)
X = 63	Lower lift position (OUT_SP_62 1-> drive lift down)
X = 66	Vacuum controller pressure
X = 70	Vacuum controller hysteresis
X = 74	Tempering medium heating bath (OUT_SP_74 0=oil, OUT_SP_74 1=water)

NAMUR Commands	Function
IN_NAME	Request designation
IN_PV_X X = 4	Read actual value
IN_SOFTWARE	Request software Id-number, date and version
IN_SP_X X = 4	Read target value input
OUT_SP_X m X = 1,60,61,62	Set target value to m
RESET	Switch to normal operation
START_X X = 4,60,61,62	Switch on appliance (remote) function
STATUS	Status output 0:Manual operation without interruption 1:Automatic operation Start (without interruption) ERROR z (z error number see table)
STOP_X X = 1,60, 61, 62	Switch off appliance function. The variable set with OUT_SP_X remains


### PC 1.1 cable

This cable is required to connect the 9-pin connector to a PC.

### USB 2.0 cable

This cable is required to connect USB interface to a PC.

## Maintenance and cleaning

The device is maintenance-free. It is only subject to the natural wear and tear of components and their statistical failure rate. The seal on the glass condenser should be checked at regular intervals and replaced if necessary.

### Cleaning

- For cleaning disconnect the mains plug!
- Use only cleaning agents which have been approved by IKA to clean IKA devices. These are water (with tenside) and isopropanol.
- Wear protective gloves during cleaning the devices.
- Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.
- Do not allow moisture to get into the device when cleaning.
- Before using another than the recommended method for cleaning or decontamination, the user must ascertain with IKA that this method does not destroy the device.

### Spare parts order

When ordering spare parts, please give:

- device type.
- serial number, see type plate.
- position number and description of spare part, see [www.ika.com](http://www.ika.com),
- software version

### Repair

**Please only send devices in for repair that have been cleaned and are free of materials which might present health hazards.** For repair, please request the "Safety Declaration (Decontamination Certificate)" from IKA or use the download printout of it from IKA website at [www.ika.com](http://www.ika.com).

If your appliance requires repair, return it in its original packaging. Storage packaging is not sufficient when sending the device - also use appropriate transport packaging.

**Note:** For glassware unit, please contact your local dealer and do not send the glassware back to our factory.

### Safety lift

The safety lift must be checked daily prior to operation!

After a long period of inactivity (approx. four weeks), use the motor to move the lift to the top or bottom position several times prior to starting distillation. (See the section entitled Safety Instructions - Safety Lift)! If the safety lift is not working, please contact the **IKA** Service department.

## Accessories

For further accessories see [www.ika.com](http://www.ika.com).

## Error messages

Any malfunctions during operation will be identified by an error message on the display.

Once a serious error message has been displayed, the lift moves to the top end position and the device can no longer be operated.

Proceed as follows in such cases:

- switch off device using the mains switch,
- carry out corrective measures,
- restart device.

Error message	Effect	Cause	Correction
No pressure change	Pressure gradient analysis does not show any deviation	Pump is not running  Vacuum hose isn't connected  RV 10.4002 valve control cable is disconnected Failure in hose system  Recipient isn't closed	Check pump control cable connection Check pump power supply Switch on pump Check vacuum hose connection (pump, recipient) Check RV 10.4002 valve control cable connection Check correct vacuum hose connection (IN / OUT/ VENT) Close recipient
System not tight	Pressure gradient analysis shows pressure deviation, but set value cannot be reached. Pressure set value cannot be reached.	Vacuum hose connection isn't tight Recipient isn't tight Vacuum pump power isn't sufficient Venting valve leakage	Check vacuum hose connection  Check recipient Check technical data of the pump  Contact service department
Pressure out of range	Vacuum is out of range Vacuum is bigger than room pressure	Pressure in the recipient is too high Sensor is broken.	Check the air flow of pump and vent recipient Contact service department
Vacuum sensor error	Vacuum sensor output is too low	Sensor is not connected Sensor is broken.	Check the sensor connection Contact service department
Vacuum calibration error	Vacuum calibration value is out of range	Input wrong vacuum actual value when doing calibration. Set wrong vacuum calibration value Sensor is broken	Check and calibration again  Check the set calibration value  Contact service department
Boiling point detecting error	Boiling point cannot be detected. Differ-temperature analysis does not show boiling point	Heating bath does not heat up  Solvent with extreme low boiling point is used	Check heating bath  Solvent distillation only manually
Venting error	Pressure gradient analysis does not show any deviation after pressing the "Venting" button.	Recipient wasn't evacuated  Failure in hose system  Venting valve doesn't work Front foil button is damaged	Evacuated recipient  Check correct vacuum hose connection (IN / OUT/ VENT) Contact service department Contact service department
Analog pump not connected	Analog pump is disconnected		Check the analog pump cable connection
No rotation	Rotary drive doesn't start or exhibits very large speed variation	Rotary drive is overload or blocked Motor or motor cable is damaged	Reduce the load  Contact service department
Temperature sensor not connected	Differential temperature shows big minus value	Temperature sensor is not connected Temperature sensor is damaged	Check the Temperature sensor connection Contact service department
No temperature difference	No differential temperature increase of cooling water	No solvent or too less solvent is distilled No water flow Water flow is wrongly connected	Check the solvent  Check flow of cooling water Check flow direction of cooling water
Temperature out of range	Differential temperature out of the measurement range	Too fast distilling Temperature sensor is damaged	Increase the set vacuum Contact service department

Error message	Effect	Cause	Correction
Temperature calibration error	Temperature calibration value is out of range	Wrong simulator resistance is chosen PCB is broken	Check the simulator resistance  Contact service department
Temperature adjustment error	Water in and water out sensor temperature have big deviation	The cooling water temperature is changing  Temperature sensor is not calibrated Temperature sensor is broken	Check the chiller and wait until the chiller reaches the set temperature Make sure there is no distilling Contact service department  Contact service department
Inside temperature too high	Analysis of internal PCB temp. Sensor. Inside temperature reaches limit value.	Room temperature > 40°C PCB is broken	Switch off device to cool down Contact service department
Flow rate out of range	Water flow rate of cooling water is out of range	Water flow rate is too high Water flow sensor is damaged	Reduce the water flow rate Contact service department
Flow rate too low	Water flow rate of cooling water is less than the minimum limit value	Water flow rate is too low Water flow sensor is damaged	Increase the water flow rate Contact service department
Condenser overload	Cooling power limit value is exceeded	Condenser limit value is too lower Distilling is too fast	Increase the condenser overload limit  Increase the set vacuum
Distillation dry	There is no more solvent distilled	Distillation is completed	Stop the distilling
Lift error	Lift can not reach the end position	Lift is overloaded or lift is blocked Position sensor failure Lift motor, cable or PCB is damaged	Check the lift  Contact service department Contact service department
Internal communication error	Internal Communication Error	Internal communication error	Switch device on/off, Contact service department
PC communication error	Error command received from PC via RS 232 or USB interface	LWS uses wrong device for RV 10 control Wrong tools are used to communicate with RV 10 control Unstable connection	Check the LWS setting  Check the PC tools  Check the connection
Heating bath communication error	Communication with heating bath is not available	Infrared interface is blocked Heating bath is switched off or detect error Heating bath with infrared interface is not present PCB is damaged	Check and clean the infrared interface Switch on heating bath  Check the heating bath  Contact service department
EEPROM error	Error while reading or writing the EEPROM	EEPROM is defective	Contact service department
Storage error	Read or write internal flash error	Read or write internal flash error	Contact service department

If the actions described fail to resolve the fault or another error message is displayed, then take one of the following steps:  
- contact **IKA** service department,  
- send the device for repair, including a short description of the fault.

## Warranty

In accordance with **IKA** warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine directly to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual.

**Note:** For claims on glassware units, please contact your local dealer and do not send the glassware back to our factory.

## Technical Data

Operating voltage range	<b>Vac</b>	(100...240) ± 10%
Rated voltage	<b>Vac</b>	100...240
Frequency	<b>Hz</b>	50/60
Power input without heating bath	<b>W</b>	100
Power input operation "standby"	<b>W</b>	3.3
Speed	<b>rpm</b>	0/5...300
Speed tolerance	<b>rpm</b>	± 1 (Set speed < 100 rpm)
	<b>%</b>	± 1 (Set speed ≥ 100 rpm)
Speed display		digital
Dimensions of visible display area (W x H)	<b>mm</b>	70 x 52
Display		TFT-Display
Multiple languages		yes
Right and left movement/ interval operation		yes
Smooth start		yes
Lift		automatic
Lifting speed	<b>mm/s</b>	50
Stroke	<b>mm</b>	140
Setting of lower end stop	<b>mm</b>	60, contactless
Head angle adjustable		0°...45°
Timer	<b>hh:mm:ss</b>	00:00:00...99:59:59
Interval	<b>mm:ss</b>	00:00...60:00
Integrated vacuum controller		yes
Measurement range vacuum	<b>mbar (hPa)</b>	1...1050
Adjustment range vacuum	<b>mbar (hPa)</b>	2...1014
Measurement accuracy vacuum	<b>mbar</b>	± 2 (± 1 Digit) if calibration is done at a constant temperature
	<b>mbar</b>	adjustable hysteresis (Two-position vacuum control)
	<b>%</b>	1 (Speed-vacuum control, set value ≥ 100 mbar)
Measurement range of differential temperature	<b>K</b>	7
Offset balancing range	<b>K</b>	± 0.5
Interface		USB, RS 232
Ramp programming		yes
Distillation processes programmable		yes
Remote control		with accessory <i>labworldsoft</i> ®
Cooling surface (Standard condenser RV 10.1, 10.10, 10.2, 10.20)	<b>cm<sup>2</sup></b>	1500
Minimum cooling water flow rate	<b>l/h</b>	30
Maximum cooling water flow rate	<b>l/h</b>	100
Cooling water pressure	<b>bar</b>	1
Cooling water temperature range	<b>°C</b>	18...22, constantly
Fuse		2 x T1, 6 A 250 V 5 x 20
Perm. On-time	<b>%</b>	100
Perm. ambient temperature	<b>°C</b>	5...40
Perm. humidity	<b>%</b>	80
Protection acc. to DIN EN 60529		IP 20
Protection class		I
Overvoltage category		II
Contamination level		2
Weight (no glassware; no heating bath)	<b>kg</b>	18.6
Dimensions (W x D x H)	<b>mm</b>	500 x 440 x 430
Operation at a terrestrial altitude	<b>m</b>	max. 2000 above sea level

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			Réglage d'usine	
Evaporator	Modes	Automatic	Bath medium.....	Water
			Start temperature.....	60 °C
		Manual.....		1013 mbar
		Pump %.....		50%
		Volume	Solvent.....	Acetic acid
			Target.....	100 ml
			Efficiency.....	80%
			Adjustment.....	-
		Program.....		-
		100%	Solvent.....	Acetic acid
		Drying.....	-	
	Cleaning	Start after.....	mm:ss.....	30:00
		Duration.....	mm:ss.....	06:00
		Pump speed.....	Value.....	20%
		Clean now.....		-
	Limits	Manual	Minimum.....	2 mbar
			Maximum.....	1013 mbar
		Pump %	Minimum.....	0%
			Maximum.....	100%
	Hysteresis	Automatic.....		✓
		Manual	Value.....	±15 mbar
	Solvent library	Add solvent		
		Add solvent		
		Add solvent	Name.....	-
		Add solvent	Formula.....	-
		Add solvent	Bath temperature.....	-
			Rotation speed.....	-
			Boiling point.....	-
			Efficiency.....	-
			Heat capacity.....	-
			Enthalpy.....	-
			Density.....	-
			i-factor.....	-
		Acetic acid		
		Acetone		
		Acetonitrile	Name.....	-
		...	Formula.....	-
			Bath temperature.....	-
			Rotation speed.....	-
			Boiling point.....	-
		Efficiency.....	-	
	Heat capacity.....	-		
	Enthalpy.....	-		
	Density.....	-		
	i-factor.....	-		
Calculator	Name.....	Acetic acid		
	Vapor temperature.....	40 °C		
	Heating bath temperature.....	60 °C		
	Pressure.....	46 mbar		
Actions after distillation	Stop heating.....		-	
	Stop rotation.....		✓	
	Lift up.....		✓	
	Close water valve.....		✓	
	Open venting valve.....		✓	
	Cleaning.....		-	

			Réglage d'usine		
Rotation	Advanced	Flow control.....		✓	
		Bath medium.....	Water		
		Cooler power limit.....	900 w		
		Max power of cooler.....	0 w		
	Interval mode	Interval	Ccw Value.....	00:10	
			Cw Value.....	00:10	
	Speed limit	Activate.....		-	
		Minimum.....	5 rpm		
		Maximum.....	300 rpm		
	Timer	Timer function	Time.....	hh:mm:ss.....	00:00:00
Beep after timeout.....				-	
Activate.....				-	
Display.....			-		
Display	Timer.....		✓		
	Flow rate.....		✓		
	Δ T.....		✓		
	Cooler power.....		✓		
Programs	Program 1-10	Select.....		-	
		Edit	Edit	Pressure.....	-
				Rotation (Speed).....	-
				hh:mm:ss.....	-
		Insert.....		-	
		Delete.....		-	
	Last Measurement	Save as.....		-	
		Edit	Edit	Pressure.....	-
				Rotation (Speed).....	-
				hh:mm:ss.....	-
		Insert.....		-	
		Delete.....		-	
Delete.....		-			
Safety	Password	Value.....		0	
		Value.....		0	
		Value.....		0	
	Continue after power failure	Activate.....		-	
Time.....		mm:ss.....	00:10		

		Réglage d'usine		
Service	Valves	Venting valve.....	-	
		Vacuum valve.....	-	
		Water valve.....	-	
	Pump.....	-		
Temperature adjustment.....	-			
Settings	Language	English.....	✓	
		Deutsch.....	-	
		Français.....	-	
		Español.....	-	
		Italiano.....	-	
		Português.....	-	
		Русский язык.....	-	
		中文.....	-	
		한국어.....	-	
		日本語.....	-	
	Units	Pressure	mbar.....	✓
			hPa.....	-
			mmHg.....	-
			Torr.....	-
		Amount	mL.....	✓
			g.....	-
	Display	Background	Black.....	✓
			White.....	-
		Brightness Value.....	100%	
	Firmware update info.....	-		
	Sound	Volume Value.....	100%	
		Key Tone.....	-	
	Factory settings	Menu values.....	-	
		Programs.....	-	
		Solvent library.....	-	
		All.....	-	
	Communication	Device name.....	RV 10 auto	
		Bluetooth.....	✓	
		Labworldsoft 5 protocol.....	✓	
	Information	Firmware update info.....	-	
		Display version.....	-	
		Logic version.....	-	
Pressure max.....		1100 mbar		
Pressure min.....		1 mbar		
Pump % max.....		100%		
Pump % min.....		0%		
Rot. speed max.....		300 rpm		
Rot. speed min.....		5 rpm		
Operating hours.....		0 hours		

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			Configuração de fábrica	
Evaporator	Modes	Automatic	Bath medium.....	Water
			Start temperature.....	60 °C
		Manual.....		1013 mbar
		Pump %.....		50%
		Volume	Solvent.....	Acetic acid
			Target.....	100 ml
			Efficiency.....	80%
			Adjustment.....	-
		Program.....		-
		100%	Solvent.....	Acetic acid
		Drying.....	-	
	Cleaning	Start after.....	mm:ss.....	30:00
		Duration.....	mm:ss.....	06:00
		Pump speed.....	Value.....	20%
		Clean now.....		-
	Limits	Manual	Minimum.....	2 mbar
			Maximum.....	1013 mbar
		Pump %	Minimum.....	0%
			Maximum.....	100%
	Hysteresis	Automatic.....		✓
		Manual	Value.....	±15 mbar
	Solvent library	Add solvent		
		Add solvent		
		Add solvent	Name.....	-
		Add solvent	Formula.....	-
		Add solvent	Bath temperature.....	-
			Rotation speed.....	-
			Boiling point.....	-
			Efficiency.....	-
			Heat capacity.....	-
			Enthalpy.....	-
			Density.....	-
			i-factor.....	-
		Acetic acid		
		Acetone		
		Acetonitrile	Name.....	-
		...	Formula.....	-
		Xylene	Bath temperature.....	-
			Rotation speed.....	-
			Boiling point.....	-
		Efficiency.....	-	
	Heat capacity.....	-		
	Enthalpy.....	-		
	Density.....	-		
	i-factor.....	-		
Calculator	Name.....	Acetic acid		
	Vapor temperature.....	40 °C		
	Heating bath temperature.....	60 °C		
	Pressure.....	46 mbar		
Actions after distillation	Stop heating.....		-	
	Stop rotation.....		✓	
	Lift up.....		✓	
	Close water valve.....		✓	
	Open venting valve.....		✓	
	Cleaning.....		-	

			Configuração de fábrica		
Advanced	Flow control.....		✓		
	Bath medium.....		Water		
	Cooler power limit.....		900 w		
	Max power of cooler.....		0 w		
Rotation	Interval mode	Interval	Ccw Value.....	00:10	
			Cw Value.....	00:10	
		Activate.....		-	
	Speed limit	Minimum.....		5 rpm	
		Maximum.....		300 rpm	
Timer	Timer function	Time.....	hh:mm:ss.....	00:00:00	
		Beep after timeout.....		-	
		Activate.....		-	
	Display.....			-	
Display	Timer.....		✓		
	Flow rate.....		✓		
	Δ T.....		✓		
	Cooler power.....		✓		
Programs	Program 1-10	Select.....		-	
		Edit	Edit	Pressure.....	-
				Rotation (Speed).....	-
				hh:mm:ss.....	-
		Insert.....		-	
		Delete.....		-	
	Last Measurement	Save as.....		-	
		Edit	Edit	Pressure.....	-
				Rotation (Speed).....	-
				hh:mm:ss.....	-
			Insert.....		-
			Delete.....		-
Delete.....			-		
Safety	Password	Value.....		0	
		Value.....		0	
		Value.....		0	
	Continue after power failure	Activate.....		-	
		Time.....	mm:ss.....	00:10	

		Configuração de fábrica		
Service	Valves	Venting valve.....	-	
		Vacuum valve.....	-	
		Water valve.....	-	
		Pump.....	-	
Settings	Language	English.....	✓	
		Deutsch.....	-	
		Français.....	-	
		Español.....	-	
		Italiano.....	-	
		Português.....	-	
		Русский язык.....	-	
		中文.....	-	
		한국어.....	-	
		日本語.....	-	
		Units	Pressure	mbar.....
	hPa.....			-
	mmHg.....			-
	Torr.....			-
	Amount		mL.....	✓
			g.....	-
	Display	Background	Black.....	✓
			White.....	-
		Brightness Value.....	100%	
	Firmware update info.....	-		
	Sound	Volume Value.....	100%	
		Key Tone.....	-	
	Factory settings	Menu values.....	-	
		Programs.....	-	
		Solvent library.....	-	
		All.....	-	
Communication	Device name.....	RV 10 auto		
	Bluetooth.....	✓		
	Labworldsoft 5 protocol.....	✓		
Information	Firmware update info.....	-		
	Display version.....	-		
	Logic version.....	-		
	Pressure max.....	1100 mbar		
	Pressure min.....	1 mbar		
	Pump % max.....	100%		
	Pump % min.....	0%		
	Rot. speed max.....	300 rpm		
	Rot. speed min.....	5 rpm		
	Operating hours.....	0 hours		

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