IKA

designed for scientists

VC 10 lite
VC 10 pro

| | | |
|---|----|-----|
| Betriebsanleitung Ursprungssprache | DE | 5 |
| Operating instructions Mode d'emploi | EN | 29 |
| Руководство пользователя | FR | 53 |
| 使用说明 | RU | 77 |
| | ZH | 101 |

| | | | | | |
|------------------------------------|----|-----|-----------------------------|----|-----|
| Breve guía | ES | 124 | Ohutusjuhised | ET | 142 |
| Brevi istruzioni | IT | 125 | Drošības norādījumi | LV | 144 |
| Veiligheidsaanwijzingen | NL | 127 | Saugos nurodymai | LT | 145 |
| Säkerhetsanvisningar | SV | 128 | Инструкции за безопасност | BG | 146 |
| Sikkerhedshenvisninger | DA | 130 | Indicații de siguranță | RO | 148 |
| Sikkerhetsinformasjon | NO | 131 | Υποδείξεις ασφαλείας | EL | 150 |
| Turvallisuusohjeet | FI | 132 | Güvenlik talimatları | TR | 151 |
| Instruções de segurança | PT | 134 | 안전 지침 | KO | 153 |
| Wskazówki dotyczące bezpieczeństwa | PL | 135 | 安全上のご注意 | JA | 154 |
| Bezpečnostní pokyny a informace | CS | 137 | คำแนะนำเกี่ยวกับความปลอดภัย | TH | 155 |
| Biztonsági utasítások | HU | 138 | Petunjuk keselamatan | ID | 157 |
| Varnostna navodila | SL | 139 | | | |
| Bezpečnostné pokyny | SK | 141 | | | |

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Fig. 1

| Pos. | Bezeichnung | (DE) | Item | Designation | (EN) |
|------|--|------|------|--|------|
| A | Display | | A | Display | |
| B | Back-Taste | | B | Back key | |
| C | Menu-Taste | | C | Menu key | |
| D | Dreh-/Druckknopf | | D | Rotating/pressing knob | |
| E | Pause-Taste | | E | Pause key | |
| F | Ein/Aus Taste | | F | On/Off key | |
| G | Tastensperre | | G | Lock key | |
| H | Belüftungs-Taste | | H | Venting key | |
| I | Belüftungs-Anschluss | | I | Venting connection | |
| J | Rezipient-Anschluss | | J | Recipient connection | |
| K | Pumpen-Anschluss | | K | Pump connection | |
| L | USB-Schnittstelle | | L | USB connection | |
| M | RS 232 Schnittstelle | | M | RS 232 connection | |
| N | Stativ-Feststellschraube | | N | Stand-locking screw | |
| O | PT 1000 Temperaturfühler-Schnittstelle (VC 10 pro) | | O | PT 1000 temperature sensor interface (VC 10 pro) | |
| P | VACSTAR Pumpen-Schnittstelle | | P | VACSTAR pump interface | |
| Q | Anschlussbuchse Netzteil 24 W | | Q | Connection socket, power supply 24 W | |
| R | Lüftungsschlitze | | R | Ventilation slots | |

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Contents

| | Page |
|--------------------------------|------|
| Device setup | 2 |
| Declaration of conformity | 29 |
| Explication of warning symbols | 29 |
| Warranty | 29 |
| Safety instructions | 30 |
| Unpacking | 31 |
| Intended use | 31 |
| Useful information | 31 |
| Setting up | 33 |
| Commissioning | 36 |
| Interfaces and outputs | 46 |
| Maintenance and cleaning | 47 |
| Error codes | 48 |
| Accessories | 49 |
| Product contact parts | 49 |
| Technical data | 50 |

Declaration of conformity

We declare under our sole responsibility that this product corresponds to the regulations 2014/35/EG, 2014/30/EG and 2011/65/EG and conforms with the standards or standardized documents EN 61010-1, EN 60529, EN 61326-1 and EN ISO 12100.

A copy of the complete EU Declaration of Conformity or further declarations of conformity can be requested at sales@ika.com.

Explication of warning symbols

Indicates an (extremely) hazardous situation, which, if not avoided, will result in death, serious injury.

Indicates a potentially hazardous situation, which, if not avoided, can result in death, serious injury.

Indicates a potentially hazardous situation, which, if not avoided, can result in injury.

Indicates practices which, if not avoided, can result in equipment damage.

Warranty

In accordance with warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual.

Safety instructions

For your protection

NOTICE

Read the operating instructions in full before starting up and follow the safety instructions.

- Keep the operating instructions in a place where it can be accessed by everyone
- Ensure that only trained staff work with the instrument.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.

DANGER

Always wear personal protective equipment in accordance with the hazard class of the media being worked with, otherwise there is a risk of splashing liquids, release of toxic or inflammable gases!

- The device and device parts must not be used on humans or animals.
- Do not expose parts of the human body or animals to the vacuum.
- Position the device as specified in the "Commissioning" section and connect the connection cables and interfaces as described.
- Never work with an incorrectly connected or defective device.
- Inhalation of or contact with media such as poisonous liquids, gases, spray mist, vapours, dusts or biological and microbiological materials can be hazardous to user. Make sure that all connections are tight and leak-free when working with such media.
- Prevent the release of the above mentioned substances. The use of suitable extraction systems e.g. fume cupboards is recommended in the place of installation. Take protective measures for personnel and the environment.
- Due to the residual leakage rate of the device, this can result in the release of medium.
- The VC 10 lite/pro vacuum controller is not designed to be installed in potentially explosive atmospheres
- The VC 10 lite/pro vacuum controller is not designed for use with self inflammable substances, substances which are inflammable without air or explosive substances.
- Take suitable precautions to prevent any formation of explosive mixtures. Use inert gas for venting and/or dilution, if necessary.
- Take into consideration interactions or possible chemical or physical reactions when working with media under reduced pressure and elevated temperature.
- Never work with the device under water or underground.
- The device must only be operated with supervision.
- Electrostatic events can occur between the medium and device and lead to a direct hazard.
- Safe operation is only guaranteed with the accessories described in the "Accessories" section.
- The instrument can only be disconnected from the mains supply by pulling out the mains plug or the connector plug.
- The socket for the mains cord must be easily accessible.
- Do not use the device if the AC adapter is defective.
- Follow the operating instructions of the additional equipment carefully (e.g. a rotary evaporator, vacuum pump) with which the VC 10 lite/pro vacuum controller is operated.

- The VC 10 lite/pro vacuum controller may only be operated under the conditions described in the "Technical data" chapter.
- The pressure at the gas inlet and gas outlet must not exceed 1100 mbar. For pressure over 1100 mbar, the total allowable pressure is exceeded and the device no longer displays the pressure accurately.
- Only use flexible hose lines.
- Flexible elements can be compressed under vacuum.
- In the event of a power failure, the integrated venting valve automatically aerates the connected recipient.
- Be aware of your emergency measures in the event of a power failure and ensure that the system is put into a safe state.
- The safety of the user cannot be guaranteed if the appliance is operated with accessories that are not supplied or recommended by the manufacturer or if the appliance is operated improperly, contrary to the manufacturer's specifications or if the instrument or the printed circuit board are modified by third parties.

For protection of the equipment

CAUTION

The voltage stated on the type plate must correspond to the mains-voltage.

- Covers or parts that can be removed from the device without tools must later be refitted to ensure safe operation, so long as no other connection is made at this point. This will prevent the infiltration of foreign objects, liquids etc.
- The device must only be operated with the original plug-in power supply unit supplied.
- Protect the appliance and accessories from bumps and impacts.
- The appliance may only be opened by experts (authorized staff).
- To ensure sufficient cooling of the VC 10 lite/pro vacuum controller, the ventilation slots on the housing must not be covered.
- Use only original spare parts for repair and maintenance in order to ensure the reliable operating condition of the device.
- Pay attention to water condensation inside and outside of the device. If the device is brought from a cold environment, allow the device to warm up.
- Never attach the vacuum controller above a heating bath.
- Ensure that no solids and/or liquids enter the VC 10 lite/pro vacuum controller via the suction power of the pump. This will result in damage to the pressure sensor and the valves.

WARNING

Check the compatibility of the substances used with the device with the materials listed under "Product contact parts"; see chapter "Product contact parts".

- Place the VC 10 lite/pro vacuum controller in an elevated position if possible, so that in the event of a fault, no condensate can collect on the pressure sensor.



Unpacking

Unpacking

- Please unpack the device carefully.
- In the case of any damage a detailed report must be sent immediately (post, rail or forwarder)

Delivery scope

| | |
|-------------------------------------|---|
| Vacuum controller VC 10 lite/pro | |
| Universal plug-in power supply unit |  Fig. 2 |
| USB A - micro B cable |  Fig. 3 |

| | |
|------------------------|---|
| Vacuum hose 1 m |  Fig. 4 |
| Y connector |  Fig. 5 |
| Operating instructions | |
| Warranty card | |

Intended use

Use

Together with the accessories recommended by , the VC 10 lite/pro vacuum controller is suitable for the controlled evacuation of air (gas) from laboratory devices and equipment (e.g. rotary evaporator or laboratory reactors, and also for classical separation or filtration or drying tasks in the laboratory) in conjunction with a suitable vacuum source (e.g. MVP 10 basic / VACSTAR). The VACSTAR can be operated automatically with the VC 10 pro

vacuum controller. Therefore, it can be used for automatic boiling point recognition, programmed pressure-time curves or programs from the solvent library, for example.

Mode of operation: Stand device

Area of use

Indoor environments similar to that a laboratory of research, teaching, trade or industry area.

The safety of the user cannot be guaranteed:

- If the device is operated with accessories that are not supplied or recommended by the manufacturer;
- If the device is operated improperly or contrary to the manufacturer's specifications;
- If the device or the printed circuit board are modified by the third parties.

WARNING

The device must not be used for:
- **vacuation of biological biotopes,**
- **evacuation of explosive, corrosive or gas volumes,**
- **throughput/use of fluids!**

Useful information

In addition to a high-resolution Al₂O₃ pressure sensor, the VC 10 lite/pro vacuum controller has an integrated vacuum valve and a venting valve. As a result, the VC 10 lite/pro vacuum controller can be used without any additional peripheral devices.

The VC 10 pro provides two basic modes of operation: Two-position control and analogue speed control.

The change between the two modes is performed automatically by plugging or unplugging the connector cable at the "VACSTAR" interface. Depending on the selection made, some submenu items may not be available (greyed out on the display).

With two-position control, parameters from the solvent library can be applied or target values specified in manual mode.

With analogue speed control, an automatic mode as well as the setting of the percentage pump speed can also be preset. Space requirement: To connect all peripheral devices you will need approximately WxDxH 150 x 200 x 155 mm³ space for the device with connected plugs.

Two-position control VC 10 pro (only)

A vacuum is made in the glassware with the help of a vacuum pump. The vacuum pump operates at constant speed which generally cannot be adjusted. Once the target value setting is reached, the integrated vacuum valve interrupts the suction line.

The device will fall slightly short of the target value setting due to the time difference in the millisecond range with the vacuum valve, the target value comparison and the switching of the vacuum valve. Natural leakage from the apparatus now causes another increase in the pressure value of the system, which in turn is detected by the integrated pressure sensor.

If the target value is exceeded, the vacuum valve is reopened and the vacuum pump in operation lowers the pressure. The pressure curve shown on the display fluctuates around the set target value. The pressure difference between activating/deactivating the valve can be set using the vacuum hysteresis value.

The quality of the control depends on the suction power of the pump (rotation speed), the tightness of the system and signal delay time of the electronics used.

Automatic boiling point recognition is not possible when using two-position control.

Schematic view of a two-position vacuum control

Analogue speed control

Imprecision of the two-position control is avoided by using the analogue speed control.

With this kind of control, the speed of the pump, and therefore its suction power, is reduced the closer the measured pressure reaches the target value. Once the target value is reached, the pump only operates according to the leakage rate of the system. This type of control offers quieter operation, more exact vacuum control and a long service life of movable pump parts.

Automatic boiling point recognition is possible with this type of control, i.e. the system reaches and keeps the solvent at its boiling point in automatic mode.


Analogue speed control with the VACSTAR:

The VACSTAR is automatically detected as soon as the connection cable is plugged into both devices. The analogue speed control is enabled, the two-position control is disabled and cannot be selected.

Menu items such as automatic boiling point recognition, programmed pressure-time curves or programs from the solvent library can be used.

Refer to the operating instructions for the VACSTAR vacuum pump.

Schematic view of speed-vacuum control

*An external temperature sensor can be connected to the "PT1000" temperature sensor interface. When the menu view is enabled, the temperature value is shown in the display.

Time sequences can be controlled using a timer.

The operating modes A, B, C and D have preset access options for adjusting and saving target values.

Display settings can be changed using the "Display" menu item.

Access to the device can be restricted by setting up a password in the "Security" menu item.

Common settings such as language, pressure and temperature units can be set via the "Settings" menu item. The display background can also be changed from black to white as well as brightness settings adjusted.

The integrated solvent library calculates the theoretical boiling point at a specified steam temperature value, or vice versa.

The upper and lower limits or the switching points of the vacuum control can be set in the submenu items "Limit values" and "Hysteresis".

*only VC 10 pro

Setting up

NOTICE

Observe the general information, always connect the recipients (load/vacuum vessel/glass cooler) at their highest points to the suction line. This will prevent the risk of fluid entering into the vacuum controller or pump.

Install a separator (e.g. Woulff bottle) in front of the intake mounting to protect against the ingress of liquid.

In the case of induced solvent vapours, a vapour condenser (accessory of pump) installed downstream of the pump helps to condense the vapours and prevents release into the atmosphere. Liquid in the pump chambers deteriorates the pump characteristics.

Attach the VC 10 lite/pro vacuum controller to a stand (d=16 mm) close to the pump.

Secure the controller with the stand locking screw (N).


Fig. 5

To do so, use a stand with diameter 16 mm, for example, stand VC 10.400, VSS 1 vacuum safety set, the stand on VACSTAR vacuum pump or attach the controller to the RV 8 rotary evaporator, or at LR 1000 laboratory reactor.

Fig. 6

VC 10 lite/pro at stand VC 10.400

Fig. 7

VC 10 lite/pro to VSS 1 with MVP 10 basic

Fig. 8

VC 10 lite/pro to VACSTAR

Fig. 9

VC 10 lite/pro to RV 8

Connecting the interfaces

Connect the vacuum hoses, pay attention to the following connections (also see **Fig. 1**):

J: Load

Hose connection d=10 mm to the system to be evacuated (recipient).

K: Pump

Hose connection d=10 mm to the vacuum source (e.g. pump).

L: Vent

Ventilation connection d=10 mm, also suitable for the connection of inert gas.

Secure the hoses with hose clamps, if necessary.

Lay the vacuum hoses according to your device configuration.

Now connect the connecting cables according to your system configuration:

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Connect the VC 10 lite/pro vacuum controller to a PC using the USB 2.0 cable (A/Micro-B).

Software updates to any existing equipment can be installed using the FUT software tool.

M: RS 232 Interface

You can connect the VC 10 lite/pro vacuum controller to a PC using an RS 232 interface cable. The pump can be operated in conjunction with other devices using labworldsoft® laboratory software. For more information see Chapter "Interfaces and outputs".

P: PT 1000 temperature sensor (VC 10 pro)

Connect the optional PT 1000.60, PT 1000.61 or PT 1000.70 temperature sensor (accessory). The measured temperature (e.g. tempering bath) is shown on the display.

The value displayed is for information purposes only and is not processed logically.

Note that the value displayed is dependent on the positioning of the sensor in the medium, its mixing, as well as the calibration of the probe and can deviate from other measured temperature values.

Incorrect or fluctuating measured values can be displayed by using a heating bath, e.g. by laying the sensor on the heating bath floor, instead of positioning in the medium.

Q: VACSTAR membrane vacuum pump (Mini DIN)

You can connect the vacuum controller and the VACSTAR vacuum pump to the analogue connection cable (accessory) for exact speed-controlled vacuum control. The vacuum controller detects the pump and switches to speed-vacuum control mode. The two-position control is disabled. The pump speed is regulated depending on the measured pressure.

R: Connector for 24 W AC adapter

Connect the universal AC adapter (included) to the socket on the rear. Observe the country-specific variants of the AC adapter!

Explanation of abbreviations used in the following schematic drawings and configurations:

- 1 Recipient (load, e.g. rotary evaporator, reactor)
- 2 Separator (Woulff bottle)
- 3 Venting valve
- 4 Pressure sensor
- 5 Vacuum valve/ball valve
- 7 Vapour condenser
- 8 Connection cable (analogue)
- 9 Vacuum source (pump, in-house vacuum)
- 10 VC 10 lite/pro.300 Check valve

Two-position control

Vacuum pump/vacuum source with VC 10 pro vacuum controller.

Once the target value is reached, the suction line is closed by the vacuum valve integrated in the VC 10 pro vacuum controller. The pump speed can be changed by manually adjusting the speed setting.

For pumps which have the speed setting option, it is recommended to work at low speed.


NOTICE

When using the VACSTAR vacuum pump: The pump must be operated in mode "B" (restart after voltage drop).

The vacuum source can also be a central in-house vacuum network.

Illustration two-position control - vacuum pump/vacuum source with vacuum controller VC 10 pro

Vacuum source with two VC 10 pro vacuum controllers.


Once the target value is reached in a recipient, the suction line is closed by the vacuum valve integrated in the VC 10 pro vacuum controller.

A sufficiently stable vacuum source (in-house vacuum, pump) can serve multiple consumers.

For pumps which have the speed setting option, it is recommended to work at low speed.

A second consumer can be controlled independent of the first using a second vacuum controller.

Illustration two-position control - vacuum source with two vacuum controllers VC 10 pro

Speed control


VACSTAR vacuum pump with VC 10 lite/pro vacuum controller

Automatic adjustment of the “analog speed control” mode if the analogue connection cable (Pos. 8) is connected to the pump with the vacuum controller.

Once the target value is reached the pump speed is automatically reduced to a minimum “0 rpm”. Once the target value is reached,

the pump operates according to the leakage rate of the vacuum system. Automatic boiling point recognition is only possible with the VC 10 pro .

Illustration speed control - Vacuum pump VACSTAR with vacuum Controller VC 10 lite/pro

Commissioning

Check whether the voltage given on the type plate corresponds to the available mains voltage.

If these conditions have been met, the machine is ready for operation when the mains plug is plugged in.

If these conditions are not met, safe operation is not guaranteed and the machine could be damaged.

Also pay attention to the ambient conditions listed in the “Technical data”.

The start screen appears for a few seconds when the device is switched using the “Standby ON/OFF” switch (F) The device name and the software/firmware versions are also displayed.


An information prompt then appears to download the firmware update tool.

The following working screen appears automatically in the display.


Display of two-position control, e.g. with MVP 10 basic.

Hysteresis display for two-position control


Display speed control: VACSTAR to analogue connection

No hysteresis display for speed control

Explanation of symbols on the working screen

The symbols displayed change depending on the status and settings of the vacuum controller:

Lock key

This symbol means that the function of the key and of the rotating/pressing knob for controlling the vacuum controller are disabled. The symbol no longer appears if the functions are enabled once again by pressing the lock key a second time (min. 1 sec.).

Operating mode

This symbol indicates the operating mode currently selected (A, B, C, D).

USB

This symbol means the vacuum controller is communicating via a USB cable. The symbol no longer appears if no USB cable is being used for communicating with the station.

Temperature sensor (VC 10 pro)

This symbol appears when the temperature display is active in the display.

Timer

This symbol appears when the timer display is active in the display.

PC control

This symbol means that the vacuum controller is connected to a computer and the device is being controlled from the computer.

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This symbol means that the vacuum controller is connected to a computer and the device is being controlled from the selected program.

Continuous mode

This symbol indicates that the vacuum controller is in continuous mode.

Vent

This symbol indicates that venting is enabled.

Vacuum

This symbol indicates that the current pressure is less than the existing atmospheric pressure.

Detection

This symbol is displayed during the automatic boiling point search, the system pressure is reduced further.

Boiling

This symbol indicates that a boiling point has been found automatically, the system pressure is kept constant.

Refrigerated circulator

This symbol indicates that refrigerated circulator is operating.

Navigation menu

Fig. 9


- Press the "Menu" key (C).
- Select the menu by turning the rotating/pressing knob (D) to the right or left to select the desired menu or sub-menu, which can then be selected by pressing the rotating/pressing knob (D).
- Press or turn the rotating/pressing knob (D) again to select the desired menu option and edit the values or settings, or activate/deactivate a function.
- Turn the rotating/pressing knob (D) to "OK" and press the "Back" (B) or "Menu" key (C) to end the procedure and return to the previous menu or working screen.
- Navigating in Menu "Graph": Press the "Back" (B) at main screen to navigate to graph menu, showing real-time pressure change on graph.

NOTICE

The menu option activated is highlighted yellow in the display. Menus which are grayed out are inactive.

Navigation menu: Press the "Menu" key (C) and turn the rotating/pressing knob (D)
 Press the "Back" key (B) or the "Menu" (C)

NOTICE


If you press the "Menu" key (C), the system skips directly back to the working screen.
 If you press the "Back" key (B), the system skips back to the previous display.

Menu structure


| Menu | Sub menu | Option | Action | Factory setting | |
|-----------------------------|---------------------------------|------------------------------|----------------------|--------------------------|-------------|
| Pump | Modes | Automatic (only VC 10 pro) | | Deactivated | |
| | | Manual | | Activated | |
| | | Pump % | | Deactivated | |
| | | Program | | Deactivated | |
| | | Start After | | - | |
| | Refrigerated circulator | Temperature | | - | |
| | | Pump speed | | - | |
| | | Cool now | | - | |
| | | Stop now | | - | |
| | | Start after | | - | |
| | Cleaning | Start after | | 05:00 [mm:ss] | |
| | | Duration | | 06:00 [mm:ss] | |
| | | Pump speed | | 20 % | |
| | | Clean now | | Deactivated | |
| | Limits | Manual | Minimum | | 1 mbar |
| | | | Maximum | | 1100 mbar |
| | | Pump % | Minimum | | 0 % |
| | | | Maximum | | 100 % |
| | Hysteresis (VC 10 pro) | Automatic | | Activated | |
| | | Manual | | Deactivated | |
| Solvent library (VC 10 pro) | Acetonitril | | - | | |
| | N-Amylalkohol | | - | | |
| | n-Pentalon | | - | | |
| | Benzol | | - | | |
| | ... | | - | | |
| Actions after stop | Open venting valve | | Activated | | |
| | Cleaning | | Deactivated | | |
| Actions after start | Stop refrigerated circulator | | Deactivated | | |
| | Turn on refrigerated circulator | | Deactivated | | |
| Leakage Test | Leakage rate | | 3 mbar | | |
| | Tolerance | | 1mbar | | |
| | Test point | | 500 mbar | | |
| | Leakage test duration | | 1 minute | | |
| | | | | - | |
| Temperature (VC 10 pro) | Sensor temperature | | - | | |
| | Display | | Deactivated | | |
| | Calibration | | 20.0 °C | | |
| Timer | Timer function | Time | | 00:00:01 [h:mm:ss] | |
| | | Beep after stop | | Activated | |
| | Display | | Deactivated | | |
| Operating mode | A | | Deactivated | | |
| | B | | Activated | | |
| | C | | Deactivated | | |
| | D | | Deactivated | | |
| Display | Timer | | Deactivated | | |
| | Temperature (VC 10 pro) | | Deactivated | | |
| Programs | Program 01 | | - | | |
| | ... | | - | | |
| | Program 10 | Select | | - | |
| | | Edit | | Edit, insert, delete | |
| | | Delete | | - | |
| Last measurement | Save as | | - | | |
| | Edit | | Edit, insert, delete | | |
| | Delete | | - | | |
| Safety | Password | | 0 0 0 | | |
| Service | Valves | Venting valve | | Open | |
| | | Vacuum valve (VC 10 pro) | | Close, open | |
| | Pump | | Close | | |
| Refrigerated circulator | | Deactivated | | | |
| Settings | Languages | English | | Activated | |
| | | German | | Deactivated | |
| | Units | Pressure | mbar | | Activated |
| | | | hPa | | Deactivated |
| | | | mmHg | | Deactivated |
| | | | Torr | | Deactivated |
| | | Temperature (VC 10 pro) | °C | | Activated |
| | | °F (not available for Japan) | | Deactivated | |
| | Display | Background | Black | | Activated |
| | | | White | | Deactivated |
| | Brightness | | 100 % | | |
| | Firmware update info | | Activated | | |
| | Sound | Volume | | 10 % | |
| | | Key tone | | Deactivated | |
| | Factory setting | Menu values | | Reset to factory setting | |
| Programs | | | - | | |
| Communication | All | | - | | |
| | Device name | | VC 10 lite/pro | | |
| | Firmware update info | | - | | |
| | Display version | | - | | |
| | Logic version | | - | | |
| | Operating mode | | B | | |
| | Pressure max | | 1100 mbar | | |
| | Pressure min | | 1 mbar | | |
| | Pump % max | | 100 % | | |
| | Pump % min | | 0 % | | |

Menu details

Pump

Modes

Automatic (VC 10 pro):


In the "Modes" menu, the user can enable automatic boiling point recognition by selecting the "Automatic" menu item. No other parameters must be set.

The boiling point is detected automatically. For rotary evaporators with heating bath, it must be ensured that the tempering medium and solvents have a constant temperature (e.g. 60 °C).

Manual:

In the "Manual" menu item, the target value can be specified (e.g. in "mbar"). The system is evacuated until the target value.

Refrigerated circulator


Start after...:

Specifies a delay before the refrigerated circulator start operating.

Temperature:

Determines the temperature, that the refrigerated circulator is operated.

Cleaning

The cleaning of the pump ensures that solvent residues are removed in the pump after the measurement. For this purpose, the pump is operated at a low level while the system is vented. The air flowing through the pump ensures that the solvent residues are discharged from the pump. If the cleaning option is selected, this starts once the measurement has finished (see "Actions after stop").

Pump %

By selecting the "Pump %" menu item, the pump can be operated continuously with a running performance of between 100 % and 1 %.

Program

Under the "Programs" menu, 10 user-defined pressure-time profiles can be created. The last measurement taken can also be viewed in this menu and saved as a program.

The target value and the running performance can only be entered within the set limits (see menu item "Limits").

Changing to the working screen and pressing the rotating/pressing knob starts the process. Pressing the knob again stops the process.

NOTICE

If no speed-controlled pump with analogue connection is connected to the "VACSTAR" interface, the system can only work in manual mode and in "two-position control" mode. The menu items "Auto" and "Pump %" are greyed out and cannot be activated.

Pump speed:

Determines the speed in rpm, that the refrigerated circulator is operated

Cool now:

If you want to perform a cooling immediately, you can select it with this setting.

Stop now:

If you want to stop the cooling immediately, you can select it with this setting.

NOTICE

The "Cleaning" menu item increases the service life of the pump, since solvent residues do not remain inside the pump.

You can also clean the pump without using the corresponding menu item. To do this, remove the hoses and operate the pump at the end of a working cycle (test, test series, end of the day ...) at idling speed.

Start after...:

Specifies after which minimum duration of a measurement the cleaning should take place. The default time is 5 minutes. This means that a measurement must be run for at least 5 minutes, so that cleaning is performed at the end of the measurement.

Duration:

Specifies how long the cleaning process should be performed.


Pump speed:

Determines at what power (in %) the pump is operated during cleaning. If the power is set too high this can lead to a pressure drop in the system.

Clean now:


If you want to perform a cleaning, you can select it with this setting.

Limits

In this menu, limits for the target value can be set.

Hysteresis (VC 10 pro)

With the "Hysteresis" option, the switching frequency and control accuracy can be influenced when in a steady state. The hysteresis value determines the upper and lower limits of the target value and therefore the closing and opening of the pump and valve (only for two-position control).

Automatic: Hysteresis is always 10 % of the actual pressure.


Manual: Manual value specification of hysteresis.

NOTICE


Limits which are too narrow shorten the service life of the devices.

Solvent Library (VC 10 pro)


In the "solvent library", the most common solvents are assigned a function to determine the steam temperature at the set pressure and vice versa. This provides assistance when setting these parameters in manual mode, for example, evaporation processes with a rotary evaporator.

Actions after stop


This menu allows the user to define actions that should be performed when a test has ended.

Actions after start

This menu allows the user to define actions that should be performed when a test has started.

Leakage test

Leakage test: User to confirm the vacuum connection. Press ok if the connection is good.

Set Value: User to input the required leakage rate with tolerance. Default leakage rate is 3mbar with tolerance 1mbar.


Test point: User to input the pressure which the leakage test will

be executed at. Default value is 500 mbar.

Duration: User to input duration of the leakage test. Default value is 1 minute.


In progress: The pump will start running to reach the test point. (Graph below: Phase 1 and 2)

The leakage measurement will start once the pressure has reached the interested test point. Meanwhile, the timer will start counting down from the pre-set duration. (Graph below: Phase 3)

Temperature (VC 10 pro)

Temperature sensor

In the "Temperature" menu the user can specify that the temperature sensor is displayed on the display/working screen. A tick (✓) shows that the option is activated.


The precondition for this is that a temperature sensor is connected to the vacuum controller. If no temperature sensor is connected, or an error is present or the temperature rises above 350 °C, three dashes are shown as the temperature value.

NOTICE Please comply with the temperature measurement range of the external temperature sensor specified in the "Technical data" section. The temperature display can be in either °C or °F (see the "Settings" menu).


Calibration

In order to calibrate the external temperature sensor, it must first be connected to the vacuum controller. The temperature sensor is then immersed in a vessel with water, for example. The water temperature is determined using a second, already calibrated temperature measuring device. The measured temperature must then be entered in the menu under "Temperature calibration". Once the entry has been confirmed with "OK", the external temperature sensor is calibrated.

Timer

Timer functions

In the "Timer" menu the user can specify that the timer is displayed on the display/working screen. A tick (✓) shows that the

option is activated. This setting allows the user to specify the actual time for the evacuation procedure.

A default time can also be set for the timer. This setting allows the user to start evacuation for a standard time. The device can emit a signal (beep) after the set target time has elapsed.

NOTICE The user can stop the evacuation procedure before expiry of the set time. In this case the countdown of the timer is interrupted.

Operating mode

Operating mode A:

In this operating mode, the set target value is not saved when the current run comes to an end or the device is switched off.

Operating mode B:

In this operating mode, the set target value is saved when the

current run comes to an end or the device is switched off, and the value can be changed.

Operating mode C:

In this operating mode, the set target value is saved when the current run comes to an end or the device is switched off, and the value cannot be changed.

Operating mode D:

In this operating mode, the target value is applied if it does not change after 3 seconds, or when the yellow background is displayed behind the target value. In operating modes A and B, the target value is applied immediately when the target value is changed by turning the knob.

In operating mode C, the target value cannot be changed.

Display


In the "Display" menu the user can specify what information will be displayed on the main screen.

Graph

In the "Graph" menu, user can observe the real-time pressure change on the shown graph.

Programs

Select:

To select a program, use the rotating/pressing knob and press "Select". Only programs that have at least one program segment can be selected.

When a program has been successfully selected, this is indicated by a tick mark (✓). The vacuum controller is then in program mode. In the main screen, the program view can then be selected by turning the knob to the right. Turning the knob to the left returns the user to the main screen.

Start:

To start a selected program, the push knob must be pressed in the main screen.

Edit:

Edit the selected program parameters. Start to edit the selected program parameters by pressing on menu option "Edit" with rotating/pressing knob. The user can edit, delete or insert one selected program segment in the program.

Delete:

Deletes the selected program. If a selected program is deleted by pressing on menu option "Delete" with rotating/pressing knob, all the program parameters will be emptied. The tick (✓) disappears.

Details for editing the program

| No. | Pressure | hh:mm:ss |
|--------------------|----------|----------|
| 01 | 600 mbar | 00:00:10 |
| 02 | 550 mbar | 00:00:20 |
| 03 | 500 mbar | 00:00:30 |
| 04 | 450 mbar | 00:00:40 |
| 05 | 400 mbar | 00:00:50 |
| 06 | 350 mbar | 00:01:00 |
| 07 | 300 mbar | 00:01:10 |
| 08 | 250 mbar | 00:01:20 |
| edit insert delete | | |

In this program, user can define up to 10 segments. The selected segment is highlighted. Then, the user can edit, delete or insert a segment in this program. The program is save automatically.

Edit:


If the background of a selected value is yellow, the user can change the setting of the pressure value or the time value.

Delete:


When a highlighted segment is deleted, the subsequent segments move upward, closing the gap created by the deleted segment.

When a program is being edited, following screen appears.


Sample for editing the program

Example of saving the last measurement


Safety


In the "Password" menu, the user can protect the vacuum controller settings using a password. (factory setting: 000)

Service

In the "Service" menu, the valves, the pump or the refrigerated circulator can be operated individually and also checked to ensure that they are operating correctly. If no VACSTAR is connected, the menu item "Pump" is greyed out and not active.

Settings

Sound:

The "Sound" option allows the user to activate/deactivate the key-press sound and to set the volume.

Factory settings:

Select the "Factory settings" option by turning and pressing the rotating/pressing knob. The system requests confirmation to re-create the factory settings. You can choose if you want to reset only the "Menu values", or only the "Programs". Alternatively, you can reset all by selecting "All". Pressing the "OK" key resets the system settings to the original standard values set at dispatch from the factory (see "Menu structure" illustration).

Communication:

The "Device name" option allows the user to input a device name. The device name appears on the opening screen when the device is switched on.

This can be useful when using multiple devices with different settings. This also facilitates the identification of each device when communicating with a PC via USB, RS232 and Bluetooth, for example.

Information:

The "Information" option offers the user an overview of the most important system settings of the vacuum controller VC 10 lite/pro.

Languages:

The "Languages" option allows the user to select the desired language by turning and pressing the rotating/pressing knob (D). A tick (✓) indicates the language that is set for the system.

Units:

The "Units" option allows the user to select the desired unit for displaying the temperature and pressure by turning and pressing the rotating/pressing knob. For the temperature, the choice is between "°C" or "°F". For the pressure, the choice is "mbar", "hPa", "mmHg" or "Torr". A tick (✓) indicates the unit that is set for the system.

Display:

The "Display" option allows the user to change the background color and brightness of the working screen.

Interfaces and outputs

The device can be operated in "Remote" mode via an RS 232 or USB interface using the laboratory software labworldsoft®.

The RS 232 interface at the back of the device is fitted with a 9-pole SUB-D port which can be connected to a PC. The pins are assigned serial signals.

NOTICE

Please comply with the system requirements together with the operating instructions and help section included with the software.

USB interface

The Universal Serial Bus (USB) is a serial bus for connecting the device to the PC. Equipped with USB devices can be connected to a PC during operation (hot plugging). Connected devices and their properties are automatically recognized.

Use the USB interface in conjunction with labworldsoft® for operation in "Remote" mode and also to update the firmware.

Installation

First, download the latest driver for devices with USB in-terface from <http://www.com/lws/download/usb-driver.zip> and install the driver by running the setup file. Then connect the device through the USB data cable to the PC.

Serial interface RS 232 (V 24)

Configuration:

- The functions of the interface connections between the device and the automation system are chosen from the signals specified in EIA standard RS 232 in accordance with DIN 66 020 Part 1.
- For the electrical characteristics of the interface and the allocation of signal status, standard RS 232 applies in accordance with DIN 66 259 Part 1.
- Transmission procedure: asynchronous character transmission in start-stop mode.
- Type of transmission: full duplex.
- Character format: character representation in accordance with data format in DIN 66 022 for start-stop mode. 1 start bit; 7 character bits; 1 parity bit (even); 1 stop bit.
- Transmission speed: 9600 bit/s
- Data flow control: none
- Access procedure: data transfer from the device to the computer takes place only at the computer's request

Command syntax and format

The following applies to the command set:

- Commands are generally sent from the computer (Master) to the stirrer machine (Slave).
- The stirrer machine sends only at the computer's request. Even fault indications cannot be sent spontaneously from the device to the computer (automation system).
- Commands are transmitted in capital letters.
- Commands and parameters including successive parameters are separated by at least one space (Code: hex 0x20).
- Each individual command (incl. parameters and data) and each response are terminated with Blank CR LF (Code: hex 0x20 hex 0x0d hex 0x20 hex 0x0A) and have a maximum length of 50 characters.
- The decimal separator in a number is a dot (Code: hex 0x2E).

The above details correspond as far as possible to the recommendations of the NAMUR working party (NAMUR recommendations for the design of electrical plug connections for analogue and digital signal transmission on individual items of laboratory control equipment, rev. 1.1).

The NAMUR commands and the additional specific commands serve only as low level commands for communication between the stirrer machine and the PC. With a suitable terminal or communications programme these commands can be transmitted directly to the stirrer equipment.

The software package, labworldsoft®, provides a convenient tool for controlling the equipment and collecting data under MS Windows, and includes graphical entry features, for motor speed ramps for example.

The following table summarises the (NAMUR) commands understood by the control equipment.

| NAMUR Commands | Function |
|----------------|---------------------------------|
| IN_NAME | Read device name |
| IN_PV_3 | Read PT 1000 value |
| IN_PV_66 | Read current pressure value |
| IN_SP_66 | Read pressure target value |
| IN_SP_70 | Read hysteresis target value |
| IN_MODE_66 | Read current evacuation mode |
| IN_ERROR | Read error status |
| OUT_SP_66 | Set pressure target value |
| OUT_SP_70 | Set hysteresis value |
| OUT_MODE_66 | Set evacuation mode |
| START_66 | Start evacuation |
| STOP_66 | Stop evacuation |
| RESET | Switch to normal operating mode |

PC 1.1 Cable (Device to PC)

Required for connecting the 9-pin socket to a PC.


Fig. 10

Connection VC 10 lite/pro - PC





Fig. 11

Connection Refrigerated circulator

PC1.3 Cable (Device - Refrigerated circulator)

Connection (Device - Refrigerated circulator)

Connection cable

This cable is required to connect the VC 10 lite/pro to the VACSTAR (speed mode).

Fig. 12

Maintenance and cleaning

The device is maintenance-free. It is only subject to the natural wear and tear of components and their statistical failure rate.

Cleaning

- For cleaning disconnect the mains plug!
- Use only cleaning agents which have been approved by to clean devices.
 - These are water (with tenside) and isopropanol.
- Wear protective gloves during cleaning the devices.
- Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.
- Do not allow moisture to get into the device when cleaning.
- Before using another than the recommended method for cleaning or decontamination, the user must ascertain with that this method does not destroy the device.

Spare parts order

When ordering spare parts, please give:

- device type.
- serial number, see type plate.
- position number and description of spare part, see www.ika.com.
- software version

Repair

Please only send devices in for repair that have been cleaned and are free of materials which might present health hazards.

For repair, please request the "Safety Declaration (Decontamination Certificate)" from or use the download printout of it from website at www.ika.com.

If your appliance requires repair, return it in its original packaging. Storage packaging is not sufficient when sending the device - also use appropriate transport packaging.

Error codes

Error:

The fault is shown by an error message in the display as following if the error occurs, e.g. Error 4.
Proceed as follows in such cases:

- Switch the device switch off,
 - carry out corrective measures,
 - switch the device switch on again.
- ⇒ Measurement is stopped.

Warning:

- Measurement is stopped.
- Warning is displayed.
- Warning can be confirmed.

Information:

- Measurement continues.
- Information is displayed.
- Information can be confirmed.

| Description | Detection | Reason | Action |
|--|---|--|--|
| Information No Pressure Change | Pressure gradient analysis does not show any deviation. | Pump is not running. Vacuum hose isn't connected. Failure in hose system. Recipient isn't closed. | Check pump control cable connection. Check pump power supply. Switch on pump. the position. Check correct vacuum hose connection (IN/ OUT/VENT). Close recipient. |
| Information System Not Tight | Pressure gradient analysis shows pressure deviation, but set value cannot be reached. | Vacuum hose connection isn't tight. Recipient isn't tight. Vacuum pump power isn't sufficient. Venting valve leakage. Boiling point of solvent is reached before. | Check vacuum hose connection. Check recipient. Check technical data of the pump. Call service department. Check set value. |
| Information Venting Error | Pressure gradient analysis does not show any deviation after pressing the "Venting" button. | Recipient wasn't evacuated at all. Failure in hose system. Venting valve leakage. Front foil button is damaged. | Evacuated recipient. Check correct vacuum hose connection (IN/ OUT/VENT). Call service department. |
| Pump not connected | Connection to Vacstar(speed mode) not present. | Damaged or no cable connected | Check connection cable is present and well connected. |
| Warning Boiling Point Error | Boiling point cannot be detected. Pressure gradient analysis does not show boiling point. | Heating bath does not heat up. Solvent with extreme low boiling point is used. | Check heating bath. Solvent distillation only manually. |
| Warning External Temperature Is Too High | Measured temperature is too high. External temperature is higher than maximum value | Medium temperature is too high. External sensor is broken. | Check the medium temperature. Call service department. |
| Warning External Temperature Is Too Low | Measured temperature is too low. External temperature is lower than minimum value. | Medium temperature is too low. External sensor is short-connected. | Check the medium temperature. Call service department. |
| Warning Control Pump State Error | Pump is disconnected during analog speed control mode. Pump is connected during 2-point control mode. | Pump is disconnected during analog speed control mode. Pump is connected during 2-point control mode | Connect the pump. Disconnect the analog pump. |
| Warning RS232 Communication Error USB Communication Error | Watchdog time elapsed. | labworldsoft® uses wrong device for VC 10 lite/pro. Wrong tools are used to communicate with VC 10 lite/pro. Unstable connection. Watchdog was set in device or via PC command. | Check the labworldsoft® setting. Check the PC tools. Check the connection. Disable watchdog when not necessary. Increase sample rate. |

| Description | Detection | Reason | Action |
|--|---|---|---|
| Warning Temperature Calibration Error | Temperature calibration value is out of range. | Wrong simulator resistance is chosen. Set wrong temperature calibration value. | Check the simulator resistance. Call service department. |
| Warning Pressure Calibration Error | Pressure calibration value is out of range. | Input wrong vacuum actual value when calibration. Set wrong vacuum calibration value. Sensor is broken. | Calibrate again. Check the set calibration value. Call service department. |
| Error 3 Device Temperature Error. | Analysis of internal PCB temp. Sensor. Inside temperature reaches limit value. | Room temperature > 40 °C. PCB or valve is broken. | Switch off device for cool down. Call service department. |
| Error 9 Logic Storage Error. Display Storage Error. | Read or write internal memory error. | Read or write internal memory error. | Call service department. |
| Error 68 Pressure Out Of Range. | Pressure is out of range. Pressure is bigger than room pressure. | Pressure in the recipient is too high. Sensor is broken. | Check the air flow of pump and vent recipient. Call service department. |
| Error 69 Pressure Sensor Error. | Pressure sensor output is too low. | Sensor is not connected. Sensor is broken. | Check the sensor connection. Call service department. |
| Error 71 Internal Communication Error | Internal Watchdog elapsed. | Communication lost between logic board and display board. | Restart device (unplug mains cable, plug in mains cable). Check communication between logic board and display board. Call service department. |

If the actions described fails to resolve the fault or another error code is displayed then take one of the following steps:

- contact the service department,
- send the instrument for repair, including a short description of the fault.

Accessories

Accessories see www.ika.com.

Product contact parts

| Designation | Material |
|--------------------|-----------------|
| Connection spigots | PP |
| Distributor | PPS |
| Pressure sensor | FPM/AL2O3 |
| Vacuum valve | PEEK; EPDM; FKM |
| Venting valve | PEEK; EPDM; FKM |

Technical data

| | Unit | VC 10 lite | VC 10 pro |
|--|------|---|-----------|
| Connection diameter suction side | mm | 8 | |
| Connection diameter pressure side | mm | 8 | |
| Connection diameter venting side | mm | 8 | |
| Input pressure min. | mbar | 1 | |
| Input pressure max. | mbar | 1050 | |
| Boiling point detection | | - | yes |
| Solvent library | | yes | |
| Two-point control | | - | yes |
| Analog speed vacuum control | | yes | |
| Display | | TFT | |
| Pressure unit / scale | | mbar, hPa, mmHg, Torr | |
| Vacuum sensor | | yes | |
| Vacuum sensor type | | ceramic Al ₂ O ₃ | |
| Pressure max. for pressure sensor | bar | 1.6 | |
| Measurement range (absolute) min. | mbar | 1 | |
| Measurement range (absolute) max. | mbar | 1100 | |
| Control range min. | mbar | 1 | |
| Control range max. | mbar | 1100 | |
| Resolution pressure | mbar | 1 | |
| Measurement uncertainty | mbar | 1 | |
| Medium temperature (gas) min. | °C | 5 | |
| Medium temperature (gas) max. | °C | 40 | |
| Vacuum valve | | - | yes |
| Venting valve | | yes | |
| Connection for ext. temperature sensor | | - | PT 1000 |
| Temperature unit | | - | °C/°F |
| Temperature measuring range min. | °C | - | -10 |
| Temperature measuring range max. | °C | - | 200 |
| Temperature measurement resolution | K | - | 1 |
| Accuracy of temperature measurement | K | - | ±1 |
| Timer | | yes | |
| Time setting min. | s | 1 | |
| Time setting max. | min | 6000 | |
| Vacuum speed control interface | | VACSTAR | |
| Material in contact with medium | | Al ₂ O ₃ , PTFE, FPM, PPS | |
| Housing material | | PBT | |
| Fastening | | stand/clamp | |
| Fastening diameter | mm | 16 | |

| | Unit | VC 10 lite | VC 10 pro |
|--|------|----------------|-----------|
| Dimensions (B x H x T) | mm | 95 x 150 x 110 | |
| Weight | kg | 1.5 | |
| Permissible ambient temperature min. | °C | 5 | |
| Permissible ambient temperature max. | °C | 40 | |
| Permissible relative humidity | % | 80 | |
| Protection class according to DIN EN 60529 | | IP 20 | |
| RS 232 interface | | yes | |
| USB interface | | yes | |
| Voltage | V | 100 - 240 | |
| Frequency | Hz | 50/60 | |
| Power input | W | 24 | |
| Power input standby | W | 2 | |
| DC Voltage | V | 24 | |
| Current consumption | mA | 1000 | |

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Phone: +90 216 394 43 43

eMail: sales.turkey@ika.com

Discover and order the fascinating products of IKA online:
www.ika.com

IKAworldwide

IKAworldwide /// #lookattheblue

@IKAworldwide

IKA

designed for scientists

IKA Vacstar digital

Operating instructions
Source language: German

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Fig. 1

Pos. Description (EN)

| | |
|---|------------------------|
| A | Handle |
| B | Display |
| C | Rotating/pressing knob |
| D | Main switch (on/off) |

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Fig. 1

Pos. Description (EN)

| | |
|---|--|
| E | Handle securing screw |
| F | Vacuum controller VC 10 stand hole |
| G | Stand securing screw |
| H | Fan / ventilation slit |
| I | Hose connection for suction line d= 8 mm (INLET) |
| J | Accessory fastening (emission condenser) |
| K | Connection for vacuum controller VC 10 (Mini DIN) |
| L | USB interface |
| M | RS 232 interface |
| N | Water valve connection |
| O | Hose connection for pressure line d= 8 mm (OUTLET) |
| P | Power supply cable connection |

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| | Page |
|--------------------------------------|------|
| Device setup | 2/3 |
| Declaration of conformity | 5 |
| Explication of warning symbols | 5 |
| Warranty | 5 |
| Safety instructions..... | 6 |
| Unpacking | 7 |
| Correct use | 8 |
| Useful information | 8 |
| Setting up..... | 9 |
| Commissioning | 12 |
| Interfaces and outputs | 13 |
| Maintenance and cleaning..... | 14 |
| Error codes | 15 |
| Accessories | 15 |
| Product contact parts..... | 15 |
| Technical data..... | 16 |

Declaration of conformity

We declare under our sole responsibility that the product to which this declaration relates is in conformity with directives 2014/35/EU, 2006/42/EC, 2014/30/EU and 2011/65/EU and conforms with the following standards or normative documents: EN 61010-1, EN 60529, EN 61326-1 and EN ISO 12100.

A copy of the complete EU Declaration of Conformity can be requested at sales@ika.com.

Explication of warning symbols

| | |
|--|--|
|  DANGER | Indicates an (extremely) hazardous situation, which, if not avoided, will result in death, serious injury. |
|  WARNING | Indicates a potentially hazardous situation, which, if not avoided, can result in death, serious injury. |
|  CAUTION | Indicates a potentially hazardous situation, which, if not avoided, can result in injury. |
|  NOTE | Indicates practices which, if not avoided, can result in equipment damage. |

Warranty

According to **IKA**'s Terms and Conditions of sale and delivery, this product is covered by a warranty for a period of 24 months. Please contact your dealer for any warranty claims. If you wish, you can send the device directly to our factory. Please provide the sales invoice and state the reasons for your guarantee claim. In this case, you are responsible for shipping and handling costs.

The warranty does not cover wearing parts, nor defects that are the result of improper use, insufficient care and maintenance or failure to follow the instructions in this operating manual.

Safety instructions

For your protection

NOTE Please read the instruction manual in full before use and follow the safety instructions.

- Keep the instruction manual in a place where it can be accessed easily.
- Ensure that only trained staff use the device.
- Be sure to comply with all safety instructions, directives and all matters of health, safety and accident prevention in the workplace.
- The device and all parts of the device must not be used on people or animals.

DANGER Always wear personal protective equipment in accordance with the hazard class of the media being worked with - otherwise there are dangers due to spraying liquids, and release of toxic or flammable gases.

- Do not expose human or animal body parts to vacuum.
- Do not work with the device underwater or underground.

Device configuration

DANGER The vacuum pump IKA Vacstar digital is not designed to be set up in potentially explosive atmospheres.

- Set up the device in accordance with the chapter "Setting up" and connect the connection lines and interfaces as described.
- Set up the device on a stable, even, non-flammable surface.

- Never work with a faulty or incorrectly connected device.
- Set up the vacuum pumps in appropriate and functioning extractor hoods, or set up an exhaust gas line to them. Ensure that the exhaust gas line cannot be kinked! The maximum permitted exhaust gas line length is 2 metres.
- Occurrence of explosive mixtures; if appropriate connect inert gas for ventilation and/or thinning.

Accessories

- Safe operation can only be ensured when working with accessories as described in the "Accessories" section.
- Carefully observe the operating instructions for additional devices (e.g. rotary evaporators, vacuum controllers), with which the vacuum pump IKA Vacstar digital is operated.
- The pressure at the gas inlet and outlet must not exceed 1100 mbar.

- Elastic elements may be pressed together under vacuum.
- Only use flexible hose lines.
- Observe your emergency measures for power failure and ensure that the device is put in a safe state (see the chapter: Commissioning, Operating mode).

Work with the device

DANGER You can prevent the release of solvent vapours into the atmosphere using a downstream emission condenser.

The vacuum pump is not suitable for use with self-igniting materials, materials that are flammable without an air supply, or explosive materials.

Do not operate the pump while it is open. Otherwise there is a danger of severe injuries as a hand may be pulled into the moving parts.

WARNING Inhaling or coming into contact with media such as poisonous liquids, gases, spray mists, vapours, dusts or biological materials can endanger the health of the user. Ensure that all connections are well sealed and free of leaks if you are working with such media.

- The vacuum pump IKA Vacstar digital must only be operated under the conditions described in the chapter "Technical data".
- Prevent release of the materials listed above. Take measures to protect staff and the environment.
- Pay attention to possible interactions or chemical or physical

reactions when working with media at reduced pressure and increased temperature.

- There can be electrostatic processes between the medium and the device which can lead to direct danger.
- Some medium may be released due to the residual leakage rate of the device.
- Before commissioning check that all the housing parts are present and fastened to the device.
- Do not lift the pump if the handle is loose and the handle securing screw (E) is loose or missing!
- Only switch the pump on if the pump is standing vertically.
- Connect the hose connections (INLET-OUTLET) and interfaces in accordance with the labelling on the device and the operating instructions.
- Please note that the temperature of the evacuated medium is below its ignition temperature. The pumping process (compression) increases the temperature of the medium additionally.
- Ensure that vapours containing solvent can be sucked into the pump.
- Do not use the pump to create pressure.
- Observe the permitted pressure at the inlet and outlet side; see the chapter "Technical data".
- The gas flow must only be regulated/throttled in the suction-side line.
- Use a solenoid valve or a check valve in the suction line if there are several load units.
- When using an emission condenser ensure that the coolant is flowing freely.

Protection of the device

WARNING Ensure that the substances used with the device are compatible with the materials of the device parts that come into contact with the product; see the chapter "Product contact parts".

CAUTION The specified settings on the rating plate must coincide with the actual power supply.

- The device is only disconnected from the power supply when the power or device plug is removed.
- The socket for the power supply cable must be easy to access.
- Covers or parts that can be removed from the device without tools must later be refitted to ensure safe operation, unless anything else is connected at this point. This will prevent the infiltration of foreign objects, liquids and other contaminants.
- Unwanted liquids that have entered can be removed by extracting air during no-load operation.
- Avoid knocks or impact to the device.


- The device must only be opened by trained, skilled personnel (who have been authorised to do so).
- The device must not be modified without authorisation from IKA.
- To ensure sufficient cooling for the vacuum pump IKA Vacstar digital, the ventilation slits on the housing must not be covered.
- Only use original replacement parts for maintenance to ensure that the device works reliably.
- Watch out for water condensation inside and outside the device. Warm up the device first, if it has been brought in from a cold environment.
- Never fasten the vacuum pump over a heating bath.
- Ensure that no solids and/or liquids can enter the diaphragm of the vacuum pump IKA Vacstar digital through the suction line of the pump. This damages the diaphragm and other internal parts of the pump. This reduces delivery capacity and the final pressure can no longer be achieved. Deposits may be produced in the interior, reducing service life and causing leaks.


Unpacking

Unpacking

- Unpack the device carefully.
- Any damage should be notified immediately to the shipping agent (post office, railway network or transport company).

Scope of delivery

| | |
|---|---|
| Diaphragm vacuum pump IKA Vacstar digital | |
| Power supply cable Laptop Euro |  |
| Power supply cable Laptop GB | |
| Power supply cable Laptop CH | |
| USB A - USB B cable | |
| Vacuum hose 1m |  |

| | |
|------------------------|---|
| Silencer |  |
| Connecting cable |  |
| Operating instructions | |
| Warranty card | |

Correct use

Use

WARNING

The device must not be used for:
 - **evacuating biological habitats,**
 - **evacuating explosive, corrosive or similar gases,**
 - **throughput/use of liquids**

The device is not suitable for:
 - **pumping up objects**
 - **pumping liquids**
 - **use in potentially explosive atmospheres**
 - **underwater**
 - **underground**

Standalone operation

The vacuum pump can be operated without a controller. The delivery capacity can be set using the adjustable speed. Generally speaking, the smaller the volume to be evacuated, the lower the speed should be. The final pressure is not controlled.

Operation with accessories

Along with the accessories recommended by **IKA**, the vacuum pump **IKA** Vacstar digital is suitable for controlled evacuation to a specified final pressure e.g. for operating an **IKA** rotary evaporator. It can also be used for classical separation, filtration or drying in the laboratory.

With an appropriate vacuum controller (e.g. VC 10) or with a rotary evaporator (e.g. RV 10 auto), the Vacstar digital can be operated automatically, for automatic boiling point detection, scheduled pressure-time curves or programs from the solvent library, for example.

Intended use: Tabletop device

The laboratory diaphragm vacuum pump **IKA** Vacstar digital is used together with other laboratory devices to create a vacuum.

For specific vacuum control, an additional vacuum controller **IKA** VC 10 or **IKA** rotary evaporator RV 10 auto is required.

Area of use (only indoors)

The device is suitable for use in residential areas and all other areas.

The safety of the user cannot be ensured:

- if the device is used in conjunctions with accessories not made or recommended by the manufacturer!

- if modifications are made to the device or parts of the device by third parties!
- if the device is operated in contravention of the safety instructions!

Useful information

DANGER

You can prevent the release of solvent vapours into the atmosphere using a downstream emission condenser.

CAUTION

Special care must be taken when working with gas mixtures that include condensable gases (e.g. solvents). If diaphragm pumps are used with such gases they condense in the pump's air flow. This leads to a build-up of pressure in the vacuum chamber, which in turn influences the effectiveness and service life of the diaphragms and valves.

If liquid or condensate enters the vacuum chambers, the specified minimum working pressure is no longer achieved.


To protect the inner valves and diaphragms against liquid, use an upstream condensate separator, such as a Woulff bottle.

CAUTION

It is generally recommended that you run the pump with no load at the end of the process (no-load operation, approx. 3 up to max. 5 minutes) to dispel any residues in the pump.

If applicable, vent the recipient. To do this, remove the suction line during manual operation.

NOTE

Observe the operating instructions of the relevant devices.

Diaphragm vacuum pumps are oscillating displacement pumps. The diaphragms expand the volume of the vacuum chamber to produce an underpressure, which sucks air into the chamber. The air is then pressed out of the vacuum chamber by subsequent compression. Changes between the inlet and outlet are controlled by diaphragm valves.

The physical unit for pressure is 1 Pascal [Pa].

However people often still use the unit millibar [mbar] colloquially.

1 mbar = 100 Pa
 1 bar = 10⁵ Pa
 1 Pa = 10⁻⁵ bar
 1 Pa = 0.01 mbar

A distinction is made between different ranges in vacuum technology.

| | | |
|-------------------|---|---|
| Rough vacuum: | 10 ⁵ to 10 ² Pa | (1000 to 1 mbar) |
| Fine vacuum: | 10 ² to 10 ⁻¹ Pa | (1 to 10 ⁻³ mbar) |
| High vacuum: | 10 ⁻¹ to 10 ⁻⁵ Pa | (10 ⁻³ to 10 ⁻⁷ mbar) |
| Ultra high vacuum | < 10 ⁻⁵ Pa | (< 10 ⁻⁷ mbar) |

The periodic movements of the diaphragms create a vacuum at the pump inlet. Diaphragm vacuum pumps work with a rough vacuum of up to 10² Pa. Lower pressures cannot be achieved with these pumps.

To create fine, high and ultra high vacuums, oil-sealed vacuum pumps, screw pumps, diffusion pumps or turbomolecular pumps must be used, for example. Diaphragm vacuum pumps are often used as pumps for the preliminary stage in these contexts.

Compressing and releasing the pressure on air creates a noise. The pump is quieter the less air is transported. This means that there is less operating noise during normal operation if the pressure is reduced.

A fitted exhaust hose with a sound absorber can also be used to reduce noise.


In automatic operation with controlled speed the pump reduces its speed when it reaches the target value. If the leakage rate in the entire system is low enough, the speed may even be set to "zero".

When combined with the vacuum controller VC 10 or the rotary evaporator with integrated vacuum controller RV 10 auto, the pump can also be used for "two-position control" (valve-controlled) or "analogue speed control" of the vacuum.

Two-position control

When the target value is reached, the suction line is interrupted. Due to the small delay between the times when the pressure value is detected, when comparison is made with the target value, when a vacuum valve is switched off and when the suction line of the vacuum-producing system are switched off, it is easy to go below the set target value. If the pressure in the system goes up again due to a natural leakage rate, the suction line valve opens.


Schematic view of two-position vacuum control

Analogue speed control

With this control, the speed of the pump, and therefore its suction power, is reduced the closer the measured pressure reaches the target value. Once the target value is reached, the pump only operates according to the leakage rate of the system. This enables quiet functioning, precise vacuum control and a long service life for the moving pump parts. Exact automatic boiling point recognition is possible with this type of control, meaning that the system reaches and keeps the solvent at its boiling point in automatic mode.

Schematic view of speed vacuum control

Setting up

NOTE

Ensure that you follow the general rule, that the recipient (load/vacuum container/glass condenser) should always be connected to the suction line at its highest point. This reduces the risk of liquids entering the vacuum pump.

Liquid in the pump chambers makes the pump properties worse.

Connect the hose connections (INLET/OUTLET) and interfaces in accordance with the labelling on the pump and the operating instructions.

The connections are barbed hose connectors. Ensure that the hoses have hose clips if applicable. Lay the hoses in such a way that no condensate can flow into the pump. Ensure that the air can enter the ventilation slit (H) unhindered.

Install a condensate separator (e.g. a Woulff bottle) in the suction line before the intake spigot of the pump to protect against water coming in.

If solvent vapours are sucked in, an emission condenser (accessory) after the pump helps to condense it, and prevents it being released into the atmosphere.

Connection of interfaces

NOTE

Observe the relevant connections (see Fig. 1).

O: Hose connection d= 8 mm OUTLET

Link this connection to the emission condenser with a vacuum hose, or fit a sound absorber to the end.

NOTE

Put the end of the hose in the extractor hood!

Check there is a free outlet on the pressurised side!

Do not use a throttle on the pressurised side and do not close the outlet! Connect the exhaust line to this connection.

K: Connection for vacuum controller VC 10 (Mini DIN)

You can connect the vacuum pump and the vacuum controller VC 10 or the rotary evaporator RV10 auto with the analogue connection cable for precise speed-controlled vacuum control. The vacuum controller detects the pump and switches to speed-vacuum control mode. Two-position control is deactivated.

L: USB interface

Connect the vacuum pump IKA Vacstar digital to a PC with a USB A - USB B cable. Any device software updates can be loaded using the IKA FUT software tool.

M: RS 232 interface

You can connect the IKA Vacstar digital vacuum pump to a PC using an RS 232 interface cable. The pump can be operated in conjunction with other devices with labworldsoft® laboratory device software. For further information see the chapter "Interfaces and outputs".

N: Water valve connection

Connect the optional water throttle valve RV 10.5001 to the diaphragm vacuum pump. The water flow to the emission condenser is controlled using the water valve. The valve is opened as soon as the pump is switched on.

I: Hose connection for suction line d= 8 mm INLET

Connect the suction line to this connection. Link this connection to the recipient (rotary evaporator condenser, laboratory reactor etc.) with a vacuum hose.

P: Power supply cable connection

Check that the voltage information on the rating plate matches your mains supply. Connect the power supply cable to supply electricity.

IKA Vacstar digital with vacuum controller VC 10

NOTE


Cannot be operated with optional water valve RV 10.5001 in conjunction with VC 2.4.

Refer to the operating instructions for the vacuum controller VC 10!

When using the optional pump cut-out unit VC 2.4 (ref. 6) the pump is switched off when the target value is reached (operating mode B).

This device configuration is also possible without the pump cut-out unit VC 2.4. The pump runs on. For good control, the pump should be operated at a low speed.


When the target value is reached, the suction line is closed by the vacuum valve integrated in the VC 10. The pump speed can be changed by manually setting the speed.

- 1 Recipient (load, e.g. rotary evaporator, reactor)
- 2 Vacuum separator (e.g. Woulff bottle)
- 3 Vent valve
- 4 Pressure sensor
- 5 Vacuum valve/ball valve
- 6 Pump cut-out unit VC 2.4
- 7 Emission condenser

Speed regulation

IKA Vacstar digital with vacuum controller VC 10 or rotary evaporator RV 10 auto


NOTE

Observe the operating instructions for the devices.

Automatic setting of the operating mode "Speed regulation" as soon as the analogue connecting cable (ref. 8) connects the pump with the vacuum controller or rotary evaporator.

When the target value is reached, the pump speed is automatically reduced to a minimum of "0 rpm". Once the target value is reached, the pump applies suction according to the leakage rate of the vacuum system. Automatic boiling point recognition is possible.

- 1 Recipient (load, e.g. rotary evaporator, reactor)
- 2 Vacuum separator (e.g. Woulff bottle)
- 7 Emission condenser
- 8 Analogue connecting cable

Commissioning

NOTE Observe the operating instructions for the devices.

Check if the voltage specified on the type plate matches the available mains voltage. If these conditions are met, the device is ready for operation when it is plugged in to the mains.


The socket used must be earthed (grounded contact).

Otherwise safe operation is not guaranteed or the device may be damaged. Observe the ambient conditions indicated under "Technical data".

Switching on

1. Switch on the device using the main switch (D). All display segments are shown to enable a visual check.

2. Software version number display (e.g. vers. 1.00).

3. Operating mode display (A; B). (Factory setting: Operating mode A)

The operating mode can be changed when you switch on the device.

Operating mode A:

In this operating mode, the set target value is not saved when the current run comes to an end or the device is switched off. No automatic restart after mains power is cut off.

Operating mode B:

In this operating mode, the set target value is saved when the current run comes to an end or the device is switched off, and the value can be changed.

Operating mode B is required particularly in combination with the pump cut-out unit VC 2.4.

Automatic restart after mains power is cut off.

Changing operating mode:


1. Press the rotating/pressing knob (C) and switch on the device at the main switch (D) at the same time to change between the two operating modes.


In standalone operation or two-position operation in combination with the vacuum controller VC 10 the target speed can be set by turning the rotating/pressing knob (C). Push to start the pump and push again to stop it. A flashing light on the display (B) indicates when the device is running.

Switching on with several interfaces

If the pump is operated via the analogue connecting cable with the vacuum controller VC 10 or the rotary evaporator RV 10 auto, "VC 10" is shown on the display.


If the vacuum controller is connected to a PC, "PC" is shown on the display (B).

Switching off

Switch off the device using the main switch (D).

Unplug the appliance from the mains socket to disconnect it from the mains.

Interfaces and outputs

NOTE Please comply with the system requirements together with the operating instructions and help section included with the software.

The device can be operated in "Remote" mode via an RS 232 or USB interface using the laboratory software labworldsoft®. The RS 232 interface (M) at the back of the device is fitted with a 9-pole SUB-D jack which can be connected to a PC. The pins are assigned serial signals.

USB Interface

The Universal Serial Bus (USB) is a serial bus system which allows the device to be connected to the PC. Devices that support USB can be connected to each other whilst they are running (hot plugging) and provide automatic recognition of the connected devices and their properties.

Use the USB interface in conjunction with labworldsoft® for operation in "Remote" mode and for updating the firmware using the "Firmware update tool".

Installation

Before the device is connected with the PC using the USB data cable, the USB driver must be installed.

The USB driver can be downloaded from the website:

Serial interface RS 232 (V 24)

Configuration:

- The functions of the interface circuit between the device and the automation system are a selection from the signals specified in the EIA standard RS232 as per DIN 66020 Part 1.
- Standard RS 232, corresponding to DIN 66259 Part 1 is valid for the electric characteristics of the interface circuits and assignment of signal states.
- Transmission process: Asynchronous character transmission in start-stop operation.
- Transmission type: Full duplex.
- Character format: Character composition according to data format in DIN 66022 for start-stop operation. 1 start bit, 7 character bits, 1 parity bit (even), 1 stop bit.
- Transmission speed: 9600 Bits/s.
- Data flow control: none
- Access method: Data transmission from the device to the computer only occurs after a request from the computer.

Command syntax and format

The following points should be noted for the instruction set:

- Commands are generally sent from the computer (master) to the device (slave).
- The device only responds to requests from the computer. Even error messages are not sent spontaneously from the device to the computer (automation system).
- The commands are transmitted in capital letters.
- Commands and parameters, as well as consecutive parameters, must be separated by at least one space (code: hex 0x20).
- Each individual command (including parameters and data) and all responses are completed with CRLF (code: hex 0x20 hex 0x0d hex 0x0a) and can have a maximum length of 50 characters.
- The dot is used for decimal separators in a floating-point value (code: hex 0x2E).

The details given above generally comply with the recommendations of NAMUR (NAMUR recommendations for the design of electrical plug-in connectors for analogue and digital signals in laboratory MSR devices. Rev. 1.1).

The NAMUR commands and the additional IKA-specific commands are only used as low-level commands for communication between the device and the PC. With an appropriate terminal or communication program, these commands can be transmitted directly to the device.

NAMUR commands Function

| NAMUR commands | Function |
|----------------|---|
| IN_NAME | Read device name |
| IN_PV_4 | Read current speed value |
| IN_SP_4 | Read rated speed value |
| OUT_SP_4 | Adjust the rated speed value |
| START_4 | Start evacuation |
| STOP_4 | Stop evacuation |
| RESET | Switch to normal operating mode |
| OUT_NAME | Set device name |
| OUT_SP_42@n | Setting WD safety limit speed with set value echo |
| OUT_WD1@m | Watchdog mode 1 |
| OUT_WD2@m | Watchdog mode 2 |

"Watchdog" functions; monitoring of the serial data flow

If, once this function has been activated (see NAMUR commands), there is no retransmission of the command from the computer within the set time ("watchdog time"), the Evacuation function is switched off in accordance with the set "watchdog" mode or is returned to previously set target values. The data transmission may be interrupted by, for example, a crash in the operating system, a power failure in the PC or an issue with the connection cable between the computer and the device.

"Watchdog" mode 1

If event WD1 should occur, the evacuation function is switched off and ER 2 is displayed. Set watchdog time to m (20 - 1,500) seconds, with watchdog time echo. This command launches the watchdog function and must be transmitted within the set watchdog time.

"Watchdog" mode 2

If there is an interruption in data communications (longer than the set watchdog time), the speed target value is changed to the set WD safety speed limit. The warning PC 2 is displayed. The WD2 event can be reset with OUT_WD2@0 - this also stops the watchdog function.

Set watchdog time to m (20 - 1,500) seconds, with watchdog time echo. This command launches the watchdog function and must be transmitted within the set watchdog time.

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
Fig. 10

USB A - USB B cable

Required for connecting the USB interface (L) to a PC.


Fig. 11

Connecting cable

Required for connecting VC 10 and RV 10 auto with the diaphragm vacuum pump IKA Vacstar digital (speed operation).


Fig. 12

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Cleaning

For cleaning disconnect the mains plug!

Use only cleaning agents which have been approved by IKA to clean IKA devices.

These are water (with tenside) and isopropanol.

Wear protective gloves during cleaning the devices.

Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.

Do not allow moisture to get into the device when cleaning.

Before using another than the recommended method for cleaning or decontamination, the user must ascertain with IKA that this method does not destroy the device.

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- device type.
- serial number, see type plate.
- position number and description of spare part, see www.ika.com.

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For repair, please request the "Safety Declaration (Decontamination Certificate)" from IKA or use the download printout of it from IKA website at www.ika.com.

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If an error occurs, this is shown in an error on the display (B), e.g. Error 4.

Proceed as follows:

- Switch the device off on the device switch,
- Take corrective measure,
- Switch device back on.

| Error | Cause | Reason | Solution |
|-------|----------------------------|-----------------------|---|
| Er 3 | Overheat | The pump does not run | Switch off the pump, disconnect it from the mains, let it cool down. |
| Er 4 | Motor blocked/overload | The pump does not run | Switch off the pump, disconnect it from the mains Let the pump cool down. Ventilate the vacuum system to normal pressure. |
| Er 8 | Hall sensor defective | The pump does not run | Service |
| Er 9 | BLP self-test/EEPROM error | The pump does not run | Service |
| Er 10 | Outlet Choked | The pump does not run | User needs to check the outlet |
| Er 11 | Membrane stuck | The pump does not run | Air the pump and run at low speed for self-cleaning |
| Er 12 | Connection problem | The pump does not run | Service: check the motor cable connection |

In the event that the fault cannot be eliminated using the measures described or another error code is displayed:

- Please contact the Service Department,
- Send the device and a brief description of the fault.

Accessories

- Vacuum controller IKA VC 10
- Vacuum control valve IKA VCV 1 and VCV 2
- Check valve IKA VC 10.300
- Vacuum hose IKA VH.SI.8
- PC 1.1 cable (RS 232)
- Maintenance kit
- Vacuum safety emission condenser VSE 1
- Woulff bottle/water valve set VSW1

For further accessories see www.ika.de

Product contact parts

| Name | Material |
|-------------------|----------|
| Connection spigot | PPS |
| Connection piece | PPS |
| Distributor | PPS |
| Head piece | PPS |
| Diaphragm valve | FFPM |
| Diaphragm | NBR/PTFE |
| Clamp | PPS |

Technical data

| | Base unit | Value |
|--|------------------------|--|
| Delivery capacity max. | m³/h | 1.32 |
| Delivery capacity max. | L/min | 22 |
| Final pressure | mbar | 2 |
| Suction steps | | 4 |
| Cylinder | | 4 |
| Connection diameter suction side | mm | 8 |
| Connection diameter pressure side | mm | 8 |
| Gas ballast valve | | no |
| Inlet pressure min. | mbar | 2 |
| Inlet pressure max. | mbar | 1030 |
| Boiling point recognition | | no |
| Solvent library | | no |
| Two-position control | | yes |
| Analogue speed-vacuum control | | yes |
| Speed setting option | | rotating/pressing knob |
| Speed min. | rpm | 285 |
| Speed max. | rpm | 1200 |
| Display | | LED |
| Vacuum sensor | | no |
| Vent valve | | no |
| Condensate separator | | no |
| Emission condenser | | no |
| Volume at minimum pressure | dB(A) | 54 |
| Product contact material | | PTFE, FFPM, PPS, NBR |
| Housing material | | Cast alu, coated |
| Appropriate for clean room | | no |
| Explosion-proof | | no |
| Width | mm | 150 |
| Depth | mm | 370 |
| Height | mm | 375 |
| Weight | kg | 11.5 |
| Permissible ambient conditions | | 5 °C to 31 °C at 80 % relative humidity. 32 °C to 40 °C decreasing linearly to a maximum 50 % relative humidity. |
| Protection class according to DIN EN 60529 | | IP 20 |
| RS 232 interface | | yes |
| USB interface | | yes |
| Analogue output | | no |
| Voltage | V | 100-240 |
| Frequency | Hz | 50/60 |
| Power consumption | W | 130 |
| Standby power consumption | W | 1.5 |

| | Base unit | Value |
|----------------------------|--|-------|
| Protection class | | I |
| Overvoltage category | | II |
| Degree of contamination | | 2 |
| Use above max. sea level | m | 2000 |
| Firmware update | | yes |
| <i>Note</i> | | |
| Protective earth conductor | The protective earth conductor (PE) is not connected to the metal housing. (Internally encapsulated DC power supply) | |

*Average values.
Subject to technical changes!*

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IKA India Private Limited
IKA England LTD.
IKA Vietnam Company Limited