IKA

designed for scientists

VC 10 lite
VC 10 pro

Betriebsanleitung Ursprungssprache	DE	5
Operating instructions Mode d'emploi	EN	29
Руководство пользователя	FR	53
使用说明	RU	77
	ZH	101

Breve guía	ES	124
Brevi istruzioni	IT	125
Veiligheidsaanwijzingen	NL	127
Säkerhetsanvisningar	SV	128
Sikkerhedshenvisninger	DA	130
Sikkerhetsinformasjon	NO	131
Turvallisuusohjeet	FI	132
Instruções de segurança	PT	134
Wskazówki dotyczące bezpieczeństwa	PL	135
Bezpečnostní pokyny a informace	CS	137
Biztonsági utasítások	HU	138
Varnostna navodila	SL	139
Bezpečnostné pokyny	SK	141

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Contents

	Page
Device setup	2
Declaration of conformity	29
Explication of warning symbols	29
Warranty	29
Safety instructions	30
Unpacking	31
Intended use	31
Useful information	31
Setting up	33
Commissioning	36
Interfaces and outputs	46
Maintenance and cleaning	47
Error codes	48
Accessories	49
Product contact parts	49
Technical data	50

Declaration of conformity

We declare under our sole responsibility that this product corresponds to the regulations 2014/35/EG, 2014/30/EG and 2011/65/EG and conforms with the standards or standardized documents EN 61010-1, EN 60529, EN 61326-1 and EN ISO 12100.

A copy of the complete EU Declaration of Conformity or further declarations of conformity can be requested at sales@ika.com.

Explication of warning symbols

Indicates an (extremely) hazardous situation, which, if not avoided, will result in death, serious injury.

Indicates a potentially hazardous situation, which, if not avoided, can result in death, serious injury.

Indicates a potentially hazardous situation, which, if not avoided, can result in injury.

Indicates practices which, if not avoided, can result in equipment damage.

Warranty

In accordance with warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual.

Safety instructions

For your protection

NOTICE

Read the operating instructions in full before starting up and follow the safety instructions.

- Keep the operating instructions in a place where it can be accessed by everyone
- Ensure that only trained staff work with the instrument.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.

DANGER

Always wear personal protective equipment in accordance with the hazard class of the media being worked with, otherwise there is a risk of splashing liquids, release of toxic or inflammable gases!

- The device and device parts must not be used on humans or animals.
- Do not expose parts of the human body or animals to the vacuum.
- Position the device as specified in the "Commissioning" section and connect the connection cables and interfaces as described.
- Never work with an incorrectly connected or defective device.
- Inhalation of or contact with media such as poisonous liquids, gases, spray mist, vapours, dusts or biological and microbiological materials can be hazardous to user. Make sure that all connections are tight and leak-free when working with such media.
- Prevent the release of the above mentioned substances. The use of suitable extraction systems e.g. fume cupboards is recommended in the place of installation. Take protective measures for personnel and the environment.
- Due to the residual leakage rate of the device, this can result in the release of medium.
- The VC 10 lite/pro vacuum controller is not designed to be installed in potentially explosive atmospheres
- The VC 10 lite/pro vacuum controller is not designed for use with self inflammable substances, substances which are inflammable without air or explosive substances.
- Take suitable precautions to prevent any formation of explosive mixtures. Use inert gas for venting and/or dilution, if necessary.
- Take into consideration interactions or possible chemical or physical reactions when working with media under reduced pressure and elevated temperature.
- Never work with the device under water or underground.
- The device must only be operated with supervision.
- Electrostatic events can occur between the medium and device and lead to a direct hazard.
- Safe operation is only guaranteed with the accessories described in the "Accessories" section.
- The instrument can only be disconnected from the mains supply by pulling out the mains plug or the connector plug.
- The socket for the mains cord must be easily accessible.
- Do not use the device if the AC adapter is defective.
- Follow the operating instructions of the additional equipment carefully (e.g. a rotary evaporator, vacuum pump) with which the VC 10 lite/pro vacuum controller is operated.

- The VC 10 lite/pro vacuum controller may only be operated under the conditions described in the "Technical data" chapter.
- The pressure at the gas inlet and gas outlet must not exceed 1100 mbar. For pressure over 1100 mbar, the total allowable pressure is exceeded and the device no longer displays the pressure accurately.
- Only use flexible hose lines.
- Flexible elements can be compressed under vacuum.
- In the event of a power failure, the integrated venting valve automatically aerates the connected recipient.
- Be aware of your emergency measures in the event of a power failure and ensure that the system is put into a safe state.
- The safety of the user cannot be guaranteed if the appliance is operated with accessories that are not supplied or recommended by the manufacturer or if the appliance is operated improperly, contrary to the manufacturer's specifications or if the instrument or the printed circuit board are modified by third parties.

For protection of the equipment

CAUTION

The voltage stated on the type plate must correspond to the mains-voltage.

- Covers or parts that can be removed from the device without tools must later be refitted to ensure safe operation, so long as no other connection is made at this point. This will prevent the infiltration of foreign objects, liquids etc.
- The device must only be operated with the original plug-in power supply unit supplied.
- Protect the appliance and accessories from bumps and impacts.
- The appliance may only be opened by experts (authorized staff).
- To ensure sufficient cooling of the VC 10 lite/pro vacuum controller, the ventilation slots on the housing must not be covered.
- Use only original spare parts for repair and maintenance in order to ensure the reliable operating condition of the device.
- Pay attention to water condensation inside and outside of the device. If the device is brought from a cold environment, allow the device to warm up.
- Never attach the vacuum controller above a heating bath.
- Ensure that no solids and/or liquids enter the VC 10 lite/pro vacuum controller via the suction power of the pump. This will result in damage to the pressure sensor and the valves.

WARNING

Check the compatibility of the substances used with the device with the materials listed under "Product contact parts"; see chapter "Product contact parts".

- Place the VC 10 lite/pro vacuum controller in an elevated position if possible, so that in the event of a fault, no condensate can collect on the pressure sensor.


Unpacking

Unpacking

- Please unpack the device carefully.
- In the case of any damage a detailed report must be sent immediately (post, rail or forwarder)

Delivery scope

Vacuum controller VC 10 lite/pro	
Universal plug-in power supply unit	 Fig. 2
USB A - micro B cable	 Fig. 3

Vacuum hose 1 m	 Fig. 4
Y connector	 Fig. 5
Operating instructions	
Warranty card	

Intended use

Use

Together with the accessories recommended by , the VC 10 lite/pro vacuum controller is suitable for the controlled evacuation of air (gas) from laboratory devices and equipment (e.g. rotary evaporator or laboratory reactors, and also for classical separation or filtration or drying tasks in the laboratory) in conjunction with a suitable vacuum source (e.g. MVP 10 basic / VACSTAR). The VACSTAR can be operated automatically with the VC 10 pro

vacuum controller. Therefore, it can be used for automatic boiling point recognition, programmed pressure-time curves or programs from the solvent library, for example.

Mode of operation: Stand device

Area of use

Indoor environments similar to that a laboratory of research, teaching, trade or industry area.

The safety of the user cannot be guaranteed:

- If the device is operated with accessories that are not supplied or recommended by the manufacturer;
- If the device is operated improperly or contrary to the manufacturer's specifications;
- If the device or the printed circuit board are modified by the third parties.

WARNING

The device must not be used for:
- **vacuation of biological biotopes,**
- **evacuation of explosive, corrosive or gas volumes,**
- **throughput/use of fluids!**

Useful information

In addition to a high-resolution Al₂O₃ pressure sensor, the VC 10 lite/pro vacuum controller has an integrated vacuum valve and a venting valve. As a result, the VC 10 lite/pro vacuum controller can be used without any additional peripheral devices.

The VC 10 pro provides two basic modes of operation: Two-position control and analogue speed control.

The change between the two modes is performed automatically by plugging or unplugging the connector cable at the "VACSTAR" interface. Depending on the selection made, some submenu items may not be available (greyed out on the display).

With two-position control, parameters from the solvent library can be applied or target values specified in manual mode.

With analogue speed control, an automatic mode as well as the setting of the percentage pump speed can also be preset. Space requirement: To connect all peripheral devices you will need approximately WxDxH 150 x 200 x 155 mm³ space for the device with connected plugs.

Two-position control VC 10 pro (only)

A vacuum is made in the glassware with the help of a vacuum pump. The vacuum pump operates at constant speed which generally cannot be adjusted. Once the target value setting is reached, the integrated vacuum valve interrupts the suction line.


The device will fall slightly short of the target value setting due to the time difference in the millisecond range with the vacuum valve, the target value comparison and the switching of the vacuum valve. Natural leakage from the apparatus now causes another increase in the pressure value of the system, which in turn is detected by the integrated pressure sensor.

If the target value is exceeded, the vacuum valve is reopened and the vacuum pump in operation lowers the pressure. The pressure curve shown on the display fluctuates around the set target value. The pressure difference between activating/deactivating the valve can be set using the vacuum hysteresis value.

The quality of the control depends on the suction power of the pump (rotation speed), the tightness of the system and signal delay time of the electronics used.

Automatic boiling point recognition is not possible when using two-position control.

Schematic view of a two-position vacuum control

Analogue speed control

Imprecision of the two-position control is avoided by using the analogue speed control.

With this kind of control, the speed of the pump, and therefore its suction power, is reduced the closer the measured pressure reaches the target value. Once the target value is reached, the pump only operates according to the leakage rate of the system. This type of control offers quieter operation, more exact vacuum control and a long service life of movable pump parts.

Automatic boiling point recognition is possible with this type of control, i.e. the system reaches and keeps the solvent at its boiling point in automatic mode.


Analogue speed control with the VACSTAR:

The VACSTAR is automatically detected as soon as the connection cable is plugged into both devices. The analogue speed control is enabled, the two-position control is disabled and cannot be selected.

Menu items such as automatic boiling point recognition, programmed pressure-time curves or programs from the solvent library can be used.

Refer to the operating instructions for the VACSTAR vacuum pump.

Schematic view of speed-vacuum control

*An external temperature sensor can be connected to the "PT1000" temperature sensor interface. When the menu view is enabled, the temperature value is shown in the display.

Time sequences can be controlled using a timer.

The operating modes A, B, C and D have preset access options for adjusting and saving target values.

Display settings can be changed using the "Display" menu item.

Access to the device can be restricted by setting up a password in the "Security" menu item.

Common settings such as language, pressure and temperature units can be set via the "Settings" menu item. The display background can also be changed from black to white as well as brightness settings adjusted.

The integrated solvent library calculates the theoretical boiling point at a specified steam temperature value, or vice versa.

The upper and lower limits or the switching points of the vacuum control can be set in the submenu items "Limit values" and "Hysteresis".

*only VC 10 pro

Setting up

NOTICE

Observe the general information, always connect the recipients (load/vacuum vessel/glass cooler) at their highest points to the suction line. This will prevent the risk of fluid entering into the vacuum controller or pump.

Install a separator (e.g. Woulff bottle) in front of the intake mounting to protect against the ingress of liquid.

In the case of induced solvent vapours, a vapour condenser (accessory of pump) installed downstream of the pump helps to condense the vapours and prevents release into the atmosphere. Liquid in the pump chambers deteriorates the pump characteristics.

Attach the VC 10 lite/pro vacuum controller to a stand (d=16 mm) close to the pump.

Secure the controller with the stand locking screw (N).


Fig. 5

To do so, use a stand with diameter 16 mm, for example, stand VC 10.400, VSS 1 vacuum safety set, the stand on VACSTAR vacuum pump or attach the controller to the RV 8 rotary evaporator, or at LR 1000 laboratory reactor.

Fig. 6

VC 10 lite/pro at stand VC 10.400

Fig. 7

VC 10 lite/pro to VSS 1 with MVP 10 basic

Fig. 8

VC 10 lite/pro to VACSTAR

Fig. 9

VC 10 lite/pro to RV 8

Connecting the interfaces

Connect the vacuum hoses, pay attention to the following connections (also see **Fig. 1**):

J: Load

Hose connection d=10 mm to the system to be evacuated (recipient).

K: Pump

Hose connection d=10 mm to the vacuum source (e.g. pump).

L: Vent

Ventilation connection d=10 mm, also suitable for the connection of inert gas.

Secure the hoses with hose clamps, if necessary.

Lay the vacuum hoses according to your device configuration.

Now connect the connecting cables according to your system configuration:

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Connect the VC 10 lite/pro vacuum controller to a PC using the USB 2.0 cable (A/Micro-B).

Software updates to any existing equipment can be installed using the FUT software tool.

M: RS 232 Interface

You can connect the VC 10 lite/pro vacuum controller to a PC using an RS 232 interface cable. The pump can be operated in conjunction with other devices using labworldsoft® laboratory software. For more information see Chapter "Interfaces and outputs".

P: PT 1000 temperature sensor (VC 10 pro)

Connect the optional PT 1000.60, PT 1000.61 or PT 1000.70 temperature sensor (accessory). The measured temperature (e.g. tempering bath) is shown on the display.

The value displayed is for information purposes only and is not processed logically.

Note that the value displayed is dependent on the positioning of the sensor in the medium, its mixing, as well as the calibration of the probe and can deviate from other measured temperature values.

Incorrect or fluctuating measured values can be displayed by using a heating bath, e.g. by laying the sensor on the heating bath floor, instead of positioning in the medium.

Q: VACSTAR membrane vacuum pump (Mini DIN)

You can connect the vacuum controller and the VACSTAR vacuum pump to the analogue connection cable (accessory) for exact speed-controlled vacuum control. The vacuum controller detects the pump and switches to speed-vacuum control mode. The two-position control is disabled. The pump speed is regulated depending on the measured pressure.

R: Connector for 24 W AC adapter

Connect the universal AC adapter (included) to the socket on the rear. Observe the country-specific variants of the AC adapter!

Explanation of abbreviations used in the following schematic drawings and configurations:

- 1 Recipient (load, e.g. rotary evaporator, reactor)
- 2 Separator (Woulff bottle)
- 3 Venting valve
- 4 Pressure sensor
- 5 Vacuum valve/ball valve
- 7 Vapour condenser
- 8 Connection cable (analogue)
- 9 Vacuum source (pump, in-house vacuum)
- 10 VC 10 lite/pro.300 Check valve

Two-position control

Vacuum pump/vacuum source with VC 10 pro vacuum controller.

Once the target value is reached, the suction line is closed by the vacuum valve integrated in the VC 10 pro vacuum controller. The pump speed can be changed by manually adjusting the speed setting.

For pumps which have the speed setting option, it is recommended to work at low speed.


NOTICE

When using the VACSTAR vacuum pump: The pump must be operated in mode "B" (restart after voltage drop).

The vacuum source can also be a central in-house vacuum network.

Illustration two-position control - vacuum pump/vacuum source with vacuum controller VC 10 pro

Vacuum source with two VC 10 pro vacuum controllers.


Once the target value is reached in a recipient, the suction line is closed by the vacuum valve integrated in the VC 10 pro vacuum controller.

A sufficiently stable vacuum source (in-house vacuum, pump) can serve multiple consumers.

For pumps which have the speed setting option, it is recommended to work at low speed.

A second consumer can be controlled independent of the first using a second vacuum controller.

Illustration two-position control - vacuum source with two vacuum controllers VC 10 pro

Speed control


VACSTAR vacuum pump with VC 10 lite/pro vacuum controller

Automatic adjustment of the “analog speed control” mode if the analogue connection cable (Pos. 8) is connected to the pump with the vacuum controller.

Once the target value is reached the pump speed is automatically reduced to a minimum “0 rpm”. Once the target value is reached,

the pump operates according to the leakage rate of the vacuum system. Automatic boiling point recognition is only possible with the VC 10 pro .

Illustration speed control - Vacuum pump VACSTAR with vacuum Controller VC 10 lite/pro

Commissioning

Check whether the voltage given on the type plate corresponds to the available mains voltage.

If these conditions have been met, the machine is ready for operation when the mains plug is plugged in.

If these conditions are not met, safe operation is not guaranteed and the machine could be damaged.

Also pay attention to the ambient conditions listed in the “Technical data”.

NOTICE

The start screen appears for a few seconds when the device is switched using the “Standby ON/OFF” switch (F) The device name and the software/firmware versions are also displayed.


An information prompt then appears to download the firmware update tool.

The following working screen appears automatically in the display.


Display of two-position control, e.g. with MVP 10 basic.

Hysteresis display for two-position control


Display speed control: VACSTAR to analogue connection

No hysteresis display for speed control

Explanation of symbols on the working screen

The symbols displayed change depending on the status and settings of the vacuum controller:

Key Lock key

This symbol means that the function of the key and of the rotating/pressing knob for controlling the vacuum controller are disabled. The symbol no longer appears if the functions are enabled once again by pressing the lock key a second time (min. 1 sec.).

A Operating mode

This symbol indicates the operating mode currently selected (A, B, C, D).

USB

This symbol means the vacuum controller is communicating via a USB cable. The symbol no longer appears if no USB cable is being used for communicating with the station.

Temperature sensor (VC 10 pro)

This symbol appears when the temperature display is active in the display.

Timer This symbol appears when the timer display is active in the display.

PC PC control

This symbol means that the vacuum controller is connected to a computer and the device is being controlled from the computer.

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This symbol means that the vacuum controller is connected to a computer and the device is being controlled from the selected program.

Continuous mode

This symbol indicates that the vacuum controller is in continuous mode.

Vent

This symbol indicates that venting is enabled.

Vacuum

This symbol indicates that the current pressure is less than the existing atmospheric pressure.

Detection

This symbol is displayed during the automatic boiling point search, the system pressure is reduced further.

Boiling

This symbol indicates that a boiling point has been found automatically, the system pressure is kept constant.

Refrigerated circulator

This symbol indicates that refrigerated circulator is operating.

Navigation menu

Fig. 9


- Press the "Menu" key (C).
- Select the menu by turning the rotating/pressing knob (D) to the right or left to select the desired menu or sub-menu, which can then be selected by pressing the rotating/pressing knob (D).
- Press or turn the rotating/pressing knob (D) again to select the desired menu option and edit the values or settings, or activate/deactivate a function.
- Turn the rotating/pressing knob (D) to "OK" and press the "Back" (B) or "Menu" key (C) to end the procedure and return to the previous menu or working screen.
- Navigating in Menu "Graph": Press the "Back" (B) at main screen to navigate to graph menu, showing real-time pressure change on graph.

NOTICE

The menu option activated is highlighted yellow in the display. Menus which are grayed out are inactive.

Navigation menu: Press the "Menu" key (C) and turn the rotating/pressing knob (D)
 Press the "Back" key (B) or the "Menu" (C)

NOTICE


If you press the "Menu" key (C), the system skips directly back to the working screen.
 If you press the "Back" key (B), the system skips back to the previous display.

Menu structure


Menu	Sub menu	Option	Action	Factory setting	
Pump	Modes	Automatic (only VC 10 pro)		Deactivated	
		Manual		Activated	
		Pump %		Deactivated	
		Program		Deactivated	
		Start After		-	
	Refrigerated circulator	Temperature		-	
		Pump speed		-	
		Cool now		-	
		Stop now		-	
		Start after		-	
	Cleaning	Start after		05:00 [mm:ss]	
		Duration		06:00 [mm:ss]	
		Pump speed		20 %	
		Clean now		Deactivated	
	Limits	Manual	Minimum		1 mbar
			Maximum		1100 mbar
		Pump %	Minimum		0 %
			Maximum		100 %
	Hysteresis (VC 10 pro)	Automatic		Activated	
		Manual		Deactivated	
	Solvent library (VC 10 pro)	Acetonitril		-	
		N-Amylalkohol		-	
		n-Pentalon		-	
		Benzol		-	
		...		-	
Actions after stop	Open venting valve		Activated		
	Cleaning		Deactivated		
Actions after start	Stop refrigerated circulator		Deactivated		
	Turn on refrigerated circulator		Deactivated		
Leakage Test	Leakage rate		3 mbar		
	Tolerance		1mbar		
	Test point		500 mbar		
	Leakage test duration		1 minute		
				-	
Temperature (VC 10 pro)	Sensor temperature		-		
	Display		Deactivated		
	Calibration		20.0 °C		
Timer	Timer function	Time		00:00:01 [h:mm:ss]	
		Beep after stop		Activated	
	Display	Activate		Deactivated	
Operating mode	A		Deactivated		
	B		Activated		
	C		Deactivated		
	D		Deactivated		
Display	Timer		Deactivated		
	Temperature (VC 10 pro)		Deactivated		
Programs	Program 01		-		
	...		-		
	Program 10	Select		-	
		Edit	Edit, insert, delete	-	
		Delete		-	
Last measurement	Save as		-		
	Edit	Edit, insert, delete	-		
	Delete		-		
Safety	Password		0 0 0		
Service	Valves	Venting valve	Close, open	Open	
		Vacuum valve (VC 10 pro)	Close, open	Close	
	Pump		Deactivated		
Refrigerated circulator		Deactivated			
Settings	Languages	English		Activated	
		German		Deactivated	
	Units	Pressure	mbar		Activated
			hPa		Deactivated
			mmHg		Deactivated
			Torr		Deactivated
		Temperature (VC 10 pro)	°C		Activated
		°F (not available for Japan)		Deactivated	
	Display	Background	Black		Activated
			White		Deactivated
	Brightness		100 %		
	Firmware update info		Activated		
	Sound	Volume		10 %	
		Key tone		Deactivated	
	Factory setting	Menu values	Reset to factory setting	-	
Programs			-		
Communication	All		-		
	Device name		VC 10 lite/pro		
Information	Firmware update info		-		
	Display version		-		
	Logic version		-		
	Operating mode		B		
	Pressure max		1100 mbar		
	Pressure min		1 mbar		
	Pump % max		100 %		
	Pump % min		0 %		

Menu details

Pump

Modes

Automatic (VC 10 pro):


In the "Modes" menu, the user can enable automatic boiling point recognition by selecting the "Automatic" menu item. No other parameters must be set.

The boiling point is detected automatically. For rotary evaporators with heating bath, it must be ensured that the tempering medium and solvents have a constant temperature (e.g. 60 °C).

Manual:

In the "Manual" menu item, the target value can be specified (e.g. in "mbar"). The system is evacuated until the target value.

Refrigerated circulator


Start after...:

Specifies a delay before the refrigerated circulator start operating.

Temperature:

Determines the temperature, that the refrigerated circulator is operated.

Cleaning

The cleaning of the pump ensures that solvent residues are removed in the pump after the measurement. For this purpose, the pump is operated at a low level while the system is vented. The air flowing through the pump ensures that the solvent residues are discharged from the pump. If the cleaning option is selected, this starts once the measurement has finished (see "Actions after stop").

Pump %

By selecting the "Pump %" menu item, the pump can be operated continuously with a running performance of between 100 % and 1 %.

Program

Under the "Programs" menu, 10 user-defined pressure-time profiles can be created. The last measurement taken can also be viewed in this menu and saved as a program.

The target value and the running performance can only be entered within the set limits (see menu item "Limits").

Changing to the working screen and pressing the rotating/pressing knob starts the process. Pressing the knob again stops the process.

NOTICE

If no speed-controlled pump with analogue connection is connected to the "VACSTAR" interface, the system can only work in manual mode and in "two-position control" mode. The menu items "Auto" and "Pump %" are greyed out and cannot be activated.

Pump speed:

Determines the speed in rpm, that the refrigerated circulator is operated

Cool now:

If you want to perform a cooling immediately, you can select it with this setting.

Stop now:

If you want to stop the cooling immediately, you can select it with this setting.

NOTICE

The "Cleaning" menu item increases the service life of the pump, since solvent residues do not remain inside the pump.

You can also clean the pump without using the corresponding menu item. To do this, remove the hoses and operate the pump at the end of a working cycle (test, test series, end of the day ...) at idling speed.

Start after...:

Specifies after which minimum duration of a measurement the cleaning should take place. The default time is 5 minutes. This means that a measurement must be run for at least 5 minutes, so that cleaning is performed at the end of the measurement.

Duration:

Specifies how long the cleaning process should be performed.


Pump speed:

Determines at what power (in %) the pump is operated during cleaning. If the power is set too high this can lead to a pressure drop in the system.

Clean now:


If you want to perform a cleaning, you can select it with this setting.

Limits

In this menu, limits for the target value can be set.

Hysteresis (VC 10 pro)

With the "Hysteresis" option, the switching frequency and control accuracy can be influenced when in a steady state. The hysteresis value determines the upper and lower limits of the target value and therefore the closing and opening of the pump and valve (only for two-position control).

Automatic: Hysteresis is always 10 % of the actual pressure.


Manual: Manual value specification of hysteresis.

NOTICE


Limits which are too narrow shorten the service life of the devices.

Solvent Library (VC 10 pro)


In the "solvent library", the most common solvents are assigned a function to determine the steam temperature at the set pressure and vice versa. This provides assistance when setting these parameters in manual mode, for example, evaporation processes with a rotary evaporator.

Actions after stop


This menu allows the user to define actions that should be performed when a test has ended.

Actions after start

This menu allows the user to define actions that should be performed when a test has started.

Leakage test

Leakage test: User to confirm the vacuum connection. Press ok if the connection is good.

Set Value: User to input the required leakage rate with tolerance. Default leakage rate is 3mbar with tolerance 1mbar.


Test point: User to input the pressure which the leakage test will

be executed at. Default value is 500 mbar.

Duration: User to input duration of the leakage test. Default value is 1 minute.


In progress: The pump will start running to reach the test point. (Graph below: Phase 1 and 2)

The leakage measurement will start once the pressure has reached the interested test point. Meanwhile, the timer will start counting down from the pre-set duration. (Graph below: Phase 3)

Temperature (VC 10 pro)

Temperature sensor

In the "Temperature" menu the user can specify that the temperature sensor is displayed on the display/working screen. A tick (✓) shows that the option is activated.


The precondition for this is that a temperature sensor is connected to the vacuum controller. If no temperature sensor is connected, or an error is present or the temperature rises above 350 °C, three dashes are shown as the temperature value.

NOTICE Please comply with the temperature measurement range of the external temperature sensor specified in the "Technical data" section. The temperature display can be in either °C or °F (see the "Settings" menu).


Calibration

In order to calibrate the external temperature sensor, it must first be connected to the vacuum controller. The temperature sensor is then immersed in a vessel with water, for example. The water temperature is determined using a second, already calibrated temperature measuring device. The measured temperature must then be entered in the menu under "Temperature calibration". Once the entry has been confirmed with "OK", the external temperature sensor is calibrated.

Timer

Timer functions


In the "Timer" menu the user can specify that the timer is displayed on the display/working screen. A tick (✓) shows that the

option is activated. This setting allows the user to specify the actual time for the evacuation procedure.

A default time can also be set for the timer. This setting allows the user to start evacuation for a standard time. The device can emit a signal (beep) after the set target time has elapsed.

NOTICE The user can stop the evacuation procedure before expiry of the set time. In this case the countdown of the timer is interrupted.

Operating mode

Operating mode A:

In this operating mode, the set target value is not saved when the current run comes to an end or the device is switched off.

Operating mode B:

In this operating mode, the set target value is saved when the

current run comes to an end or the device is switched off, and the value can be changed.

Operating mode C:


In this operating mode, the set target value is saved when the current run comes to an end or the device is switched off, and the value cannot be changed.

Operating mode D:

In this operating mode, the target value is applied if it does not change after 3 seconds, or when the yellow background is displayed behind the target value. In operating modes A and B, the target value is applied immediately when the target value is changed by turning the knob.


In operating mode C, the target value cannot be changed.

Display


In the "Display" menu the user can specify what information will be displayed on the main screen.

Graph

In the "Graph" menu, user can observe the real-time pressure change on the shown graph.

Programs

Select:

To select a program, use the rotating/pressing knob and press "Select". Only programs that have at least one program segment can be selected.

When a program has been successfully selected, this is indicated by a tick mark (✓). The vacuum controller is then in program mode. In the main screen, the program view can then be selected by turning the knob to the right. Turning the knob to the left returns the user to the main screen.

Start:

To start a selected program, the push knob must be pressed in the main screen.

Edit:

Edit the selected program parameters. Start to edit the selected program parameters by pressing on menu option "Edit" with rotating/pressing knob. The user can edit, delete or insert one selected program segment in the program.

Delete:

Deletes the selected program. If a selected program is deleted by pressing on menu option "Delete" with rotating/pressing knob, all the program parameters will be emptied. The tick (✓) disappears.

Details for editing the program

No.	Pressure	hh:mm:ss
01	600 mbar	00:00:10
02	550 mbar	00:00:20
03	500 mbar	00:00:30
04	450 mbar	00:00:40
05	400 mbar	00:00:50
06	350 mbar	00:01:00
07	300 mbar	00:01:10
08	250 mbar	00:01:20
edit insert delete		

In this program, user can define up to 10 segments. The selected segment is highlighted. Then, the user can edit, delete or insert a segment in this program. The program is save automatically.

Edit:

If the background of a selected value is yellow, the user can change the setting of the pressure value or the time value.

Delete:


When a highlighted segment is deleted, the subsequent segments move upward, closing the gap created by the deleted segment.

When a program is being edited, following screen appears.

Sample for editing the program

Example of saving the last measurement


Safety


In the "Password" menu, the user can protect the vacuum controller settings using a password. (factory setting: 000)

Service

In the "Service" menu, the valves, the pump or the refrigerated circulator can be operated individually and also checked to ensure that they are operating correctly. If no VACSTAR is connected, the menu item "Pump" is greyed out and not active.

Settings

Sound:

The "Sound" option allows the user to activate/deactivate the key-press sound and to set the volume.

Factory settings:

Select the "Factory settings" option by turning and pressing the rotating/pressing knob. The system requests confirmation to re-create the factory settings. You can choose if you want to reset only the "Menu values", or only the "Programs". Alternatively, you can reset all by selecting "All". Pressing the "OK" key resets the system settings to the original standard values set at dispatch from the factory (see "Menu structure" illustration).

Communication:

The "Device name" option allows the user to input a device name. The device name appears on the opening screen when the device is switched on.

This can be useful when using multiple devices with different settings. This also facilitates the identification of each device when communicating with a PC via USB, RS232 and Bluetooth, for example.

Information:

The "Information" option offers the user an overview of the most important system settings of the vacuum controller VC 10 lite/pro.

Languages:

The "Languages" option allows the user to select the desired language by turning and pressing the rotating/pressing knob (D). A tick (✓) indicates the language that is set for the system.

Units:

The "Units" option allows the user to select the desired unit for displaying the temperature and pressure by turning and pressing the rotating/pressing knob. For the temperature, the choice is between "°C" or "°F". For the pressure, the choice is "mbar", "hPa", "mmHg" or "Torr". A tick (✓) indicates the unit that is set for the system.

Display:

The "Display" option allows the user to change the background color and brightness of the working screen.

Interfaces and outputs

The device can be operated in "Remote" mode via an RS 232 or USB interface using the laboratory software labworldsoft®.

The RS 232 interface at the back of the device is fitted with a 9-pole SUB-D port which can be connected to a PC. The pins are assigned serial signals.

NOTICE

Please comply with the system requirements together with the operating instructions and help section included with the software.

USB interface

The Universal Serial Bus (USB) is a serial bus for connecting the device to the PC. Equipped with USB devices can be connected to a PC during operation (hot plugging). Connected devices and their properties are automatically recognized.

Use the USB interface in conjunction with labworldsoft® for operation in "Remote" mode and also to update the firmware.

Installation

First, download the latest driver for devices with USB in-terface from <http://www.com/lws/download/usb-driver.zip> and install the driver by running the setup file. Then connect the device through the USB data cable to the PC.

Serial interface RS 232 (V 24)

Configuration:

- The functions of the interface connections between the device and the automation system are chosen from the signals specified in EIA standard RS 232 in accordance with DIN 66 020 Part 1.
- For the electrical characteristics of the interface and the allocation of signal status, standard RS 232 applies in accordance with DIN 66 259 Part 1.
- Transmission procedure: asynchronous character transmission in start-stop mode.
- Type of transmission: full duplex.
- Character format: character representation in accordance with data format in DIN 66 022 for start-stop mode. 1 start bit; 7 character bits; 1 parity bit (even); 1 stop bit.
- Transmission speed: 9600 bit/s
- Data flow control: none
- Access procedure: data transfer from the device to the computer takes place only at the computer's request

Command syntax and format

The following applies to the command set:

- Commands are generally sent from the computer (Master) to the stirrer machine (Slave).
- The stirrer machine sends only at the computer's request. Even fault indications cannot be sent spontaneously from the device to the computer (automation system).
- Commands are transmitted in capital letters.
- Commands and parameters including successive parameters are separated by at least one space (Code: hex 0x20).
- Each individual command (incl. parameters and data) and each response are terminated with Blank CR LF (Code: hex 0x20 hex 0x0d hex 0x20 hex 0x0A) and have a maximum length of 50 characters.
- The decimal separator in a number is a dot (Code: hex 0x2E).

The above details correspond as far as possible to the recommendations of the NAMUR working party (NAMUR recommendations for the design of electrical plug connections for analogue and digital signal transmission on individual items of laboratory control equipment, rev. 1.1).

The NAMUR commands and the additional specific commands serve only as low level commands for communication between the stirrer machine and the PC. With a suitable terminal or communications programme these commands can be transmitted directly to the stirrer equipment.

The software package, labworldsoft®, provides a convenient tool for controlling the equipment and collecting data under MS Windows, and includes graphical entry features, for motor speed ramps for example.

The following table summarises the (NAMUR) commands understood by the control equipment.

NAMUR Commands	Function
IN_NAME	Read device name
IN_PV_3	Read PT 1000 value
IN_PV_66	Read current pressure value
IN_SP_66	Read pressure target value
IN_SP_70	Read hysteresis target value
IN_MODE_66	Read current evacuation mode
IN_ERROR	Read error status
OUT_SP_66	Set pressure target value
OUT_SP_70	Set hysteresis value
OUT_MODE_66	Set evacuation mode
START_66	Start evacuation
STOP_66	Stop evacuation
RESET	Switch to normal operating mode

PC 1.1 Cable (Device to PC)

Required for connecting the 9-pin socket to a PC.


Fig. 10

Connection VC 10 lite/pro - PC





Fig. 11

Connection Refrigerated circulator

PC1.3 Cable (Device - Refrigerated circulator)

Connection (Device - Refrigerated circulator)

Connection cable

This cable is required to connect the VC 10 lite/pro to the VACSTAR (speed mode).


Fig. 12

Maintenance and cleaning

The device is maintenance-free. It is only subject to the natural wear and tear of components and their statistical failure rate.

Cleaning

- For cleaning disconnect the mains plug!
- Use only cleaning agents which have been approved by to clean devices.
 - These are water (with tenside) and isopropanol.
- Wear protective gloves during cleaning the devices.
- Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.
- Do not allow moisture to get into the device when cleaning.
- Before using another than the recommended method for cleaning or decontamination, the user must ascertain with that this method does not destroy the device.

Spare parts order

When ordering spare parts, please give:

- device type.
- serial number, see type plate.
- position number and description of spare part, see www.ika.com.
- software version

Repair

Please only send devices in for repair that have been cleaned and are free of materials which might present health hazards.

For repair, please request the "Safety Declaration (Decontamination Certificate)" from or use the download printout of it from website at www.ika.com.

If your appliance requires repair, return it in its original packaging. Storage packaging is not sufficient when sending the device - also use appropriate transport packaging.

Error codes

Error:

The fault is shown by an error message in the display as following if the error occurs, e.g. Error 4. Proceed as follows in such cases:

- Switch the device switch off,
 - carry out corrective measures,
 - switch the device switch on again.
- ⇒ Measurement is stopped.

Warning:

- Measurement is stopped.
- Warning is displayed.
- Warning can be confirmed.

Information:

- Measurement continues.
- Information is displayed.
- Information can be confirmed.

Description	Detection	Reason	Action
Information No Pressure Change	Pressure gradient analysis does not show any deviation.	Pump is not running. Vacuum hose isn't connected. Failure in hose system. Recipient isn't closed.	Check pump control cable connection. Check pump power supply. Switch on pump. the position. Check correct vacuum hose connection (IN/ OUT/VENT). Close recipient.
Information System Not Tight	Pressure gradient analysis shows pressure deviation, but set value cannot be reached.	Vacuum hose connection isn't tight. Recipient isn't tight. Vacuum pump power isn't sufficient. Venting valve leakage. Boiling point of solvent is reached before.	Check vacuum hose connection. Check recipient. Check technical data of the pump. Call service department. Check set value.
Information Venting Error	Pressure gradient analysis does not show any deviation after pressing the "Venting" button.	Recipient wasn't evacuated at all. Failure in hose system. Venting valve leakage. Front foil button is damaged.	Evacuated recipient. Check correct vacuum hose connection (IN/ OUT/VENT). Call service department.
Pump not connected	Connection to Vacstar(speed mode) not present.	Damaged or no cable connected	Check connection cable is present and well connected.
Warning Boiling Point Error	Boiling point cannot be detected. Pressure gradient analysis does not show boiling point.	Heating bath does not heat up. Solvent with extreme low boiling point is used.	Check heating bath. Solvent distillation only manually.
Warning External Temperature Is Too High	Measured temperature is too high. External temperature is higher than maximum value	Medium temperature is too high. External sensor is broken.	Check the medium temperature. Call service department.
Warning External Temperature Is Too Low	Measured temperature is too low. External temperature is lower than minimum value.	Medium temperature is too low. External sensor is short-connected.	Check the medium temperature. Call service department.
Warning Control Pump State Error	Pump is disconnected during analog speed control mode. Pump is connected during 2-point control mode.	Pump is disconnected during analog speed control mode. Pump is connected during 2-point control mode	Connect the pump. Disconnect the analog pump.
Warning RS232 Communication Error USB Communication Error	Watchdog time elapsed.	labworldsoft® uses wrong device for VC 10 lite/pro. Wrong tools are used to communicate with VC 10 lite/pro. Unstable connection. Watchdog was set in device or via PC command.	Check the labworldsoft® setting. Check the PC tools. Check the connection. Disable watchdog when not necessary. Increase sample rate.

Description	Detection	Reason	Action
Warning Temperature Calibration Error	Temperature calibration value is out of range.	Wrong simulator resistance is chosen. Set wrong temperature calibration value.	Check the simulator resistance. Call service department.
Warning Pressure Calibration Error	Pressure calibration value is out of range.	Input wrong vacuum actual value when calibration. Set wrong vacuum calibration value. Sensor is broken.	Calibrate again. Check the set calibration value. Call service department.
Error 3 Device Temperature Error.	Analysis of internal PCB temp. Sensor. Inside temperature reaches limit value.	Room temperature > 40 °C. PCB or valve is broken.	Switch off device for cool down. Call service department.
Error 9 Logic Storage Error. Display Storage Error.	Read or write internal memory error.	Read or write internal memory error.	Call service department.
Error 68 Pressure Out Of Range.	Pressure is out of range. Pressure is bigger than room pressure.	Pressure in the recipient is too high. Sensor is broken.	Check the air flow of pump and vent recipient. Call service department.
Error 69 Pressure Sensor Error.	Pressure sensor output is too low.	Sensor is not connected. Sensor is broken.	Check the sensor connection. Call service department.
Error 71 Internal Communication Error	Internal Watchdog elapsed.	Communication lost between logic board and display board.	Restart device (unplug mains cable, plug in mains cable). Check communication between logic board and display board. Call service department.

If the actions described fails to resolve the fault or another error code is displayed then take one of the following steps:

- contact the service department,
- send the instrument for repair, including a short description of the fault.

Accessories

Accessories see www.ika.com.

Product contact parts

Designation	Material
Connection spigots	PP
Distributor	PPS
Pressure sensor	FPM/AL2O3
Vacuum valve	PEEK; EPDM; FKM
Venting valve	PEEK; EPDM; FKM

Technical data

	Unit	VC 10 lite	VC 10 pro
Connection diameter suction side	mm	8	
Connection diameter pressure side	mm	8	
Connection diameter venting side	mm	8	
Input pressure min.	mbar	1	
Input pressure max.	mbar	1050	
Boiling point detection		-	yes
Solvent library		yes	
Two-point control		-	yes
Analog speed vacuum control		yes	
Display		TFT	
Pressure unit / scale		mbar, hPa, mmHg, Torr	
Vacuum sensor		yes	
Vacuum sensor type		ceramic Al ₂ O ₃	
Pressure max. for pressure sensor	bar	1.6	
Measurement range (absolute) min.	mbar	1	
Measurement range (absolute) max.	mbar	1100	
Control range min.	mbar	1	
Control range max.	mbar	1100	
Resolution pressure	mbar	1	
Measurement uncertainty	mbar	1	
Medium temperature (gas) min.	°C	5	
Medium temperature (gas) max.	°C	40	
Vacuum valve		-	yes
Venting valve		yes	
Connection for ext. temperature sensor		-	PT 1000
Temperature unit		-	°C/°F
Temperature measuring range min.	°C	-	-10
Temperature measuring range max.	°C	-	200
Temperature measurement resolution	K	-	1
Accuracy of temperature measurement	K	-	±1
Timer		yes	
Time setting min.	s	1	
Time setting max.	min	6000	
Vacuum speed control interface		VACSTAR	
Material in contact with medium		Al ₂ O ₃ , PTFE, FPM, PPS	
Housing material		PBT	
Fastening		stand/clamp	
Fastening diameter	mm	16	

	Unit	VC 10 lite	VC 10 pro
Dimensions (B x H x T)	mm	95 x 150 x 110	
Weight	kg	1.5	
Permissible ambient temperature min.	°C	5	
Permissible ambient temperature max.	°C	40	
Permissible relative humidity	%	80	
Protection class according to DIN EN 60529		IP 20	
RS 232 interface		yes	
USB interface		yes	
Voltage	V	100 - 240	
Frequency	Hz	50/60	
Power input	W	24	
Power input standby	W	2	
DC Voltage	V	24	
Current consumption	mA	1000	

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IKA Japan K.K.
Phone: +81 6 6730 6781
eMail: info_japan@ika.ne.jp

INDIA

IKA India Private Limited
Phone: +91 80 26253 900
eMail: info@ika.in

UNITED KINGDOM

IKA England LTD.
Phone: +44 1865 986 162
eMail: sales.England@ika.com

VIETNAM

IKA Vietnam Company Limited
Phone: +84 28 38202142
eMail: sales.lab-vietnam@ika.com

THAILAND

IKA Works (Thailand) Co. Ltd.
Phone: +66 2059 4690
eMail: sales.lab-thailand@ika.com

TURKEY

IKA Turkey A.Ş.
Phone: +90 216 394 43 43
eMail: sales.turkey@ika.com

Discover and order the fascinating products of IKA online:
www.ika.com

IKA

designed for scientists

IKA Vacstar digital

Operating instructions
Source language: German

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Fig. 1

Pos. Description (EN)

A	Handle
B	Display
C	Rotating/pressing knob
D	Main switch (on/off)

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Fig. 1

Pos. Description (EN)

E	Handle securing screw
F	Vacuum controller VC 10 stand hole
G	Stand securing screw
H	Fan / ventilation slit
I	Hose connection for suction line d= 8 mm (INLET)
J	Accessory fastening (emission condenser)
K	Connection for vacuum controller VC 10 (Mini DIN)
L	USB interface
M	RS 232 interface
N	Water valve connection
O	Hose connection for pressure line d= 8 mm (OUTLET)
P	Power supply cable connection

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	Page
Device setup	2/3
Declaration of conformity	5
Explication of warning symbols	5
Warranty	5
Safety instructions.....	6
Unpacking	7
Correct use	8
Useful information	8
Setting up.....	9
Commissioning	12
Interfaces and outputs	13
Maintenance and cleaning.....	14
Error codes	15
Accessories	15
Product contact parts.....	15
Technical data.....	16

Declaration of conformity

We declare under our sole responsibility that the product to which this declaration relates is in conformity with directives 2014/35/EU, 2006/42/EC, 2014/30/EU and 2011/65/EU and conforms with the following standards or normative documents: EN 61010-1, EN 60529, EN 61326-1 and EN ISO 12100.

A copy of the complete EU Declaration of Conformity can be requested at sales@ika.com.

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Warranty

According to **IKA**'s Terms and Conditions of sale and delivery, this product is covered by a warranty for a period of 24 months. Please contact your dealer for any warranty claims. If you wish, you can send the device directly to our factory. Please provide the sales invoice and state the reasons for your guarantee claim. In this case, you are responsible for shipping and handling costs.

The warranty does not cover wearing parts, nor defects that are the result of improper use, insufficient care and maintenance or failure to follow the instructions in this operating manual.

Safety instructions

For your protection

NOTE Please read the instruction manual in full before use and follow the safety instructions.

- Keep the instruction manual in a place where it can be accessed easily.
- Ensure that only trained staff use the device.
- Be sure to comply with all safety instructions, directives and all matters of health, safety and accident prevention in the workplace.
- The device and all parts of the device must not be used on people or animals.

DANGER Always wear personal protective equipment in accordance with the hazard class of the media being worked with - otherwise there are dangers due to spraying liquids, and release of toxic or flammable gases.

- Do not expose human or animal body parts to vacuum.
- Do not work with the device underwater or underground.

Device configuration

DANGER The vacuum pump IKA Vacstar digital is not designed to be set up in potentially explosive atmospheres.

- Set up the device in accordance with the chapter "Setting up" and connect the connection lines and interfaces as described.
- Set up the device on a stable, even, non-flammable surface.

- Never work with a faulty or incorrectly connected device.
- Set up the vacuum pumps in appropriate and functioning extractor hoods, or set up an exhaust gas line to them. Ensure that the exhaust gas line cannot be kinked! The maximum permitted exhaust gas line length is 2 metres.
- Occurrence of explosive mixtures; if appropriate connect inert gas for ventilation and/or thinning.

Accessories

- Safe operation can only be ensured when working with accessories as described in the "Accessories" section.
- Carefully observe the operating instructions for additional devices (e.g. rotary evaporators, vacuum controllers), with which the vacuum pump IKA Vacstar digital is operated.
- The pressure at the gas inlet and outlet must not exceed 1100 mbar.

- Elastic elements may be pressed together under vacuum.
- Only use flexible hose lines.
- Observe your emergency measures for power failure and ensure that the device is put in a safe state (see the chapter: Commissioning, Operating mode).

Work with the device

DANGER You can prevent the release of solvent vapours into the atmosphere using a downstream emission condenser.

The vacuum pump is not suitable for use with self-igniting materials, materials that are flammable without an air supply, or explosive materials.

Do not operate the pump while it is open. Otherwise there is a danger of severe injuries as a hand may be pulled into the moving parts.

WARNING Inhaling or coming into contact with media such as poisonous liquids, gases, spray mists, vapours, dusts or biological materials can endanger the health of the user. Ensure that all connections are well sealed and free of leaks if you are working with such media.

- The vacuum pump IKA Vacstar digital must only be operated under the conditions described in the chapter "Technical data".
- Prevent release of the materials listed above. Take measures to protect staff and the environment.
- Pay attention to possible interactions or chemical or physical

reactions when working with media at reduced pressure and increased temperature.

- There can be electrostatic processes between the medium and the device which can lead to direct danger.
- Some medium may be released due to the residual leakage rate of the device.
- Before commissioning check that all the housing parts are present and fastened to the device.
- Do not lift the pump if the handle is loose and the handle securing screw (E) is loose or missing!
- Only switch the pump on if the pump is standing vertically.
- Connect the hose connections (INLET-OUTLET) and interfaces in accordance with the labelling on the device and the operating instructions.
- Please note that the temperature of the evacuated medium is below its ignition temperature. The pumping process (compression) increases the temperature of the medium additionally.
- Ensure that vapours containing solvent can be sucked into the pump.
- Do not use the pump to create pressure.
- Observe the permitted pressure at the inlet and outlet side; see the chapter "Technical data".
- The gas flow must only be regulated/throttled in the suction-side line.
- Use a solenoid valve or a check valve in the suction line if there are several load units.
- When using an emission condenser ensure that the coolant is flowing freely.

Protection of the device

WARNING Ensure that the substances used with the device are compatible with the materials of the device parts that come into contact with the product; see the chapter "Product contact parts".

CAUTION The specified settings on the rating plate must coincide with the actual power supply.

- The device is only disconnected from the power supply when the power or device plug is removed.
- The socket for the power supply cable must be easy to access.
- Covers or parts that can be removed from the device without tools must later be refitted to ensure safe operation, unless anything else is connected at this point. This will prevent the infiltration of foreign objects, liquids and other contaminants.
- Unwanted liquids that have entered can be removed by extracting air during no-load operation.
- Avoid knocks or impact to the device.



- The device must only be opened by trained, skilled personnel (who have been authorised to do so).
- The device must not be modified without authorisation from IKA.
- To ensure sufficient cooling for the vacuum pump IKA Vacstar digital, the ventilation slits on the housing must not be covered.
- Only use original replacement parts for maintenance to ensure that the device works reliably.
- Watch out for water condensation inside and outside the device. Warm up the device first, if it has been brought in from a cold environment.
- Never fasten the vacuum pump over a heating bath.
- Ensure that no solids and/or liquids can enter the diaphragm of the vacuum pump IKA Vacstar digital through the suction line of the pump. This damages the diaphragm and other internal parts of the pump. This reduces delivery capacity and the final pressure can no longer be achieved. Deposits may be produced in the interior, reducing service life and causing leaks.


Unpacking

Unpacking

- Unpack the device carefully.
- Any damage should be notified immediately to the shipping agent (post office, railway network or transport company).

Scope of delivery

Diaphragm vacuum pump IKA Vacstar digital	
Power supply cable Laptop Euro	
Power supply cable Laptop GB	
Power supply cable Laptop CH	
USB A - USB B cable	
Vacuum hose 1m	

Silencer	
Connecting cable	
Operating instructions	
Warranty card	

Correct use

Use

WARNING

The device must not be used for:
 - **evacuating biological habitats,**
 - **evacuating explosive, corrosive or similar gases,**
 - **throughput/use of liquids**

The device is not suitable for:
 - **pumping up objects**
 - **pumping liquids**
 - **use in potentially explosive atmospheres**
 - **underwater**
 - **underground**

Standalone operation

The vacuum pump can be operated without a controller. The delivery capacity can be set using the adjustable speed. Generally speaking, the smaller the volume to be evacuated, the lower the speed should be. The final pressure is not controlled.

Operation with accessories

Along with the accessories recommended by **IKA**, the vacuum pump **IKA** Vacstar digital is suitable for controlled evacuation to a specified final pressure e.g. for operating an **IKA** rotary evaporator. It can also be used for classical separation, filtration or drying in the laboratory.

With an appropriate vacuum controller (e.g. VC 10) or with a rotary evaporator (e.g. RV 10 auto), the Vacstar digital can be operated automatically, for automatic boiling point detection, scheduled pressure-time curves or programs from the solvent library, for example.

Intended use: Tabletop device

The laboratory diaphragm vacuum pump **IKA** Vacstar digital is used together with other laboratory devices to create a vacuum.

For specific vacuum control, an additional vacuum controller **IKA** VC 10 or **IKA** rotary evaporator RV 10 auto is required.

Area of use (only indoors)

The device is suitable for use in residential areas and all other areas.

The safety of the user cannot be ensured:

- if the device is used in conjunctions with accessories not made or recommended by the manufacturer!

- if modifications are made to the device or parts of the device by third parties!
- if the device is operated in contravention of the safety instructions!

Useful information

DANGER

You can prevent the release of solvent vapours into the atmosphere using a downstream emission condenser.

CAUTION

Special care must be taken when working with gas mixtures that include condensable gases (e.g. solvents). If diaphragm pumps are used with such gases they condense in the pump's air flow. This leads to a build-up of pressure in the vacuum chamber, which in turn influences the effectiveness and service life of the diaphragms and valves.

If liquid or condensate enters the vacuum chambers, the specified minimum working pressure is no longer achieved.

To protect the inner valves and diaphragms against liquid, use an upstream condensate separator, such as a Woulff bottle.

CAUTION

It is generally recommended that you run the pump with no load at the end of the process (no-load operation, approx. 3 up to max. 5 minutes) to dispel any residues in the pump.

If applicable, vent the recipient. To do this, remove the suction line during manual operation.

NOTE

Observe the operating instructions of the relevant devices.

Diaphragm vacuum pumps are oscillating displacement pumps. The diaphragms expand the volume of the vacuum chamber to produce an underpressure, which sucks air into the chamber. The air is then pressed out of the vacuum chamber by subsequent compression. Changes between the inlet and outlet are controlled by diaphragm valves.

The physical unit for pressure is 1 Pascal [Pa].

However people often still use the unit millibar [mbar] colloquially.

1 mbar = 100 Pa
 1 bar = 10⁵ Pa
 1 Pa = 10⁻⁵ bar
 1 Pa = 0.01 mbar

A distinction is made between different ranges in vacuum technology.

Rough vacuum:	10 ⁵ to 10 ² Pa	(1000 to 1 mbar)
Fine vacuum:	10 ² to 10 ⁻¹ Pa	(1 to 10 ⁻³ mbar)
High vacuum:	10 ⁻¹ to 10 ⁻⁵ Pa	(10 ⁻³ to 10 ⁻⁷ mbar)
Ultra high vacuum	< 10 ⁻⁵ Pa	(< 10 ⁻⁷ mbar)

The periodic movements of the diaphragms create a vacuum at the pump inlet. Diaphragm vacuum pumps work with a rough vacuum of up to 10² Pa. Lower pressures cannot be achieved with these pumps.

To create fine, high and ultra high vacuums, oil-sealed vacuum pumps, screw pumps, diffusion pumps or turbomolecular pumps must be used, for example. Diaphragm vacuum pumps are often used as pumps for the preliminary stage in these contexts.

Compressing and releasing the pressure on air creates a noise. The pump is quieter the less air is transported. This means that there is less operating noise during normal operation if the pressure is reduced.

A fitted exhaust hose with a sound absorber can also be used to reduce noise.


In automatic operation with controlled speed the pump reduces its speed when it reaches the target value. If the leakage rate in the entire system is low enough, the speed may even be set to "zero".

When combined with the vacuum controller VC 10 or the rotary evaporator with integrated vacuum controller RV 10 auto, the pump can also be used for "two-position control" (valve-controlled) or "analogue speed control" of the vacuum.


Two-position control

When the target value is reached, the suction line is interrupted. Due to the small delay between the times when the pressure value is detected, when comparison is made with the target value, when a vacuum valve is switched off and when the suction line of the vacuum-producing system are switched off, it is easy to go below the set target value. If the pressure in the system goes up again due to a natural leakage rate, the suction line valve opens.

Schematic view of two-position vacuum control

Schematic view of speed vacuum control

Setting up

NOTE

Ensure that you follow the general rule, that the recipient (load/vacuum container/glass condenser) should always be connected to the suction line at its highest point. This reduces the risk of liquids entering the vacuum pump.

Liquid in the pump chambers makes the pump properties worse.

Connect the hose connections (INLET/OUTLET) and interfaces in accordance with the labelling on the pump and the operating instructions.

The connections are barbed hose connectors. Ensure that the hoses have hose clips if applicable. Lay the hoses in such a way that no condensate can flow into the pump.

Ensure that the air can enter the ventilation slit (H) unhindered.

Install a condensate separator (e.g. a Woulff bottle) in the suction line before the intake spigot of the pump to protect against water coming in.

If solvent vapours are sucked in, an emission condenser (accessory) after the pump helps to condense it, and prevents it being released into the atmosphere.

Connection of interfaces

NOTE

Observe the relevant connections (see Fig. 1).

O: Hose connection d= 8 mm OUTLET

Link this connection to the emission condenser with a vacuum hose, or fit a sound absorber to the end.

NOTE

Put the end of the hose in the extractor hood!

Check there is a free outlet on the pressurised side!

Do not use a throttle on the pressurised side and do not close the outlet! Connect the exhaust line to this connection.

K: Connection for vacuum controller VC 10 (Mini DIN)

You can connect the vacuum pump and the vacuum controller VC 10 or the rotary evaporator RV10 auto with the analogue connection cable for precise speed-controlled vacuum control. The vacuum controller detects the pump and switches to speed-vacuum control mode. Two-position control is deactivated.

L: USB interface

Connect the vacuum pump IKA Vacstar digital to a PC with a USB A - USB B cable. Any device software updates can be loaded using the IKA FUT software tool.

M: RS 232 interface

You can connect the IKA Vacstar digital vacuum pump to a PC using an RS 232 interface cable. The pump can be operated in conjunction with other devices with labworldsoft® laboratory device software. For further information see the chapter "Interfaces and outputs".

N: Water valve connection

Connect the optional water throttle valve RV 10.5001 to the diaphragm vacuum pump. The water flow to the emission condenser is controlled using the water valve. The valve is opened as soon as the pump is switched on.

I: Hose connection for suction line d= 8 mm INLET

Connect the suction line to this connection.

Link this connection to the recipient (rotary evaporator condenser, laboratory reactor etc.) with a vacuum hose.

P: Power supply cable connection

Check that the voltage information on the rating plate matches your mains supply. Connect the power supply cable to supply electricity.

IKA Vacstar digital with vacuum controller VC 10

NOTE

Cannot be operated with optional water valve RV 10.5001 in conjunction with VC 2.4.


Refer to the operating instructions for the vacuum controller VC 10!

When using the optional pump cut-out unit VC 2.4 (ref. 6) the pump is switched off when the target value is reached (operating mode B).

This device configuration is also possible without the pump cut-out unit VC 2.4. The pump runs on. For good control, the pump should be operated at a low speed.


When the target value is reached, the suction line is closed by the vacuum valve integrated in the VC 10. The pump speed can be changed by manually setting the speed.

- 1 Recipient (load, e.g. rotary evaporator, reactor)
- 2 Vacuum separator (e.g. Woulff bottle)
- 3 Vent valve
- 4 Pressure sensor
- 5 Vacuum valve/ball valve
- 6 Pump cut-out unit VC 2.4
- 7 Emission condenser

Speed regulation

IKA Vacstar digital with vacuum controller VC 10 or rotary evaporator RV 10 auto


NOTE

Observe the operating instructions for the devices.

Automatic setting of the operating mode "Speed regulation" as soon as the analogue connecting cable (ref. 8) connects the pump with the vacuum controller or rotary evaporator.

When the target value is reached, the pump speed is automatically reduced to a minimum of "0 rpm". Once the target value is reached, the pump applies suction according to the leakage rate of the vacuum system. Automatic boiling point recognition is possible.


- 1 Recipient (load, e.g. rotary evaporator, reactor)
- 2 Vacuum separator (e.g. Woulff bottle)
- 7 Emission condenser
- 8 Analogue connecting cable

Commissioning

NOTE Observe the operating instructions for the devices.

Check if the voltage specified on the type plate matches the available mains voltage. If these conditions are met, the device is ready for operation when it is plugged in to the mains.

 The socket used must be earthed (grounded contact).

Otherwise safe operation is not guaranteed or the device may be damaged. Observe the ambient conditions indicated under "Technical data".

Switching on

- Switch on the device using the main switch (D). All display segments are shown to enable a visual check.

8.8.8.8.

- Software version number display (e.g. vers. 1.00).

0.03

- Operating mode display (A; B). (Factory setting: Operating mode A)

A

The operating mode can be changed when you switch on the device.

Operating mode A:

In this operating mode, the set target value is not saved when the current run comes to an end or the device is switched off. No automatic restart after mains power is cut off.

Operating mode B:

In this operating mode, the set target value is saved when the current run comes to an end or the device is switched off, and the value can be changed.

Operating mode B is required particularly in combination with the pump cut-out unit VC 2.4.

Automatic restart after mains power is cut off.

Changing operating mode:

- Press the rotating/pressing knob (C) and switch on the device at the main switch (D) at the same time to change between the two operating modes.

A

b

In standalone operation or two-position operation in combination with the vacuum controller VC 10 the target speed can be set by turning the rotating/pressing knob (C). Push to start the pump and push again to stop it. A flashing light on the display (B) indicates when the device is running.

Switching on with several interfaces

If the pump is operated via the analogue connecting cable with the vacuum controller VC 10 or the rotary evaporator RV 10 auto, "VC 10" is shown on the display.


VC 10

If the vacuum controller is connected to a PC, "PC" is shown on the display (B).

PC

Switching off

Switch off the device using the main switch (D).

 Unplug the appliance from the mains socket to disconnect it from the mains.

Interfaces and outputs

NOTE Please comply with the system requirements together with the operating instructions and help section included with the software.

The device can be operated in "Remote" mode via an RS 232 or USB interface using the laboratory software labworldsoft®. The RS 232 interface (M) at the back of the device is fitted with a 9-pole SUB-D jack which can be connected to a PC. The pins are assigned serial signals.

USB Interface

The Universal Serial Bus (USB) is a serial bus system which allows the device to be connected to the PC. Devices that support USB can be connected to each other whilst they are running (hot plugging) and provide automatic recognition of the connected devices and their properties.

Use the USB interface in conjunction with labworldsoft® for operation in "Remote" mode and for updating the firmware using the "Firmware update tool".

Installation

Before the device is connected with the PC using the USB data cable, the USB driver must be installed.

The USB driver can be downloaded from the website:

Serial interface RS 232 (V 24)

Configuration:

- The functions of the interface circuit between the device and the automation system are a selection from the signals specified in the EIA standard RS232 as per DIN 66020 Part 1.
- Standard RS 232, corresponding to DIN 66259 Part 1 is valid for the electric characteristics of the interface circuits and assignment of signal states.
- Transmission process: Asynchronous character transmission in start-stop operation.
- Transmission type: Full duplex.
- Character format: Character composition according to data format in DIN 66022 for start-stop operation. 1 start bit, 7 character bits, 1 parity bit (even), 1 stop bit.
- Transmission speed: 9600 Bits/s.
- Data flow control: none
- Access method: Data transmission from the device to the computer only occurs after a request from the computer.

Command syntax and format

The following points should be noted for the instruction set:

- Commands are generally sent from the computer (master) to the device (slave).
- The device only responds to requests from the computer. Even error messages are not sent spontaneously from the device to the computer (automation system).
- The commands are transmitted in capital letters.
- Commands and parameters, as well as consecutive parameters, must be separated by at least one space (code: hex 0x20).
- Each individual command (including parameters and data) and all responses are completed with CRLF (code: hex 0x20 hex 0x0d hex 0x0a) and can have a maximum length of 50 characters.
- The dot is used for decimal separators in a floating-point value (code: hex 0x2E).

The details given above generally comply with the recommendations of NAMUR (NAMUR recommendations for the design of electrical plug-in connectors for analogue and digital signals in laboratory MSR devices. Rev. 1.1).

The NAMUR commands and the additional IKA-specific commands are only used as low-level commands for communication between the device and the PC. With an appropriate terminal or communication program, these commands can be transmitted directly to the device.

NAMUR commands Function

NAMUR commands	Function
IN_NAME	Read device name
IN_PV_4	Read current speed value
IN_SP_4	Read rated speed value
OUT_SP_4	Adjust the rated speed value
START_4	Start evacuation
STOP_4	Stop evacuation
RESET	Switch to normal operating mode
OUT_NAME	Set device name
OUT_SP_42@n	Setting WD safety limit speed with set value echo
OUT_WD1@m	Watchdog mode 1
OUT_WD2@m	Watchdog mode 2

"Watchdog" functions; monitoring of the serial data flow

If, once this function has been activated (see NAMUR commands), there is no retransmission of the command from the computer within the set time ("watchdog time"), the Evacuation function is switched off in accordance with the set "watchdog" mode or is returned to previously set target values. The data transmission may be interrupted by, for example, a crash in the operating system, a power failure in the PC or an issue with the connection cable between the computer and the device.

"Watchdog" mode 1

If event WD1 should occur, the evacuation function is switched off and ER 2 is displayed. Set watchdog time to m (20 - 1,500) seconds, with watchdog time echo. This command launches the watchdog function and must be transmitted within the set watchdog time.

"Watchdog" mode 2

If there is an interruption in data communications (longer than the set watchdog time), the speed target value is changed to the set WD safety speed limit. The warning PC 2 is displayed. The WD2 event can be reset with OUT_WD2@0 - this also stops the watchdog function.

Set watchdog time to m (20 - 1,500) seconds, with watchdog time echo. This command launches the watchdog function and must be transmitted within the set watchdog time.

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
Fig. 10

USB A - USB B cable

Required for connecting the USB interface (L) to a PC.


Fig. 11

Connecting cable

Required for connecting VC 10 and RV 10 auto with the diaphragm vacuum pump IKA Vacstar digital (speed operation).


Fig. 12

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Cleaning

For cleaning disconnect the mains plug!

Use only cleaning agents which have been approved by IKA to clean IKA devices.

These are water (with tenside) and isopropanol.

Wear protective gloves during cleaning the devices.

Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.

Do not allow moisture to get into the device when cleaning.

Before using another than the recommended method for cleaning or decontamination, the user must ascertain with IKA that this method does not destroy the device.

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- device type.
- serial number, see type plate.
- position number and description of spare part, see www.ika.com.

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For repair, please request the "Safety Declaration (Decontamination Certificate)" from IKA or use the download printout of it from IKA website at www.ika.com.

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If an error occurs, this is shown in an error on the display (B), e.g. Error 4.

Proceed as follows:

- Switch the device off on the device switch,
- Take corrective measure,
- Switch device back on.

Error	Cause	Reason	Solution
Er 3	Overheat	The pump does not run	Switch off the pump, disconnect it from the mains, let it cool down.
Er 4	Motor blocked/overload	The pump does not run	Switch off the pump, disconnect it from the mains Let the pump cool down. Ventilate the vacuum system to normal pressure.
Er 8	Hall sensor defective	The pump does not run	Service
Er 9	BLP self-test/EEPROM error	The pump does not run	Service
Er 10	Outlet Choked	The pump does not run	User needs to check the outlet
Er 11	Membrane stuck	The pump does not run	Air the pump and run at low speed for self-cleaning
Er 12	Connection problem	The pump does not run	Service: check the motor cable connection

In the event that the fault cannot be eliminated using the measures described or another error code is displayed:

- Please contact the Service Department,
- Send the device and a brief description of the fault.

Accessories

- Vacuum controller IKA VC 10
- Vacuum control valve IKA VCV 1 and VCV 2
- Check valve IKA VC 10.300
- Vacuum hose IKA VH.SI.8
- PC 1.1 cable (RS 232)
- Maintenance kit
- Vacuum safety emission condenser VSE 1
- Woulff bottle/water valve set VSW1

For further accessories see www.ika.de

Product contact parts

Name	Material
Connection spigot	PPS
Connection piece	PPS
Distributor	PPS
Head piece	PPS
Diaphragm valve	FFPM
Diaphragm	NBR/PTFE
Clamp	PPS

Technical data

	Base unit	Value
Delivery capacity max.	m³/h	1.32
Delivery capacity max.	L/min	22
Final pressure	mbar	2
Suction steps		4
Cylinder		4
Connection diameter suction side	mm	8
Connection diameter pressure side	mm	8
Gas ballast valve		no
Inlet pressure min.	mbar	2
Inlet pressure max.	mbar	1030
Boiling point recognition		no
Solvent library		no
Two-position control		yes
Analogue speed-vacuum control		yes
Speed setting option		rotating/pressing knob
Speed min.	rpm	285
Speed max.	rpm	1200
Display		LED
Vacuum sensor		no
Vent valve		no
Condensate separator		no
Emission condenser		no
Volume at minimum pressure	dB(A)	54
Product contact material		PTFE, FFPM, PPS, NBR
Housing material		Cast alu, coated
Appropriate for clean room		no
Explosion-proof		no
Width	mm	150
Depth	mm	370
Height	mm	375
Weight	kg	11.5
Permissible ambient conditions		5 °C to 31 °C at 80 % relative humidity. 32 °C to 40 °C decreasing linearly to a maximum 50 % relative humidity.
Protection class according to DIN EN 60529		IP 20
RS 232 interface		yes
USB interface		yes
Analogue output		no
Voltage	V	100-240
Frequency	Hz	50/60
Power consumption	W	130
Standby power consumption	W	1.5

	Base unit	Value
Protection class		I
Overvoltage category		II
Degree of contamination		2
Use above max. sea level	m	2000
Firmware update		yes
<i>Note</i>		
Protective earth conductor	The protective earth conductor (PE) is not connected to the metal housing. (Internally encapsulated DC power supply)	

*Average values.
Subject to technical changes!*

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EN

Contents

	Page		
Declaration of conformity	19	Commissioning	29
Explication of warning symbols	19	Assembling the stand pillar RV 8.1	30
Safety instructions	19	Maintenance and cleaning	31
Correct use	23	Accessories	31
Unpacking	23	Error Codes	31
Useful information	24	Technical data	32
Installation	24	Solvent table (excerpt)	33
Drive RV 8		Warranty	33
Glassware			
Hose system			
Heating bath			

Declaration of conformity

We declare under our sole responsibility that this product corresponds to the regulations 2014/35/EU, 2006/42/EC, 2014/30/EU and 2011/65/EU and conforms with the standards or standardized documents EN 61010-1, EN 61010-2-051, EN 61326-1, EN 60529 and EN 12100.

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Explication of warning symbols

General hazard.

This symbol identifies information **that is of vital importance for protecting your health and safety**. Disregarding this information may lead to health impairment and injuries.

This symbol identifies information **that is of vital importance for the technically correct functioning of the system**.

Disregarding this information may result in damage to the instrument or to system components.

This symbol indicates information **which is important for proper use and ensuring that the operations of the instrument are performed efficiently**.

Failure to observe this information may result in inaccurate results.

Danger - Reference to exposure to a hot surface!

Safety instructions

For your protection

The safety lift operates at power failure and when the device is switched off. It is designed to raise itself with the glassware attached. Therefore, please note that the lift without glass assemblies goes up quickly due to the reduced weight.

- Read the operating instructions completely before starting up and follow the safety instructions.
- Keep the operating instructions in a place where they can be accessed by everyone.

- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations. **When working under a vacuum in particular!**
- Wear your personal protective equipment in accordance with the hazard category of the medium to be processed. There may be a risk of:
 - splashing liquids,
 - body parts, hair, clothing and jewellery getting caught,
 - injury as a result of glass breakage.

Inhalation of or contact with media such as poisonous liquids, gases, spray mist, vapours, dusts or biological and microbiological materials can be hazardous to user.

- Set up the device in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.
- Ensure that there is sufficient space above the device as the glass assembly may exceed the height of the device.
- Prior to each use, always check the device, accessories and especially the glass parts for damage. Do not use damaged components.
- Ensure that the glass assembly is tension-free! Danger of cracking as a result of:
 - stress due to incorrect assembly,
 - external mechanical hazards,
 - local temperature peaks.
- Ensure that the stand does not start to move due to vibrations respectively unbalance.
- Beware of hazards due to:
 - flammable materials,
 - combustible media with a low boiling temperature,
 - glass breakage.

Only process and heat media that has a flash point higher than the adjusted safe temperature limit of the heating bath that has been set. The safe temperature limit of the heating bath must always be set to at least 25 °C lower than the fire point of the media used. (acc. to EN 61010-2-010)

- Do **not** operate the appliance in explosive atmospheres, with hazardous substances or under water.
- Only process media that will not react dangerously to the extra energy produced through processing. This also applies to any extra energy produced in other ways, e.g. through light irradiation.
- Tasks with the device must only be performed when operation is monitored.
- Operation with excess pressure is not permitted (for cooling water pressure see "Technical Data").
- Do not cover the ventilation slots of the device in order to ensure adequate cooling of the drive.
- There may be electrostatic discharges between the medium and the drive which could pose a direct danger.
- The appliance is not suitable for manual operation (except lift movement).
- Safe operation is only guaranteed with the accessories described in the "Accessories" chapter.
- Refer to the operating instructions for the heating bath.
- Refer to the operating instructions for the accessories, e.g. vacuum pump.
- Position the positive pressure outlet of the vacuum pump under a fume hood.
- Only use the device under an all side-closed exhaust, or a comparable protective device.
- Adapt the quantity and the type of distill and to the size of the distillation equipment. The condenser must work properly. Monitor the cooling water flow rate at the condenser outlet.
- The glass equipment must always be ventilated when working under normal pressure (e.g. open outlet at condenser) in order to prevent a pressure build-up.

- Please note that dangerous concentrations of gases, vapors or particulate matter can escape through the outlet at the condenser. Take appropriate action to avoid this risk, for example, downstream cold traps, gas wash bottles or an effective extraction system.
- Evacuated glass vessels must not be heated on one side; the evaporating flask must rotate during the heating phase.
- The glassware is designed for operation under a vacuum of down to 2 mbar. The equipment must be evacuated prior to heating (see chapter "Commissioning"). The equipment must only be aired again after cooling. When carrying out vacuum distillation, uncondensed vapours must be condensed out or safely dissipated. If there is a risk that the distillation residue could disintegrate in the presence of oxygen, only inert gas must be admitted for stress relief.

Avoid peroxide formation. Organic peroxides can accumulate in distillation and exhaust residues and explode while decomposing! Keep liquids that tend to form organic peroxides away from light, in particular from UV rays and check them prior to distillation and exhaust for the presence of peroxides. Any existing peroxides must be eliminated. Many organic compounds are prone to the formation of peroxides, e.g. dekalin, diethyl ether, dioxane, tetrahydrofuran, as well as unsaturated hydrocarbons, such as tetralin, diene, cumene and aldehydes, ketones and solutions of these substances.

The heating bath, tempering medium, evaporation flask and glass assembly can become hot during operation and remain so for a long time afterwards! Let the components cool off before continuing work with the device.

Avoid delayed boiling! Never heat the evaporating flask in the heating bath without switching on the rotary drive! Sudden foaming or exhaust gases indicate that flask content is beginning to decompose. Switch off heating immediately. Use the lifting mechanism to lift the evaporation flask out of the heating bath. Evacuate the danger zone and warn those in the surrounding area!

Safety lift

When the device is switched off or the power supply disconnected, the internal safety lift removes the evaporating flask from the heating bath.

The safety lift at loss of power is designed for a maximum total weight (glassware and solvent) of 3.1 kg.

Sample calculation of the maximum load with vertical glassware and a 1 litre flask:

Condenser + receiving flask + evaporating flask + fittings =
1200 gr + 400 gr + 280 gr + 100 gr = 1980 gr

Maximum loading of solvent = 3100 gr – 1980 gr = 1120 gr

Due to the design, operation of the safety lift cannot be guaranteed for higher loads!

When using other types of condensers such as dry ice or intensive condensers, and also when using return distillation distributors with slip-on condensers, it may be necessary to reduce the load by the amount of the added weight of the glass apparatus.

Thus, prior to distillation, check whether the lift, laden with the glassware and distillation material, rises when power is lost.

The safety lift must be checked daily prior to operation!

Attach the maximum total weight of 3.1 kg to the lift. Manually move the lift to the lowest position and press the "Power" key on the front panel or the main switch on the back right side of the device.

The evaporating flask is lifted out of the heating bath.

If the safety lift does not work, please contact the **IKA** Service department.

For the evaporating equipment (evaporating flask plus contents), the maximum permissible weight is 3.0 kg! Loads greater than this risk breakage of glass at the steam pipe!

Ensure that the safety lift has been powered down.

When working with large loads, always use low speeds.

Unbalanced loads can result in breakage of the steam pipe!

- A vacuum may be formed inside the glassware in the case of power outage. The glassware must be vented manually.


Never operate the device when the evaporation flask is rotating and the lift is raised. Always lower the evaporation flask into the heating bath first before starting the rotation drive. Otherwise hot tempering medium may be sprayed out!

- Set the speed of the drive so no tempering medium is sprayed out as a result of the evaporation flask rotating in the heating bath. If necessary reduce the speed.
- Do not touch rotating parts during operation.
- Imbalance may result in uncontrolled resonance behavior of the device or assembly. Glass apparatus may be damaged or destroyed. In the event of unbalance or unusual noises, switch off the appliance immediately or reduce the speed.
- The appliance does not start up again automatically following a cut in the power supply.
- The device is only disconnected from the power supply network if the device power switch is off or the plug is pulled out.
- The socket for the mains cord must be easily accessible.


For protection of the equipment

- The voltage stated on the type plate must correspond to the mains voltage.
- Socket must be earthed (protective ground contact).
- Removable parts must be refitted to the appliance to prevent the infiltration of foreign objects, liquids, etc.
- Protect the appliance and accessories from bumps and impacts.
- The appliance may only be opened by experts.

DANGER!

- ☞ Take care! Fast lift movement without glassware because of spring inside!
- ☞ Shut down device only in upper position!

Note:

- ☞ Safe lift movement!
- ☞ Device can be switched off at any time!

Intended use

• Use

Together with the accessories recommended by IKA, the device is suitable for:

- quick and gentle distillation of liquids,
- evaporation of solutions and suspensions,
- crystallization, synthesis or cleaning of fine chemicals,
- drying of powder and granulate material,
- recycling of solvents.

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- if the device is operated with accessories that are not supplied or recommended by the manufacturer,
- if the device is operated improperly or contrary to the manufacturer's specifications,
- if the device or the printed circuit board are modified by third parties.

• Further remarks

Legal bases and country-specific characteristics, especially for medicine and food:

The operator is responsible for compliance with the laws and regulations applicable to him.

Unpacking

• Unpacking

- Please unpack the device carefully
- In the case of any damage a report must be sent immediately (post, rail or forwarder)

• Contents of package

	Drive RV 8	Heating bath HB digital	Vertical glassware RV 10.1	Vertical glassware RV 10.10 coated	Serrated washer M6	Cylindrical M6x25	Ring spanner	Desktop switching RV 10.8001 Seal	Foot	Right-angle T20	HEXAGON DIN911 SW5	Holding bracket	Condensate flask	Screw Plug D	Screwed -Joint Cap	Tube clip D=9mm	Operating instructions	Base	Vapour tube	Power cable
Package 1 / RV 8 V	1	1	1		4	4	1	1	1	1	1	1	1	1	2	2	1	1	1	2
Package 2 / RV 8 VC	1	1		1	4	4	1	1	1	1	1	1	1	1	2	2	1	1	1	2
Package 3 / RV 8 flex	1	1			4	4	1	1	1	1	1	1	1	1	2	2	1	1	1	2

Contents of package RV 8


Fig. 6

Useful information

Distillation is a thermal separating process for liquid compounds based on substance-specific, pressure-dependent boiling points through evaporation and subsequent condensation.

The boiling point temperature decreases with decreasing external pressure which means that work is usually done under reduced pressure. In this way the heating bath can be maintained at a constant temperature (e.g. 60 °C). Using the vacuum, the boiling point is set with a steam temperature of approx. 40 °C. The cooling water for the condensation condenser should not be warmer than 20 °C (60-40-20 rule).

A chemical resistant membrane pump with a vacuum controller should be used to create the vacuum. The pump is protected from solvent residue by the addition of a Woulff bottle and/or a vacuum separator.

Working with a jet pump to create a vacuum can only be recommended to a limited extent as the solvents may contaminate the environment when using these systems.

Speed, temperature, flask size and system pressure all affect the evaporator air capacity.

The optimum capacity of the flow-through condenser is approx. 60%.

This corresponds to condensation on approx. 2/3 of the cooling coil. With larger capacities there is the risk that the uncondensed solvent vapor will be extracted.

The device is equipped with a Lift safety unit.

If the power cuts out, the evaporation flask is automatically lifted out of the heating bath by an integrated gas spring.

“CAUTION! The safety lift must be checked daily before use. See safety notes / safety lift ! “

The glass apparatus can be evacuated after power outage!

When using other types of condensers such as dry ice or intensive condensers as well as when using return distillation distributors with slip-on condensers, it may be necessary to reduce the mass of the media to compensate for the mass of this additional glass-ware. Thus, prior to distillation, check whether the lift goes up without power when laden with the glass and distillation material.

The device is designed for operation with a cooling water supply system (e.g. laboratory thermostat), but can also be run off a water supply line. Please refer to the Technical Data for information on cooling water pressure, temperature constancy, and flow rate.

Installation

Drive RV 8

1.) Mount the base to the lift. (Fig. 7)


Fig. 7

- Apply Cylindrical M6x25 (F) (4x) and Serrated washer M6 (G) (4x) to mount base (H) with Lift Unit (I) after you unpack the package. (Fig. 7)

2.) Remove transportation lock (Fig. 8)


Fig. 8

- Hold the lift with your hand on the height position and remove the thumb screw (J) by turning counterclockwise on the back of the appliance.

Take care! DANGER!

Once the transportation lock has been removed, the lift moves fast to its upper end position. The distance is approx. 120 mm.

3.) Fix the Desktop switching. (Fig. 9)


Fig. 9

- Attach the switching adapter (K) on the rear side. Secure it with the Velcro, by contracting both tabs on the switching adapter. The Velcro must not be below the switching adapter.
- Connect the Low-voltage connector (M) to the connection socket (N). (Fig. 9)
- Plug the power cord (L) in the switching adapter and connect it to the power source.
- Observe the valid supply voltage.

4.) Fix the Holding bracket (Fig.10)

Fig.10

5.) Insert the bottle and attach the supplied hose connectors to the bottle. (Fig. 11)

Fig.11

6.) Adjust the angle of the rotation head. (Fig. 12)

- Remove the clamping device for the angle setting of the rotation drive on the right side of the lift by rotating the knob counterclockwise.

Fig.12

7.) Set the drive at an angle of approx. 30°. (Fig. 13)

Fig.13

- Then, secure the rotary drive against rotation by tightening the adjustable knob clockwise.

Glassware

NOTE: Please read the operating instructions of the glassware for the safe handling of laboratory glassware !

- There is a locking knob (P) for locking or unlocking the locking device (O).

If a red mark (Q) is shown, the locking device is unlocked, otherwise it is locked.

To lock or unlock the locking device (O), push the locking knob (P) to the end position.


Fig. 14a

- Open the locking device (O) on the drive head by turning it 60° counterclockwise. Then you can see a red mark (Q) .


Fig. 14b

- Feed the vapour tube (R) in until it stops.
- Then, lock the locking device (O) by turning it clockwise by 60°.


Fig. 14c

- Pushing the lock knob (P) to the end position, the red mark (Q) should be covered and invisible.
- The vapour tube is not allowed to be pulled out!
- Check the correct axial locking device on the vapour tube.
- Keep the red mark (Q) invisible.

Fit & unfit the evaporation flask / push off mechanism


Fig. 15a


Fig. 15b

Turn clockwise (↻) for tight-fitting
Turn counterclockwise (↺) for loosen tight-fitting

Mount the evaporation flask:

- Place the evaporation flask on the Vapour tube (R). (Fig. 15a)
- Turn the plastic nut (S) with the clip (T) counterclockwise until the clip can be put on the flange. (Fig. 15b)
- Now turn the plastic nut clockwise until the clamp is hard up against the flange. When doing so, hold the locking device (O) on the steam pipe.

Remove the evaporation flask:

- Hold the locking device (O) and loosen the plastic nut by turning it counterclockwise. This releases a tightly clamped evaporating flask.
- Hold the evaporating flask at the flange and open the jaws of the clamp.
- Remove the evaporating flask.
- Check that the locking of vapour tube is still closed!

- Place the RV 10.8001 condenser seal (U) in the condenser receptacle and fit the glassware to the device according to the assembly instructions. (Fig. 16, 17 and 18)

Fig. 16

Mount the condenser

- Place the cap nut into the condenser, and then put the spring ring to the flange of the condenser (V).
- Place the condenser on the rotary drive (W) and tighten the cap nut by hand. (Fig. 17)


Fig. 17

- Fix the receiving flask and hose connectors as shown. See also operating manual of glassware. (Fig. 18)

Fig. 18

Mounting the washer


Fig. 19a

Drainage washer

Fig. 19b

Note: Pay attention to the correct position of the washer.

Removing the condenser

- Use the ring spanner provided to loosen union nuts that are tightly fitted.
- Loosen the union nut by turning anticlockwise.

Description of special condensers

(see products/accessories at www.ika.com)

• RV 10.3 Vertical-intensive condenser with manifold

This vertical-intensive condenser features a double jacket design for particularly efficient condensation. Also available with coating (RV 10.30)

• RV 10.4 Dry ice condenser

Dry ice condenser for distilling low-boiling solvents. Cooling by dry ice, no cooling water required. Maximum condensation thanks to low temperatures.

Also available with coating (RV 10.40)

- **RV 10.5 Vertical-condenser with manifold and cut-off valve for reflux distillation**

Also available with coating (RV 10.50)

- **RV 10.6 Vertical-intensive condenser with manifold and cut-off valve for reflux distillation**

This vertical-intensive action condenser features a double jacket design for particularly efficient reflux distillation.

Also available with coating (RV 10.60)

Hose system

- Connect the water hoses (**H₂O**) to the condenser according to the counter-flow principle. (Fig. 20)
- Install the vacuum connections to the condenser, Wouff bottle, vacuum controller with valve and vacuum pump.
- Always connect the vacuum hose (**vac**) to the condenser at the highest point to minimize solvent losses during suctioning.
- Use standard laboratory vacuum hoses with an inside diameter of 8 mm and a wall thickness of 5 mm (see accessories).


Fig. 20

Heating bath

Refer to the chapter entitled "Operation" in the operating instructions of heating bath!

- **Assemble the foot of the heating bath**

For operation of the heating bath in combination with the RV 8, the foot plate supplied must be fitted to the heating bath.

Caution: It is not permitted to use heating bath fitted with a foot plate in combination with any other type of rotary evaporator.


Fig. 21

- Place the heating bath on the stand of the rotation drive and push it into the left position. (Fig. 22)

Fig. 22

Commissioning

The unit is ready for service when the mains plug has been plugged in.

Fig.5

Filling the evaporation flask:

You can fill the evaporation flask manually prior to creating the vacuum. The evaporation flask should not be filled more than half its volume.

The maximum allowed load (evaporation flask and contents) is 3 kg.

Filling medium by vacuum: prior to filling the evaporation flask, a vacuum controller is used to regulate the glass apparatus to the target pressure.

- Now fill the evaporation flask using the backfeed line.
- Due to the vacuum present, the solvent is suctioned into the evaporation flask. This enables you to keep solvent loss due to suctioning to a minimum.

Setting up the heating bath:

Also refer to the operating instructions for the heating bath!

- Move the lift to the bottom position and check the position of the heating bath in relation to the evaporation flask. When using larger evaporation flasks (2 or 3 liters) or depending on the angle of the rotation drive, you can move the heating bath to the right.
- Fill the heating bath with the tempering medium until the evaporation flask is surrounded by tempering medium to 2/3 of its volume.
- Switch on the rotation drive and slowly increase the speed.

Note: Avoid creating waves.

- Switch on the heating bath using the main on/off switch.
- Note: Avoid stress on the glass due to different evaporation flask and heating bath temperatures when lowering the evaporation flask into the heating bath!*

Function keys (Fig.3)

Fig.3

- **"Power" key (1)**

When the power switch (5) is switched on, the device will be in "activated" mode.

Touch the "Power" key (1) for about 1 second, the status of the device will be changed from "activated" mode to "stand-by" mode.

Touch the "Power" key again shortly, the status of the device will be changed back to "activated" mode immediately.

- **Display window/speed (2)**

Display (2) shows "8.8.8."

then Softwareversion "X.X".

Display shows "set value".

Note: The device is not deenergized until the mains power supply has been interrupted, for instance by unplugging the power cable.

- **Rotating knob (3)**

To start the rotation, press the rotary knob.


The indicator flashes until the set speed in rpm (revolution per minute) is reached.

By turning the rotary knob, the speed can be changed.

To stop the rotation, press the rotary knob again.

The last displayed speed value is stored.

Lift position key (4&6) (Fig. 23a and Fig. 23b)

To adjust the lift position, grasp the handle and press the left or right lock button (4).
The hand must rest on the sensor (6).
Note: Hand detection sensor (6), only to be touched by hand, do not press!
The lift lock is unlocked and you can move the lift up and down. Release the lock button or remove the hand from the lift, the lift is locked in the adjusted position.
When locking and unlocking, a clear "click" sound is heard.

Assembling the stand pillar RV 8.1 (Accessories)

No stand is required if the condenser has been properly fitted and secured by tightening the condenser union nut on the rotary drive.
The purpose of the stand is solely to prevent the condenser twisting.

Caution: If the RV 8.1 stand is fitted incorrectly, the strong forces at the Velcro strap can cause stresses within the glass, which in turn can cause damage to the glass condenser.

After the condenser has been fitted, the RV 8.1 stand is attached to the rotary drive. Make sure that the condenser is installed parallel to the lift body.

- Guide the stand rod through the hole on the underside of the rotary drive,
- Screw the knurled screw into the tapped hole but do not tighten it,
- Move the stand until the upper rubber sleeve rests against the glass condenser,
- Tighten the knurled screw finger-tight,
- Secure the condenser to the stand using the pre-fitted Velcro strap.

Fig.24

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Cleaning

For cleaning disconnect the mains plug!
Use only cleaning agents which have been approved by IKA to clean IKA devices.

These are water (with tenside) and isopropanol.
Wear protective gloves during cleaning the devices.
Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.
Do not allow moisture to get into the device when cleaning.
Before using another than the recommended method for cleaning or decontamination, the user must ascertain with IKA that this method does not destroy the device.

Spare parts order

When ordering spare parts, please give:
- device type.
- serial number, see type plate.
- position number and description of spare part, see www.ika.com.
- software version.

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For accessories see www.ika.com.

Error Codes

Any malfunctions during operation will be identified by an error message on the display.
Once a serious error message has been displayed, the lift moves to the top end position and the device can no longer be operated.

Proceed as follows in such cases:

- switch off device using the main switch,
- carry out corrective measures,
- restart device.

Error code	Cause	Effect	Correction
E03	Temperature inside the device is too high	Room temperature > 40 °C	- Switch off the device - Let the device cool down - Check the location (see Technical data)
E04	Target speed is not reached Motor blocked	Load too large Cable break	- Reduce the volume in the evaporator
E09	Flash memory storage error	Transmission error	- Turn the unit off and on again

If the actions described fail to resolve the fault or another error code is displayed, then take one of the following steps:

- contact the service department,
- send the device for repair, including a short description of the fault.

Technical data

Operating voltage range	Vac	(100 – 240) ± 10%
Rated voltage	Vac	100 – 240
Frequency	Hz	50 / 60
Rated power of RV 8 device (without heating bath)	W	75
Rated power (max.) of switching adaptor	W	90
Speed	rpm	5 – 300
Speed tolerance	Set speed: < 100 rpm Set speed: ≥ 100 rpm	rpm % ± 1 ± 1
Speed display	mm	digital
Dimensions of visible display area (W x H)		37 x 18
Display		7 segments display
Smooth start		yes
Lift		manual
Stroke	mm	120
Head angle adjustable		0° – 45°
Cooling surface	cm ²	1500
Cooling water flow rate min.	l/h	30
Cooling water flow rate max.	l/h	100
Cooling water pressure max.	bar	1
Perm. On-time	%	100
Perm. ambient temperature	°C	5 – 40
Perm. relative humidity	%	80
Protection acc. to DIN EN 60529		IP 20
Protection class		I
Contamination level		2
Weight (no glassware; no heating bath)	kg	12.5
Dimensions (W x D x H)	mm	595 x 390 x 615
Operation at a terrestrial altitude	m	max. 2000

Subject to technical changes!

Solvent table (excerpt)

Solvent	Formula	Pressure for boiling point 40 °C in mbar (For HB approx. 60 °C)
Acetic acid	C ₂ H ₄ O ₂	44
Acetone	C ₃ H ₆ O	556
Acetonitrile	C ₂ H ₃ N	226
N-Amylalcohol	C ₅ H ₁₂ O	11
n-Pentanol	C ₅ H ₁₀ O	11
n-Butanol	C ₄ H ₁₀ O	25
tert. Butanol	C ₄ H ₁₀ O	130
2-Methyl-2-Propanol	C ₄ H ₁₀ O	130
Butylacetate	C ₆ H ₁₂ O ₂	39
Chlorobenzene	C ₆ H ₅ Cl	36
Chloroform	CHCl ₃	474
Cyclohexane	C ₆ H ₁₂	235
Dichloromethane	CH ₂ Cl ₂	atm. press.
Methylenechloride	CH ₂ Cl ₂	atm. Press.
Diethylether	C ₄ H ₁₀ O	atm. press.
1,2,-Dichloroethylene (trans)	C ₂ H ₂ Cl ₂	751
Diisopropylether	C ₆ H ₁₄ O	375
Dioxane	C ₄ H ₈ O ₂	107
Dimethylformamide (DMF)	C ₃ H ₇ NO	11
Ethanol	C ₂ H ₆ O	175
Ethylacetate	C ₄ H ₈ O ₂	240
Ethylmethylketone	C ₄ H ₈ O	243
Heptane	C ₇ H ₁₆	120
Hexane	C ₆ H ₁₄	335
Isopropylalcohol	C ₃ H ₈ O	137
Isoamylalcohol	C ₅ H ₁₂ O	14
3-Methyl-1-Butanol	C ₅ H ₁₂ O	14
Methanol	CH ₄ O	337
Pentane	C ₅ H ₁₂	atm. press.
n-Propylalcohol	C ₃ H ₈ O	67
Pentachloroethane	C ₂ HCl ₅	13
1, 1, 2, 2, -Tetrachloroethane	C ₂ H ₂ Cl ₄	35
1, 1, 1, -Trichloroethane	C ₂ H ₃ Cl ₃	300
Tetrachloroethylene	C ₂ Cl ₄	53
Tetrachloromethane	CCl ₄	271
Tetrahydrofurane (THF)	C ₄ H ₈ O	357
Toluene	C ₇ H ₈	77
Trichloroethylene	C ₂ HCl ₃	183
Water	H ₂ O	72
Xylene	C ₈ H ₁₀	25

Warranty

In accordance with **IKA** warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine directly to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

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