

S25D- ... / T25 basic, DI25 basic


S18D- ... / T18 basic, DI18 basic

# IKA

designed for scientists

**Kunststoff - Dispergierwerkzeuge**  
**Plastic dispersion tools**  
**Outils dispersants**

S18/25-ET50

# IKA

designed for scientists

**IKA-Werke GmbH & Co. KG**  
 Janke & Kunkel-Straße 10, 79219 Staufen, Germany  
 Phone: +49 7633 831-0, Fax: +49 7633 831-98  
 eMail: sales@ika.de

**USA**  
**IKA Works, Inc.**  
 Phone: +1 910 452-7059  
 eMail: sales@ika.net

**KOREA**  
**IKA Korea Ltd.**  
 Phone: +82 2 2136 6800  
 eMail: sales-lab@ika.kr

**BRAZIL**  
**IKA Brazil**  
 Phone: +55 19 3772 9600  
 eMail: sales@ika.net.br

**MALAYSIA**  
**IKA Works (Asia) Sdn Bhd**  
 Phone: +60 3 6099-5666  
 eMail: sales.lab@ika.my

**CHINA**  
**IKA Works Guangzhou**  
 Phone: +86 20 8222 6771  
 eMail: info@ika.cn

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Source language: German	<span><span></span></span> <b>EN</b>
Safety instructions	
<ul style="list-style-type: none"><li>No liquid may get into the drive unit of the disperser.</li> <li>Abrasion of the dispersion equipment or the rotating accessories can get into the medium you are working on.</li> <li>While working with the disperser, the user must wear his personal protective equipment, in particular safety goggles.</li> <li>Follow the health and safety regulations accident prevention regulations.</li> <li>Do not use any damaged dispersion tools e.g. with hairline cracks or splits.</li> <li>Note the operating instructions of the disperser unit.</li> <li>When using disposable tube <b>S18/25-ET50</b>, ensure that the dispersing element has fully engaged in the cover.</li> <li>The maximum capacity is 30 ml.</li> <li>Only operate in a vertical working position.</li> <li>The distance between the stand and the dispersing instrument must not exceed 80<span> </span>mm.</li></ul>	
Product information	
<ul style="list-style-type: none"><li>The dispersion tool comprises shaft tube, rotor shaft and two bearing shells.</li> <li>Dispersion tools made of plastics are not sterile and not free of DNA or RNA.</li> <li>The tools are ideal for applications, with which no cross contamination is permissible.</li> <li>The tools are ideally appropriate for a single - use. However they can be used, depending upon the application, several times too.</li> <li>In this case please consider the cleaning references.</li> <li>Use disposable tube <b>S18/25-ET50</b> for working with the plastic dispersing element in a closed system.</li></ul>	
Application instruction	

Application instruction
<ul style="list-style-type: none"><li><i><b>Attention!</b></i> Never run the dispersion tool dry. The bearing will be destroyed without cooling by the medium.</li> <li>Always ensure that the dispersing element does not run dry in some cases should spouts develop.</li> <li>The fill level of the medium must be between the minimum and maximum marks.</li> <li>Medium may leak out of the overflow hole around the maximum mark.</li> <li>The dispersion tools are not suitable for continuous operation.</li> <li>Dispersing causes the medium to heat up.</li> <li>The dispersing element must not be operated at the maximum permitted speed of 25000 rpm for more than one minute.</li> <li>Store the dispersing elements in a dry place, in their packaging and at room temperature.</li> <li>The solid bodies in samples must not be frozen or too hard.</li> <li>Likewise the optimal dispersion duration and rotating frequency must be determined by attempts. Usually a few seconds are sufficient, maximum duration is one minute.</li> <li>Longer application times bring no improvements, increase however the sample temperature substantially.</li> <li>Dispose of the tools after use according to the regulations. Observe the laboratory regulations.</li></ul>
Materials and chemical resistance

- The shaft tubes are made of polycarbonate (PC), the rotor shafts and bearing shells are made of polyetheretherketon (PEEK).
- The plastic materials used conform to the FDA (Food and Drug Administration Approval) - rules.
- The parts have a good stability against weak acids, chlorides, hypochlorides and many other chemicals.

Cleaning, disinfection and sterilization		
The dispersion tools are not sterile, but they can be sterilised before use.		
	Allowed procedures	Not allowed procedures
Method of sterilisation	Autoclave up to 125 <span> </span> °C with 2 bar <p>Ethylene oxides Formaldehyde (g) Gluteralhyde (2%) Gamma radiation up to 5 MRAD</p>	UV - irradiation <p>Hot-air sterilisation 160-190<span> </span>°C</p>
Biological disinfection	Alcohol (70%) <p>Hypochlorides (5%) Formaldehyde (40%)</p>	Phenol derivatives <p>Expanded use of guanidine</p>
Radioactive decontamination	Radiacwash <p>Count - Off Water / ethyl alcohol / SDS</p>	

Technical data					
		<b>S18D–10G–KS</b>	<b>S18D–14G–KS</b>	<b>S25D–10G–KS</b>	<b>S25D–14G–KS</b>
Suitable for disperser		T 18 digital	T 18 digital	T 25 digital	T 25 digital
Working range	<b>ml</b>	10–100	10–500	10–100	10–500
Stator / Rotor (Ø)	<b>mm</b>	10 / 6.8	14 / 9.8	10 / 6.75	14 / 9.5
Max. allowable speed	<b>rpm</b>	25000			
Max. circumferential speed	<b>m/s</b>	8.9	12.8	8.5	12
Min. / max. immersion depth	<b>mm</b>	15 / 85			
Materials in contact with medium		Polycarbonat (PC), Polyetheretherketon (PEEK)			
Max. working temperature	<b>°C</b>	100			
Sterilization methods		autoclavable			

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