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Warranty

In accordance with IKA warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover wearing parts, nor does it apply to faults resulting from improper use or insufficient care and maintenance contrary to the instructions in this operating manual.

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
Assembling of the dispersion tool:

- Open the knurled screw (pos. 10) so that the thread does not protrude into the hole.
- Insert the dispersing element (pos. 9) from the bottom into the drive unit as far as the stop. Ensure that the dispersion tool is aligned correctly when inserted (refer to picture: the tabs must project sideways when viewed from the front of the drive unit/dispersion tool). The shaft will audibly click into place when a small amount of force is applied. The dispersion tool is correctly mounted when the edge of the projecting part is aligned with the lower edge of the drive unit (see fig. 2).
- Now secure the dispersing element by screwing the knurled screw tight (pos. 10).

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
IKA Works Guangzhou
Phone: +86 20 8222 6771
eMail: info@ika.cn

INDIA

IKA India Private Limited
Phone: +91 80 26253 900
eMail: info@ika.in

BRAZIL

IKA Brazil
Phone: +55 19 3772 9600
eMail: sales@ika.net.br

POLAND

IKA Poland Sp. z o.o.
Phone: +48 22 201 99 79
eMail: sales.poland@ika.com

UNITED KINGDOM

IKA England LTD.
Phone: +44 1865 986 162
eMail: sales.england@ika.com

Discover and order the fascinating products of IKA online:

www.ika.com

IKAworldwide

IKAworldwide /// #lookattheblue

@IKAworldwide

Technical specifications may be changed without prior notice.

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

Safety instructions

- The dispersion tool may heats up during operation. In particular in the processing of hot medium, the tool must be cooled before it can be removed from the drive.
- When in operation the dispersion tool can become hot. If the dispersion tool is not inserted into the drive flange correctly or the medium is hot, it can become extremely hot.

- No liquid may get into the drive unit of the disperser.
- While working with the disperser, the user must wear his personal protective equipment, in particular safety goggles.
- Follow the health and safety regulations accident prevention regulations.
- The work has to be interrupted immediately if you notice unusual noise and/ or emission of liquid from the upper rinsing drilling. The seals must be controlled, and replaced if necessary.
- Note the operating instructions of the disperser unit.
- Abrasion of the dispersion equipment or the rotating accessories can get into the medium you are working on.
- The gasket and bearings are made of PTFE and rustproof steel; the following points should therefore be noted: *Chemical reactions of PTFE occur in contact with molten or solute alkali metals and alkaline earth metals, as well as with fine powders of metals in groups 2 and 3 of the periodic system at temperatures above 300 °C - 400 °C. Only elementary fluorine, chlorotrifluoride and alkali metals attack it; halogenated hydrocarbons have a reversible swelling effect.* (Source: Römpps Chemie-Lexikon and "Ulmann", Volume 19)

Product information

Rotor-laterally the shaft bearing of the dispersion tool is a sleeve bearing bushing (material PTFE). Output - laterally the shaft of the rotor is inserted in the coupling of the drive unit.

The plastic materials used are approved by FDA (Food and Drug Administration Approval).


Application instructions

- Attention!** Never run the dispersion tool dry. Without cooling by the medium the bearing will be destroyed. Always ensure that the dispersing element does not run dry in some cases should spouts develop.
- Dispersing causes the medium to heat up.
- The fill level of the medium must be at least 10 mm above the lower overflow hole of the shaft.
- While working with the disperser, the user must wear his personal protective equipment.
- Longer application times bring hardly any improvements, increase however the sample temperature substantially.
- The function of the dispersion tools depends on the condition of the sharp edges of rotor and stator. Abrasive media can round off these edges fast, whereby the disperse effect diminishes.

Spare parts list

Item	Designation	Item	Designation	Spare parts diagram see front page.
1	Shaft tube	1	Shaft tube	For orders of spare parts, please specify the type of disper-sion tool and the designation of the spare part.
2	Shaft	2	Grooved piston ring bearing	
3	Rotor	3	Shaft	
4	Grooved piston ring bearing	4	Stator	
		5	Rotar	
		6	Driveshaft ring	

Maintenance and cleaning

Allowed procedures	Sterilisation techniques	Flat key Shaft key
Damp heat	Autoclave to 121°C at 2 bar positive pressure	
Chemical procedures	By germ - killing solvents alcohol, phenol,formalin…., disinfectant remainders are to be removed with germ - free water	
Hot air	Germ - killing by hot-air with 160 to 190°C (approx. 30 min)	Fig. 1

For cleaning the dispersion tool, operate it in a dissolvent for dissolving the residuaries of the substance. Due to the large flow rate, rotor and stator are cleaned to a large extent. The dispersion tool must be divided and cleaned immediately after working, so that the adhering residuaries of the substance don't cultivate unwanted bacterial cultures.

For cleaning purposes the dispersion tools can be disassembled by the flat key and the shaft key, which are parts of the tool sets of the drive, how described below:

- S 18 N-10 G:**

Unscrew the rotor (3) from the shaft (2) by the flat key (locking of the rotor (3)) and the shaft key (turning of the shaft (2) clutch-laterally). Pull the shaft (2) out upwards of the shaft tube (1). Now press the grooved piston ring bearing (4) out upwards of the shaft tube (1).

- S 18 N-19 G:**

Unscrew the stator (4) from the shaft tube (1) by the flat key (***Attention!*** Left-handed thread) and pull it out of the shaft tube (1), completely with shaft (1) and rotor (5). Unscrew the rotor (5) from the shaft (3): Put the shaft key on the end of the shaft and turn the shaft key, while locking the rotor (5) with the flat key. Now pull the grooved piston ring bearing (2) out of the stator.

The assembly of the dispersion tool takes place in reverse order. ***Attention!*** The rotor must not be tightened too strongly, otherwise the thread / rotor teeth will be damaged.

Technical Data

		S 18 N - 10 G	S 18 N - 19 G
Max. allowable speed	rpm	25000	
Working range	ml	1–100	10–1500
Stator / Rotor (Ø)	mm	10 / 7.5	19 / 12.7
Max. circumferential speed	m/s	9.8	16.6
Min. / max. immersion depth	mm	25 / 70	35 / 170
Material in contact with medium		PTFE / AISI 316L	
Suitable for solvents		Yes	
Max. working temperature	°C	180	
Sterilization methods		all methods	
Ultimate fineness, suspensions	µm	10–50	
Ultimate fineness, emulsions	µm	1–10	

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