20000007206

IKA®

C 1.20_052016

IKA® Heizung C 1.20

Betriebsanleitung	DE	3
<i>Ursprungssprache</i>		
Operating instructions	EN	11
Mode d'emploi	FR	19
Instrucciones de manejo	ES	27
Instruções de serviço	PT	35
Руководство по эксплуатации	RU	43
操作说明	ZH	51
取扱説明書	JA	59
사용 설명서	KO	67

IKA®-Werke, Germany
Reg. No. 004343

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
Fig. 1

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Contents

	Page
Device setup	2
Declaration of conformity	11
Explication of warning symbols	11
Safety instructions.....	12
Correct use	12
Unpacking	13
System components.....	13
Commissioning	14
Error codes	16
Accessories	17
Maintenance and cleaning.....	17
Warranty	17
Technical Data	18

Declaration of conformity

We declare under our sole responsibility that this product is in compliance with the regulations 2014/35/EU, 2014/30/EU and 2011/65/EU and conforms to the standards or standardized documents: EN 61010-1, EN 61010-2-010, EN 61326-1, EN 60529 and EN ISO 12100.

Explication of warning symbols

Indicates an (extremely) hazardous situation, which, if not avoided, will result in death, serious injury.

Indicates a potentially hazardous situation, which, if not avoided, can result in death, serious injury.

Indicates a potentially hazardous situation, which, if not avoided, can result in injury.

Indicates practices which, if not avoided, can result in equipment damage.

Safety instructions

Read the operating instructions in full before starting up and follow the safety instructions.

- Keep the operating instructions in a place where they can be accessed by everyone.
 - Ensure that only trained staff work with the appliance.
 - Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
 - Set up the device in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.
 - Do not operate the appliance in explosive atmospheres, with hazardous substances.
 - Protect the appliance and accessories from bumps and impacts.
 - Check the appliance and accessories beforehand for damage each time you use them. Do not use damaged components.
 - Safe operation is only guaranteed with the accessories described in the "Accessories" chapter.
 - The socket for the mains cord must be easily accessible.
 - The plug to be used must be earthed (grounded contact).
 - The voltage stated on the type plate must correspond to the mains voltage.
 - The appliance can only be disconnected from the mains supply by pulling out the mains plug or the connector plug.
 - Disconnect the power cord before attaching or replacing accessories.
 - Disconnect the power cord before cleaning, maintaining, or transporting the device.
 - The device must only be opened by trained, skilled personnel, even when repairs are to be made. Unplug the device before opening. Current-carrying parts in the interior of the device may remain live for some time after the device has been unplugged.
- The C 1 calorimeter's power supply unit must never be positioned beneath the heater or beneath water hoses.
 - Connect the C 1.20 heater to the C 1 in accordance with these operating instructions.
 - Use only drinking-quality tap water for operation.
 - The permissible operating pressure must not exceed the value given in the technical data.
 - Check the hose system between the devices for tightness before commissioning.
 - Shut off the water supply before disassembling the hose system.
 - The device may only be used for its intended purpose and as described in these operating instructions. This applies to operation by trained and qualified staff as well.
 - Transport the device carefully.

WARNING

Users must make sure that the on/off switch of the IKA® device can be accessed immediately, directly and without risk at any time. If installation or positioning cannot ensure this access at all times, an additional EMERGENCY OFF switch that can be easily accessed must be installed in the work area.

- Ensure that the hoses and tubes do not slip or bend.
- Regularly check hoses, tubes, and heater for possible material fatigue (cracks/leakage).

DANGER

Do not start up the device if:

- it is damaged or leaking,
- cables (not just the power cable) are damaged.

Correct use

• Use:

The C 1.20 heater may only be operated in conjunction with the C 1 calorimeter. The C 1.20 heater serves to maintain the correct temperature of drinking-quality tap water, which must be connected to the heater at a maximum pressure of 6 bar/0.6 MPa. The C 1 control connection cable supplied is used to control the C 1 calorimeter.

The conditions for operating the C 1 calorimeter with the C 1.20 heater on a water supply line are:

- tap water: drinking-quality tap water is recommended,
- temperature range: see "Technical data" section; the measuring mode must be selected according to the water temperature:
water temperature < 21 °C : measuring mode 22 °C,
water temperature > 21 °C : measuring mode 30 °C.
- water pressure: see "Technical data" section (for high or irregular pressure, we recommend using the C 25 water pressure control valve, available in our accessories range).

Operating mode: tabletop device

• Range of use (interior areas only):

- Laboratories
- Schools
- Universities

Use the C 1.20 heater in closed spaces only.

The unit is suitable for use in all areas apart from domestic areas and areas directly connected to a low-voltage supply which also serves domestic areas.

The safety of the user cannot be guaranteed if the appliance is operated with accessories that are not supplied or recommended by the manufacturer or if the appliance is operated improperly contrary to the manufacturer's specifications or if the appliance or the printed circuit board are modified by third parties.

Unpacking

- **Unpacking:**

- Please unpack the device carefully,
- In the case of any damage a detailed report must be sent immediately (post, rail or forwarder).

- **Delivery scope Heater C 1.20**

- Heater **C 1.20**
- C 1 control connecting cable, 8-pin
- Power supply cable
- Hose connector set
- "Handle" tool for disconnecting the inflow pipe
- Operating instructions
- Warranty card

System components

Heater C 1.20

Fig. 2

C 1 control connecting cable, 8-pin.


Fig. 3

Power supply cable (not shown)

"Handle" tool for disconnecting the inflow pipe

Fig. 4

Hose connector set



Fig. 5

C 1 inflow pipe

(supplied with the C 1 calorimeter)

Commissioning

• Place of installation:

Set up the device in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.

The installation area must be large and well-ventilated enough to ensure that the heat given off by the device does not overheat the room. There should be 10-15 cm of clearance on the front and back sides.

• Installation:

☞ Remove the C 1.20 heater from the packaging.

Fig. 7

☞ Remove the water inlet hose on the C 1.

☞ Set the C 1.20 heater behind the C 1 on a level, stable surface.

Fig. 8

• Connecting the water hose

☞ Connect the C 1 water inlet hose to the C 1.20 heater as follows (IN).

Preparing the C 1 inflow pipe for the C 1.20 heater


Fig. 9

- Separate the C 1 water inlet hose using the supplied "Handle" tool into part **A** (IN C 1.20) and Part **B** (OUT C 1.20).


Fig. 10

- Attach hose connector **C** from the supplied hose connector set to the filter from Part **A**.


Fig. 11

- Attach hose connector **D** from the supplied hose connector set to Part **B**.


Fig. 12

- This results in


Fig. 13

Connecting parts X/Y C 1 to the C 1.20 heater (see Fig. 14)

- Connect Part **X** to the C 1.20 heater's input.
- Connect Part **Y** to the C 1.20 heater's output.

Connecting the water supply

NOTICE Shut off the water supply before disassembling the hose system.

- Connect Part **X** with your water supply (note the maximum pressure in accordance with the "Technical Data" section).
- If you are connecting a water tap: Ensure that the C 1 return and overflow hoses are positioned over a drain (see Fig. 14).
- Before startup, check the hoses between the devices for leaks. Gradually turn on the water tap. After the test, turn off the water tap.


Fig. 14

• Connection to C 1 calorimeter and mains power cable

- ☞ Connect the heater to the C 1 calorimeter with the supplied C 1 control cable.




Fig. 15

• Starting the C 1.20 heater

Verify that the voltage on the rating plate corresponds to the available power supply.

   The plug to be used must be earthed (grounded contact).

The power supply unit must never be positioned beneath the heater or beneath water hoses.

If these conditions are met the unit is ready for service when the mains plug has been plugged in.

Otherwise, safe operation cannot be ensured, and the device may suffer damage.

Also pay attention to the ambient conditions listed in the chapter "Technical Data".

Check the hose system between the devices for tightness before commissioning.

- ☞ Turn on the heater using the on/off switch.
- The switch lights up green.


Fig. 16

• C 1 calorimeter settings

- ☞ Turn on the C 1 calorimeter and commission it according to the operating instructions.

NOTICE Oxygen must be connected and the Pressure regulating valver must be open.

- ☞ On the C 1 calorimeter, select "Heater on" in the "Accessories" menu and exit the menu by saving.

ACCESSORIES	
Heater:	on
Combustible crucible:	no
Printer:	off
Service info:	no
Balance:	off
CANCEL SAVE	

Fig. 17

☞ Turn the C 1 calorimeter off and on again; at this point, the heater will be recognized by the C 1 calorimeter.

NOTICE The error message "Heater not found" should not be displayed (see "Error codes" section).

☞ Open the water tap slowly.

CAUTION Check connections for tightness.

☞ Check the cooling water temperature on the C 1 calorimeter in the "Accessories" menu under the "Condenser test" item.

☞ Then choose the corresponding measuring mode for the C 1 calorimeter in the "settings" menu (see "Use for the intended purposes" section).

NOTICE If the water inlet temperature is too high, select a higher temperature range for the C 1 to allow a measurement to be performed.

NOTICE If the water inlet temperature is low, you should always select the lower operating temperature for the C 1 to save time and energy.

☞ Perform a system test in accordance with the C 1 calorimeter operating instructions or a measurement.


Fig. 18


Fig. 19

Error Codes

Proceed as follows in the case of error code:

- ☞ Switch the device off,
- ☞ Perform corrective measures as described below,
- ☞ Perform a restart.

Error	Cause	Effect	Solution
Heater not found	C 1.20 heater is not recognized by the C 1 calorimeter.	No heater	Always turn the heater on first. Check the heating connection cable plug contact. Switch the C 1 calorimeter off and on again.


Fig. 20

If the actions described fail to resolve the fault or another error code is displayed then take one of the following steps:

- Contact the service department,
- Send the device for repair, including a short description of the fault..

Accessories


C 25 Pressure regulating valve

Fig. 21

Maintenance and cleaning

Cleaning:

For cleaning disconnect the main plug.

Only clean **IKA®** appliances using these **IKA®** approved cleaning agents:

Dirt	Cleaning agent
Dyes	isopropanol
Building materials	water containing detergent/ isopropanol
Cosmetics	water containing detergent/ isopropanol
Food	water containing detergent
Fuels	water containing detergent

For other materials, please consult our Technical Application Laboratory.

Wear protective gloves when cleaning the devices.

Do not place electrical appliances into the cleaning agents for cleaning purposes.

Do not allow moisture to get into the appliance when cleaning.

Before using a cleaning or decontamination method other than that recommended by the manufacturer, check with the manufacturer that the intended method will not destroy the device.

Ordering spare parts

When ordering spare parts, please give:

- Device type
- Manufacturing number, see type plate
- Item number and designation of the spare part, see www.ika.com.
- Software version

Repair

Please send instrument in for repair only after it has been cleaned and is free from any materials which may constitute a health hazard.

For this you should request the **“Decontamination Clearance Certificate”** from **IKA®**, or use the download printout of it from the **IKA®** website www.ika.com.

Return the instrument in its original packaging. Storage packaging is not sufficient. Also, please use suitable shipping package materials.

Warranty

In accordance with **IKA®** warranty conditions, the warranty period is 12 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual.

Technical Data

Rated voltage	VAC	230 ± 10 %
Rated voltage	VAC	100 - 115 ± 10 %
Frequence	Hz	50/60
Max. input power	W	1700
Medium		drinking-quality tap water
Max. permissible medium operating pressure	bar	6
Min. medium temperature	°C	5
Max. medium temperature	°C	25
Min. flow rate at 1 bar (0.1 MPa)	l/h	40
Max. flow rate at 1 bar (0.1 MPa)	l/h	80
Recommended flow rate	l/h	60
Water usage per measurement	l	4
Fuse		2 x T6,25 A 250 V 2 x T12 A 115 V
Perm. duration of operation	%	100
Protection acc. DIN EN 60529		IP 20
Perm. ambient temperature	°C	+ 5 ... + 40
Perm. rel. humidity	%	80
Dimensions (W x D x H)	mm	265 x 115 x 250
Weight	kg	6
Operation at a terrestrial altitude	m	max. 2000 above sea level

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IKA®-Werke GmbH & Co.KG