IKA

designed for scientists

IKA C 248

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Contents

EN

	Page
EU Declaration of Conformity	7
Warranty	7
Warning symbols	7
Safety instructions	8
Correct use and description of the device	8
Unpacking	8
Installation	9
Commissioning	9
Maintenance and cleaning	9
Accessories	10
Technical data	10
List of spare parts-Illustration of spar parts	2

EU Declaration of conformity

We declare under our sole responsibility that this product corresponds to the directives 2014/68/EU and 2011/65/EU and conforms with the following standards or normative documents: EN ISO 12100.

A copy of the complete EU Declaration of Conformity or further declarations of conformity can be requested at sales@ika.com.

Warranty

In accordance with IKA Terms and Conditions of Sale, the warranty period is 12 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual.

Warning symbols

DANGER

Indicates an (extremely) hazardous situation, which, if not avoided, will result in death, serious injury.


WARNING

Indicates a potentially hazardous situation, which, if not avoided, can result in death, serious injury.

CAUTION

Indicates a potentially hazardous situation, which, if not avoided, can result in injury.

NOTICE

Indicates practices which, if not avoided, can result in equipment damage.

Safety instructions

For your protection

Read the operating instructions in full before starting up and follow the safety instructions.

- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.

Work with the device

Before use, check the outlet pressure on the pressure reducer and compare it with the pressure displayed on the pressure gauge.

Oxygen as a compressed gas is oxidising; intensively aids combustion; can react violently with flammable materials.

- The filling time should not be set to less than 30 seconds. This will prevent combustible substances from being whirled around in the crucible.
- Insufficient ventilation can cause formation of explosive mixtures. Therefore, only use the device in well-ventilated areas.
- Check the seals for damage before each use.
- After finishing work, the oxygen supply on the pressure reducer must be closed.

For protection of the device

Maintenance tasks should be performed without any excess pressure on the system.

- Keep lines and fittings through which oxygen flows free of grease.

Correct use and description of the device

Use

The **IKA C 248** oxygen filling station is only suitable for filling decomposition vessels **C 5010** with oxygen.

Description of the device

- The 2 m pressure hose, included with delivery, can easily be laid to connect to the oxygen bottle.
- Therefore, the location of the **C 248** filling station does not depend on the oxygen bottle being in the immediate vicinity.
- The filling pressure that is set and the pressure in the decomposition vessel can be read directly from the manometer of the **C 248** (11), thus preventing insufficient filling quantities because of not observing the minimum filling period.
- Constant opening and reasoning on the oxygen bottle are void because the automatic non-return valve (bulb compl. (2)) only opens if a decomposition vessel is being filled.

Unpacking

Unpacking

- Please unpack the device carefully.
- In the case of any damage a detailed report must be sent immediately (post, rail or forwarder).

Scope of delivery

- **C 248** oxygen filling station
- 3 x O-rings
- Screw fitting
- Gasket
- Operating instruction


Installation

Connecting the pressure hose

The pressure hose, included with delivery, is designed for a maximum pressure of 40 bar. It is 2 m long and easily be laid to connect to the oxygen bottle or another oxygen inlet with pressure reducing valve (smallest permissible bending radius 80 mm).

There is an R1/4inch threading (SW 17) on the **C 29** reducing valve (accessory) of the oxygen bottle to which the other end of the pressure hose can be connected.

An appropriate adapter is included for US pressure reducers.

Decomposition vessels

The **IKA** decomposition vessel **C 5010** is provided with an identification on the base or on the cap screw that contains information including the year of manufacturing number.

Commissioning

The required filling pressure -usually 30 bar- can be set on the pressure reducer valve of the oxygen bottle. The adjusted pressure can also be read on the pressure gauge (11).

The bulb compl. (2) integrated into the oxygen filling station prevents oxygen from flowing out when this is not intended.

When the task is complete, the main valve on the oxygen supply should be closed.

The decomposition vessel should be placed on the marked position.

All you need to do then to fill up the decomposition vessel is to swing the rocking arm downwards. The increase in pressure (ca. 2 bar) can be observed on the pressure gauge.

The filling time should not be less than 30 seconds.

Maintenance and cleaning

The device is maintenance-free. It is only subject to the natural wear and tear of components and their statistical failure rate.

Cleaning

- For cleaning disconnect the mains plug!
Use only cleaning agents which have been approved by IKA to clean the devices:
Water containing surfactant / isopropyl alcohol.
- Wear protective gloves during cleaning the devices.
- Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.
- Do not allow moisture to get into the device when cleaning.
- Before using another than the recommended method for cleaning or decontamination, the user must ascertain with IKA that this method does not destroy the device.

It should be expected that frequent use will result in wear on the O rings in the bulb compl. (2). The filling head (3) is screwed in place and can be unscrewed after the pressure hose is removed. The lower gasket (3007), which is subject to particularly high wear and tear, can be reached without disassembling the pressure hose. Pull with a suitable tool (tweezers) the gasket (3007) down. Put the new gasket into the recess of the decomposition vessel. Place the decomposition vessel under the filling head (3) and fill you the decomposition vessel. The gasket remains at the bulb compl. (2) back.

Spare parts order

When ordering spare parts, please give:

- device type.
- serial number, see type plate.
- position number and description of spare part, see www.ika.de.

Repair

Please only send devices in for repair that have been cleaned and are free of materials which might present health hazards.

For repair, please request the "**Safety Declaration (Decontamination Certificate)**" from IKA or use the download printout of it from IKA website at www.ika.com.

If your appliance requires repair, return it in its original packaging. Storage packaging is not sufficient when sending the device - also use appropriate transport packaging.

Accessories

See more accessories on www.ika.com.

Technical data

Permissible operating pressure	bar	max. 40
Pressure gauge	bar	0-50
Typical filling pressure		30
Pressure hose	m	2
Filling speed	sec	> 30
Current filling pressure		readable
Oxygen valve		automatic non-return valve
Dimensions (B x T x H)	mm	130 x 180 x 330
Weight	kg	6

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IKA India Private Limited
IKA England LTD.
IKA Vietnam Company Limited
IKA Works (Thailand) Co. Ltd.