Sicherheits-Laborgasbrenner

Safety Laboratory Gas Burner

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*with the enclosed description fulfils the requirements of the Low Voltage Directive 2006/95/EC (valid to 19. April 2016) the Low Voltage Directive 2014/35/EU (valid from 20. April 2016), the regulations of the EMC Directive 2014/30/EC and the eco design Directive 2009/125/EC.*

*Hereby, we certify that our products, regardless of the production location, RoHS compliant and fulfill the directive 2011/65/EC (revised version: directive 2002/95/EC).*

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
*A. Wegener*



Firmenstempel / Company stamp

Armin Wegener  
 Vice President Research & Development


- 1 - Function knob
- 2, 2a - Dual knob: 2 - Gas adjustment  
2a - Air adjustment
- 3 - LED Standard
- 4 - LED Start-Stop
- 5 - LED Sensor
- 6 - IR-Sensor
- 7 - LED-Button
- 8 - LED Burner head HOT / BHC
- 9 - Burner head
- 9a - Flame orifice

- 10 - Connector for foot pedal
- 11 - Power connector for 9 V DC
- 12 - Gas inlet R 1/4" L
- 13 - Holding device for inoculation loop holder

- 14 - Burner head screw
- 15 - Monitor electrode
- 16 - Ignition electrode

- 17 - Active nozzle
- 18 - Nozzle holder for alternative gas
- 19 - Cover of the burner shaft
- 20 - Retaining screw for cover of the burner shaft
- 20a - Position screw for the cover
- 21 - Guide slots for tilt adjustment

## Table of contents

<b>Use</b> .....	22
<b>Safety Precautions</b> .....	22
<b>The range</b> .....	24
<b>1. Set up</b> .....	25
1.1 Choice of nozzle.....	25
1.2 Gas connection.....	25
1.3 Electrical connection.....	26
1.4 Foot pedal connection.....	26
<b>2. Operation</b> .....	26
2.1 On-Off switch, operating mode.....	26
2.2 Flame regulation.....	26
2.3 Switch-off the unit.....	26
<b>3. Application programs and operations</b> .....	27
3.1 Button Start-Stop.....	27
3.2 Pedal Standard.....	27
3.3 Pedal Start-Stop.....	27
3.4 Sensor Start-Stop.....	27
3.5 Sensor Auto-Off.....	28
3.5.1 Learn function Sensor Auto-Off burning time.....	28
3.6 Adjusting the IR-Sensor detection range.....	28
<b>4. Safety symbols and safety functions</b> .....	29
4.1 Residual heat display.....	29
4.2 Burner Head Control (BHC).....	29
4.3 Automatic unit switch off.....	29
4.4 DoubleClick IR-Sensor.....	29
<b>5. Error displays</b> .....	30
5.1 Ignition failure.....	30
5.2 Flame failure.....	30
5.3 Overtemperature.....	30
5.4 Burner head assembly monitor.....	30
5.5 Burner Head Control (BHC).....	30
<b>6. Cleaning and sterilizing</b> .....	31
6.1 Burner head disassembly and cleaning.....	31
6.2 Burner shaft cleaning.....	31
<b>7. Turbo flame</b> .....	32
<b>8. Tilt adjustment</b> .....	32
<b>9. Warranty</b> .....	32
<b>Technical data</b> .....	33
<b>Troubleshooting guide</b> .....	34
<b>Declaration of Conformity</b> .....	36

Read these instructions carefully to familiarize yourself with the product. Please retain these operating instruction for future reference.

## Use

**Safety laboratory gas burners for heating and flame sterilizing.  
Ideal for use in cleanroom workbenches and laboratory.**

**All users who have been assigned to use this device must have read and understood these operating instructions or have been instructed by an expert user so that this device can be used safely without causing danger.**

**ATTENTION: OPERATE THIS LABORATORY GAS BURNER UNDER CONSTANT SUPERVISION ONLY!**

## Safety Precautions

- On unpacking the unit, check for possible transportation damages. Do not operate the unit if damages are visible.
- After use or for any longer period of time without attendance, turn the main gas supply off and turn off the gas burner at the function knob (1).
- All gas connections must be adequately tightened (left-hand thread). Ensure gas proofness with suitable test equipment. **DO NOT** seal up the thread of the gas connection (12) of the laboratory gas burner with Teflon tape, etc.
- **BEFORE** using the device carefully check the gas feed tube for leaks. Check this even if the device has been installed by your distributor. To do this, carry out all the procedures mentioned in these operating instructions (see paragraph 1.).
- In the event that gas can be smelled: immediately turn off the gas supply to the device. Extinguish any open flames. Pull out the mains plug. Check all gas connections for gas proofness. If the smell of gas persists, the appropriate authorities must be notified (janitor, gas utility company, Fire Brigade).  
**LEAKING GAS CAN CAUSE A FIRE OR AN EXPLOSION. THIS MAY RESULT IN SEVERE INJURIES, FATAL ACCIDENTS AND DAMAGE TO PROPERTY.**
- The device can be dangerous if operated or used in an incorrect manner by untrained staff.
- An incorrect gas connection may create a hazard. Observe the installation instructions in the manual.


- Do not store spare or unconnected gas cartridges / gas bottles in the vicinity of this device.
- Even in an apparently empty gas cartridge / gas bottle, some gas may still remain. Gas cartridges / gas bottles should be transported and stored accordingly. Empty gas cartridges should be properly disposed.
- Do not use the device if there is a smell of gas or if there is a leak.
- NEVER try to loosen or unfasten gas connections while the gas supply is turned on and the device is in operation.
- Pay attention to your relevant rules for using liquid gas.
- Only use DVGW safety tubings with thread or tubing connectors. Check the condition of the tube / hose frequently. Depending upon type of tube / hose, hose clamps are required.
- NEVER use an open flame to look for leaks.
- Do not smoke if you are searching for leaks.
- Keep hands or other parts of the body away from the burner orifice (9a).
- Do not operate the unit near flammable liquids or hazardous materials.
- Unattended operation of the unit is not permissible.
- Always work in a well-ventilated area.
- Note that the burner orifice (9, 9a) remains hot after the flame has been extinguished. Do not touch. Can cause burns.
- Allow sufficient time for flame orifice (9a) to cool down prior to cleaning, disinfecting, servicing or transport. Ensure that the unit and the gas supply are turned off.
- Because of the connectors at the back of the unit the backside should not be sterilized with a flame.
- Allow sufficient time for burner head (9) to cool down prior to disassembling.
- Operate the unit with assembled burner head (9) only.
- After cleaning the burner head (9) allow sufficient time to dry before assembling again.
- Do not allow anything to fall into the flame orifice (9a).
- Before mounting a nozzle check the O-Ring (26). Replace the O-Ring if damaged or worn.
- When working with this device, always wear protective glasses.

## The range

### gasprofi 1 micro SCS

**Art.-No. 6.004.000**

- with DoubleClick IR-Sensor and stainless steel foot pedal
- adjustable IR-Sensor reaction distance
- 5 standard programs for button (function knob), IR-Sensor and foot pedal
- SCS (Safety Control System), BHC (Burner Head Control)
- Removable and decomposable burner head
- Tilt mechanism, right / left **(24)**
- Holding device for 3 inoculation loop holders **(28)**
- Nozzles for natural gas, propane/butane gas
- Turbo flame
- Wrench 17 mm **(22)** for gas connection
- Screwdriver **(23)** for burner head and cover of the burner shaft
- Tubing connector with swivel nut for 10 mm pipe spout tubing **(25)**
- Power connection with 4 adapter **(27)**, global
- Instruction manual
- 2-year warranty

## 1. Set up

Initial operation is to be carefully carried out as described in the following paragraphs.

Failure to observe the instruction manual **may give rise to hazards from leaks and / or bursts of flame.**



### 1.1 Choice of nozzle

**The unit is shipped with the nozzle for natural gas (17) (N, Ø 0.8 mm) installed. The nozzle for propane / butane (P, Ø 0.5 mm) can be found in the nozzle holder (18).**

The nozzle on the device must be fitted or the gas type in use, in which case proceed as follows: On the bottom of the device, remove the nozzle **P** for propane / butane by unscrewing it from the nozzle holder (18) with a coin (1 Cent coin) or the screwdriver (22). Then remove the active nozzle **N** for natural gas (17) (preinstalled) and exchange the nozzles. Tighten both nozzles again after changing over. The optional nozzle **B** (Ø 0.45 mm) is only for usage with the adapter CV 360.

Each time before mounting a nozzle, check the O-ring (26) for damage. **Replace O-ring if damaged, Art. No.: 8.000.010**

EN

### 1.2 Gas connection


Connect the device via the gas inlet (12) to the gas supply, if the operating pressure for natural gas is 18 - 25 mbar, for propane / butane 47.5 - 57.5 mbar. Only use DVGW or other gas approved safety tubing with thread or tubing connectors.

Check the condition of the tube / hose regularly.

Depending on the type of hose, hose clips and/or the hose clamps and the swivel nut included (25) are to be used.

A DVGW tested or other gas approved pressure regulator (50mbar) must be used for liquid gas. Pay attention to your relevant rules for using liquid gas.

If using gas cartridges, use only an original WLD-TEC gas cartridge adapter! Observe the operating instructions in question for the gas cartridge adapter.

**DO NOT seal up the thread of the gas inlet (12) of the laboratory gas burner with Teflon tape, etc.**

All gas connections must be adequately tightened with a wrench SW17 **(22)** (included in delivery). Note: left-hand thread. Check gas proofness with testing equipment or soapy water. To do this with the gas supply turned on, wet the gas connections with test solution or soapy water and look for bubbles. If bubbles form or if one bubble becomes larger, there is a leak. Immediately turn off the gas supply. Check all gas supply connections or gas joints (e.g.: at the gas cartridge / gas bottle, at the laboratory gas burner or at an additional gas leak protection). In the event of leaks, tighten up the gas joints with a wrench and test again for gas proofness.

**Never use an open flame to look for leaks.**

### 1.3 Electrical connection

Insert the power cord into the socket **(11)** on the back panel of the unit, or into the socket of the foot pedal. The mains power supply unit included is designed for a voltage from 100 - 240 V, 50/60 Hz. The mains power supply unit may only be connected if these values comply with the intended electricity supply. Before using it, push the corresponding socket adapter onto the mains power supply plug.

### 1.4 Foot pedal connection

Insert the connection cable of the foot pedal into the socket **(10)** at the back of the unit.

## 2. Operation

This section describes the operation of the unit with the function knob and the gas / air adjustment.

### 2.1 On-Off switch, operating mode

Switch the unit on by a short push on the function knob **(1)**. It can be turned off by a long push ( 2 seconds +) on the function knob. By turning the function knob all programs can be selected and the IR sensor detection range can be adjusted.

### 2.2 Flame regulation

The flame can be varied in size and intensity by turning the gas knob **(2)** and adjusting the air knob **(2a)** to suit all requirements.

When operating the unit for the first time or after changing the nozzle, turn the gas adjustment knob **(2)** two revolutions to the left and turn the air adjustment knob **(2a)** 3-4 revolutions to the left, too.

### 2.3 Switch-off the unit

The unit can be turned off by pushing the function knob **(1)** for more than 2 seconds.

**To depressurize gas hose, shut off the gas supply and activate the burner again until the flame extinguishes in order to burn the residual gas. Afterwards turn off the gas burner at the function knob (1).**

### 3. Application programs and operations

This chapter explains the various application programs and explains how to adjust the burn times and how to set the IR-sensor range.

#### 3.1 BUTTON StartStop

(LED Button is on)

The flame is ignited by operation of the function knob **(1)**. The flame is extinguished after renewed actuation of the function knob **(1)**. Additionally the flame is automatically extinguished when the burning timer has expired after 60 min.


#### 3.2 PEDAL Standard

(LED Standard is on)

The flame is ignited by operation of the foot pedal. The foot pedal remains depressed for the duration of use. The flame is extinguished upon release of the pedal.

#### 3.3 PEDAL Start-Stop

(LED StartStop is on)

The flame is ignited by operation of the foot pedal. The flame is extinguished after renewed actuation of the foot pedal. Additionally the flame is automatically extinguished when the burning timer has expired after 60 min. Alternatively the flame can be extinguished by a short push on the function knob **(1)**.

#### 3.4 SENSOR Start-Stop

(LED StartStop and LED Sensor is on)

For flame ignition, activate the DoubleClick IR-Sensor **(6)** by passing your hand over it twice. **(see paragraph 4.4)**. The flame is extinguished after renewed activation of the IR-Sensor. Additionally the flame is automatically extinguished when the burning timer has expired after 60 min. Alternatively the flame can be extinguished by a short push on the function knob **(1)**.

### 3.5 SENSOR Auto-Off

(LED Standard and LED Sensor are on)

For flame ignition, activate the DoubleClick IR-Sensor **(6)** by passing your hand over it twice. **(see paragraph 4.4)** The flame is extinguished automatically when the adjusted burning time expires (factory setting: 5 sec). Alternatively, the flame can be extinguished by a short push on the function knob **(1)**.

The burning time can be adjusted with the learn function “SENSOR Auto-Off burning time“ **(see paragraph 3.5.1 learn function “SENSOR Auto-Off burning time“)**. To activate the learn function, press the function knob **(1)** while the application program SENSOR Auto-Off is selected.

#### 3.5.1 Learn function “SENSOR Auto-Off burning time“

(LED Standard is on, LED Sensor blinks slowly)

This learn function is used to adjust the burning time for the application program SENSOR Auto-Off. To activate the learn function, press the function knob **(1)** while the application program SENSOR Auto-Off is selected.

Learning the burning time: For flame ignition and to start the time measurement, activate the DoubleClick IR-Sensor **(6)** by passing your hand over it twice.

When the needed burning time is reached (1 sec - 60 min) press the function knob **(1)** to stop the flame. Thereby, the reached burning time is memorized and the application program SENSOR Auto-Off is selected again. Now, the memorized time is available for the application program SENSOR Auto-Off.

#### 3.6 Adjusting the IR-Sensor detection range

(LED Sensor blinks fast)

Turn the function knob **(1)** clockwise till the LED Sensor **(5)** blinks fast to select the function for adjusting the IR-Sensor detection range. To modify the IR-Sensor detection range, press the function knob **(1)** briefly. (IR-Sensor coverage, 35 - 70 mm) Now, the IR-Sensor range can be adjusted by turning the function knob **(1)**. The frequency of the blinking LED Sensor **(5)** indicates the adjusted detection range.

Clockwise rotation:

The IR-Sensor detection range increases and the blinking of the LED Sensor **(5)** becomes faster.

Counter-clockwise rotation:

The IR-Sensor detection range decreases and the blinking of the LED Sensor **(5)** becomes slower.

Hold the hand in front of the IR-Sensor **(6)** to test the adjusted IR-Sensor detection range. When the IR-Sensor is activated by the hand, the LED StartStop **(4)** lights up. If the required detection range is not yet reached, modify the IR-Sensor detection range again by turning the function knob **(1)** clockwise or counter-clockwise.

When the required IR-Sensor detection range is reached press the function knob **(1)**. The adjusted IR-Sensor detection range is now saved and the application program returns to SENSOR StartStop.

#### 4. Safety symbols and safety functions

The safety LED's appear during usage to warn the user of potential hazards.

##### 4.1 Residual heat display

The LED „Burner Head Hot / BHC“  **(8)** indicates a hot burner head. As long as the LED shines, DO NOT TOUCH the burner head. This can cause serious burns. Even after switching-off the unit, the residual heat LED remains "on" till the burner head is cooled down.

**Disconnecting the power supply or removing the power cord will clear the residual heat display even if the burner head is still hot.**

##### 4.2 Burner Head Control (BHC)

If the burner head is clogged, the amber LED "Burner Head HOT / BHC" **(8)** will flash you have to clean the burner head (**see paragraph 5.5**), otherwise the maximum burning time in the operating mode "Button" and Pedal "Start-Stop" is limited to 30 seconds (**see paragraph 6.1**). If burning times longer than 30 seconds are required in case of a clogged burner head, the operating mode "Pedal Standard" can be used without time limit.

##### 4.3 Automatic unit switch off

The unit switches itself off automatically after 4 hours if it is not used during this time. The button & IR-Sensor programs switch off after 10 minutes. Even in the event of a fault, the device shuts off after 4 hours and the error display goes out (**see paragraph 5**). For further operation, switch the unit on again.

##### 4.4 DoubleClick IR-Sensor

For flame ignition, activate the IR-Sensor **(6)** by passing your hand over it twice within a time range of 1 sec. After the first activation the LED Sensor blinks fast for 1 sec. Within this time range the IR-Sensor **(6)** needs to be activated a second time to ignite the flame.

## 5. Error displays

If there is an error during operation, eventually the gas valve of the burner automatically closed and different flashing LEDs indicate the error.

All error displays can be reset by a long push (> 2 seconds) on the function knob (1). In case of overtemperature the unit needs to be cooled down and in case of burner head assembly monitor, the burner head needs to be reinstalled before a reset is possible.

### 5.1 Ignition failure

#### **Green LED "Button", "Standard" or "Start-Stop" blinks 2x**

This signal appears and indicates a malfunction if the flame fails to ignite after 7 seconds. In case of ignition failure check the burner head (9) for possible clogging, check the correct input pressure of the gas supply and verify that the correct nozzle is installed. In case of this malfunction the gas valve of the burner will be shut off automatically.

Nozzle N: natural gas, 18 - 25 mbar

Nozzle P, B: propane / butane gas, 47.5 - 57.5 mbar

### 5.2 Flame failure

#### **Green LED "Button", "Standard" or "Start-Stop" blinks 3x**

This signal indicates a malfunction if the flame is extinguished by external factors and fails to reignite within 5 sec. In case of flame failure check the burner head (9) for possible clogging and verify the correct input pressure of the gas supply. In case of this malfunction the gas valve of the burner will be shut off automatically.

### 5.3 Overtemperature

#### **Green LED "Button", "Standard" or "Start-Stop" blinks 4x**

This signal indicates a malfunction if the interior temperature has exceeded 70 °C. At a normal room temperature with normal air circulation the unit is suited for continuous operation. In case of overtemperature increase the air ventilation or change the operation site. In case of this malfunction the gas valve of the burner will be shut off automatically.

### 5.4 Burner head assembly monitor

#### **Green LED "Button", "Standard" or "Start-Stop" blinks 5x**

This message indicates that the burner head is removed. Further operation is possible after the burner head is reinstalled.

### 5.5 Burner Head Control (BHC)

#### **Orange LED "Burner Head HOT / BHC" (8) flashes**

This signal indicates that the time limit (30 sec) is turned on in operating mode "Start-Stop" and "Button" due to a clogged burner head. For cleaning the burner head **see paragraph 5.1.**

## 6. Cleaning and sterilizing

Allow sufficient time for burner orifice (**9, 9a**) to cool down before disassembling or cleaning the burner head. Check the unit is disconnected and that the gas supply is turned off at the mains. The burner can be cleaned with customary commercial disinfectants. Additionally, it is possible to remove the burner head and to clean it separately. The stainless steel and glass construction allow 100% UV-radiation sterilization and short time surface flame sterilization.

Because of the connectors at the back of the unit the backside should not be sterilized with a flame!


### 6.1 Burner head disassembly and cleaning

Allow sufficient time for burner orifice (**9,9a**) to cool down before disassembling or cleaning the burner head. Check the unit is turned off, that the gas supply is turned off at the mains. Clean the burner head with customary commercial disinfectants, sterilize it in an autoclave or wash it in a dishwasher. To remove the burner head proceed as follows:

Unscrew the burner head screw (**14**) completely with the included screwdriver. Turn approx. 8 revolutions to the left.

Now remove the burner head from the device by pulling it upwards. Reinstallation is performed in the reverse sequence. The dismantled burner head can be even dismantled into the individual components for in-depth cleaning: Unscrew both screws (**14a**) with the included screwdriver and take off the base plate (**14b**) of the burner head which was fixed by the two screws (**14a**). After the base plate is removed both electrodes can be pulled out for separate cleaning. Reinstallation is performed in the reverse sequence.

**i** When dismantling the burner head completely the sealing ring placed around the burner head screw (**14**) could dropout. Ensure that the sealing ring is placed around the burner head screw (**14**) when reassembling.


### 6.2 Burner shaft cleaning

Unscrew the screw (**20**) completely with the included screwdriver. Take off the cover (**19**) of the burner shaft.

Now the burner shaft can be cleaned or solid substances which have fallen into the unit can be removed. Reinstallation is performed in the reverse sequence.

Take care that the notch of the cover fits to the screw (**20a**). Screw in the screw (**20**) completely with the screwdriver.


## 7. Turbo flame

If the cover of the burner shaft (19) is removed the flame is extremely firm and consistent. To take off the cover of the burner shaft unscrew the screw (20) completely with the included screwdriver. With an open burner shaft the intensity of the flame cannot be adjusted by the air knob (2a) any longer. During the use of the turbo flame most of the needed air is taken inside through the open burner shaft. Remounting the cover of burner shaft (see paragraph 6.2).

## 8. Tilt adjustment

Insert the tilt adjustment (24) into the slots (21) at the bottom of the unit. The tilt adjustment can be used to the left or right side to protect the burning chamber from contamination when working with liquids.

## 9. Warranty

All WLD-TEC gas burners are covered under our two-year manufacturer warranty against any manufacture defects in material and workmanship. The WLD-TEC warranty guarantees all gasburners under normal usage conditions and does not cover any damages as a direct result of user misuse or/and abuse. The warranty is void upon any unauthorized servicing, disassembly or modifications.

## Technical data

Technology:	Microprocessor
<b>Programs</b>	
Foot pedal:	Standard (flame during pressed foot pedal) Start-Stop with timer, 60 min
Button:	Start-Stop with timer, 60 min
IR Sensor:	Start-Stop with timer, 60 min Auto-Off with learn function, 1 sec - 60 min

## Safety features

Safety Control System (SCS)  
with gas safety cut off:

ignition and flame control, temperature monitor with gas safety cut off  
burner head clogging and assembly monitor (BHC)  
foot pedal 4 h, Button & IR-Sensor 10 min  
indicates a hot burner head

Automatic unit switch off:  
Residual heat display:

## Gas supply and consumption

Gas supply:	1/4" left + filter
Gas types:	natural gas E/LL, 18 - 25 mbar, liquid gas 47.5 - 57.5 mbar, II <sub>2ELL3B/P</sub>
Connected load:	95 l/h natural gas (LL), 70 g/h liquid gas
Continuous cartridge operation: (approx.)	<i>Campingaz</i> : CV 360 - 40 min, C 206 - 170 min, CP 250 - 210 min, CV 470 - 370 min, CG 1750 - 150 min, <i>Guilbert</i> : Express 444 / 445 - 50 min <i>Coleman</i> : C 250 - 160 min, C 500 - 320 min

## Temperatures

Flame temperature:	1350 °C on liquid gas / 1300 °C on natural gas (E)
Temperature threshold level:	1 kW liquid gas, 1 kW natural gas

## Electrical

Power consumption:	2 VA
Power connection:	100 - 240V / 50/60Hz / max. 0.3 A; 9V DC / 1.3A Level 6 (stand by power consumption max. 0.1W)
IR-Sensor coverage:	35 - 70 mm, adjustable
DoubleClick IR-Sensor:	Time range 1 sec

## Mechanical

Casing and operating controls:	stainless steel / glass, UV and solvent resistant
Burner head:	removable and decomposable, stainless steel
Cover of the burner shaft:	Ø 23 mm, with drains
Measurements (w x h x d):	85 x 49 x 86 mm
Weight:	700 g

## Licences

DIN-DVGW Reg.-No.:	NG-2211AS0167
CE:	EN 61326-1, EN61010-1, EN61010-2-010
EU guidelines:	2014/30/EU, 2014/35/EU, 2011/65/EU

## Troubleshooting guide

- **The green LED's do not light up**

Check for correct connection and specification of the power adapter.  
Ensure that the original power adapter is used.

Specifications: 9 V / DC, 1.3 A

Polarity: 

---

- **The foot pedal does not work**

Check for correct connection of the foot pedal. Ensure that the foot pedal socket and plug is not twisted or broken.

---

- **No Flame**


In case of ignition or flame failure check if the burner head is clogged.  
Verify the input pressure of the used gas. Ensure that the correct nozzle is installed in the unit.

Nozzle N: natural gas, 18 - 25 mbar

Nozzle P, B: propane / butane gas, 47.5 - 57.5 mbar

---

- **LED “Burner head HOT / BHC“ flashes / Inspection of the burner head (clogging)**

Take care that there are no liquids or other substances at the Bypass (area between the inner and the outer ring). Contaminants can be removed, for example, with a brush. If there are contaminants in the area of the electrodes, the flame cannot encircle the electrodes correctly. For this reason pay particular attention to clogging in the areas marked in light colour and remove. The burner head can be cleaned with customary commercial disinfectants, or it can be sterilized in an autoclave or washed in a dishwasher.

---

- **In operating mode Pedal “StartStop“, SENSOR Auto-Off, SENSOR StartStop and „Button“ the flame burns 30 seconds, only**

BHC time limit is active, LED “Burner head hot / BHC“ is flashing.  
The burner head is clogged and must be cleaned.  
(see paragraph 4.2 & 6.1).

---

- **Flame too small / large / soft**

Check the position of the air and gas adjustment.  
Check if the correct nozzle is installed.

Nozzle N: natural gas, 18 - 25 mbar

Nozzle P, B: propane / butane gas, 47.5 - 57.5 mbar

Check if the drilling of the active nozzle is blocked (see paragraph 1.1). Unscrew the active nozzle. If the drilling is blocked clean with a brush or compressed air.

- **No ignition spark / LED Burner head "HOT / BHC" flashes but the burner head is clean**

Check if the ceramic electrodes are in good condition. In some cases the electrodes may break. To check move the ends of the electrodes. If they are not moving they should be okay. If they are moving more than 0.5 mm the electrodes are broken. The electrodes can be dismantled and changed by the user (**see paragraph 6.1**).

- **The burner shuts-off due to overtemperature frequently**

In case of overtemperature increase the air ventilation or change the operation site.

- **The IR-Sensor cannot be activated**

It is possible to adjust the sensor range to 0 mm or more than 70 mm. Then the IR-Sensor is out of range and cannot be activated. In that case, increase or decrease the IR-Sensor range at the function for adjusting the IR-Sensor detection range (**see paragraph 3.6**).

- **Green LED "Button", "Standard" or "Start-Stop" blinks 2x**

Ignition failure (**see paragraph 5.1**).

- **Green LED "Button", "Standard" or "Start-Stop" blinks 3x**

Flame failure (**see paragraph 5.2**).

- **Green LED "Button", "Standard" or "Start-Stop" blinks 4x**

Overtemperature (**see paragraph 5.3**).

- **Green LED "Button", "Standard" or "StartStop" blinks 5x**

This message indicates that the burner head is removed or not mounted correctly. Further operation is possible after the burner head is reinstalled (**see paragraph 5.4**).

- **Amber LED "Burner head HOT / BHC" (8) is on permanently**

Residual heat display is active. Attention: DO NOT TOUCH the burner head (**see paragraph 4.1**).

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