**Vidrio Industrial POBEL, S.A.**

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## 1.2. General information

- 1) Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please, notify to your supplier immediately.
- 2) Do not install or use the equipment without reading this handbook before.
- 3) This handbook must always be attached to the equipment and it must be available for all users.
- 4) **IMPORTANT! WE WILL NOT ACCEPT ANY APPARATUS TO BE REPAIRED IF IT IS NOT DULY CLEANED**
- 5) If any modification, elimination or lack of maintenance of any device of the equipment by the user transgress the directive 89/655/CEE, the manufacturer is not responsible for the damage that can occur.
- 6) Do not use the apparatus with liquids which can give off vapours capable of making explosive mixtures.

## 1.3. Packing list

The standard equipment consists of the following components:

Description	Qty.
Heater	1
Filter	1
Feed tube	1
Set of silicone tubes	1

## 1.4. Recommendations

- Electric:
  - Electrical wiring must be, at least, of 4.0 mm<sup>2</sup> section.
  - It must have ground.
  - Power point must let intensity of 13 A pass.
- Outlet:
  - The tap water should have screw mouth of 3/4".

### 1.5. Technical features

<i>Code</i>	<i>Length (cm)</i>	<i>Width (cm)</i>	<i>Height (cm)</i>	<i>Weight (Kg)</i>
720 720	39	58	17	5,2

<i>Code</i>	<i>Output (lts/h)</i>	<i>Conductivity at 20 °C*</i>	<i>Water consumption (l/min)</i>	<i>Consumption (W)</i>	<i>Voltage (V)</i>	<i>Intensity (A)</i>
720 720	4.2	1 $\mu$ S / cm	1,0	1 x 3000	230	13

*\*As distillates are open to the atmosphere before measurements can be made, the value of conductivity is really much less (down to 10 times) than it is really measured. Any ultrapure water in contact with the atmosphere has a conductivity of 1-2  $\mu$ S/cm and a pH of 5, due to the small amount of CO<sub>2</sub> (0.5 ppm) absorbed from the atmosphere.*


Fig. 1


Fig. 2

### 1.6. Assembly

Steps to follow:


1. Be sure that heater elements (washer of silicone rubber, guide fitting and PTPB nut) are in the order and position that the fig. 1 shows.

2. Insert the heater through the screw mouth (see fig. 2); hold the heater of the fitting (fig. 1,

piece 2) – and do not of the cap (fig. 1, piece 1) – and turn the PTPB nut (fig. 1, piece 3)

3. Be sure of putting the filter, of tightening and connecting the other end to the tap.

- Once you have assembled, check:
  - The silicone tubes hold suitable position to let water pass and to avoid water pressures.

### 1.7. Operating


Steps to follow:

1. Connect the end wiring of the heater to the socket, which is on the metallic base (piece 12).
2. Turn the tap on.
3. Turn the tap off when you see that the water level is on the heater.

### 1.8. Switch on

Steps to follow:

1. Switch on the green switch (fig. 2, piece 10). This will light showing that the still is connecting to the power point. The warning light will also light (fig. 2, piece 9) showing that the water supply is not enough.
2. Turn the tap on slowly until the pilot (fig. 2, piece 9) turn off.
  - Before few minutes, the heater will be getting red and the boiling will start. The stilled water goes out through the outlet (fig. 2, piece 2).
  - If water steam went out through steam hole (fig. 2, piece 3), you would have to increase water supply.

**1.9.Maintenance**

- You should clean the still regularly. So, in this way, the results will be better and you can get a still more durable.
- The cleaning periodicity depends on the salts which water supply contains.
- You should clean the filters of the tubes whenever you clean the still.

Steps to follow:

1. Put acetic acid, hydrochloric acid or whichever another product into the boiler through cleaning mouth (piece 5). The water should be hot.
2. You have to wait for the attack to the residues and empty the boiler through the same cleaning mouth (piece 5).
3. Rinse with clean water once or twice.
4. Fill with water until you see the water level is on the heater.

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