



designed for scientists



AS 130.1 Universal attachment

/// Hoja de Datos

Para el uso con diversos tipos de recipientes mediante rodillos de sujeción universales infinitamente variables.

Dimensiones (An x Pr x Al): 325 x 234 x 88 mm

Placa de montaje: 220 x 220 mm

Capacidad de recipientes: 15 x AS 2.1 (25 ml), 12 x AS 2.2 (50 ml), 6 x AS 2.3 (100 ml), 2 x AS 2.4 (250 ml), 2 x AS 2.5 (500 ml), 1 x AS 2.6 (1000 ml), 1 x AS 2.7 (2000 ml)

Peso: 850 g



designed for scientists

Información Técnica

Ancho [mm]	64
Profundidad [mm]	39.5
Peso [kg]	0.844
Temperatura ambiental permitida [°C]	5 - 40

