



designed for scientists



AS 260.1 Universal attachment

/// Hoja de Datos

Para diversos tipos de recipientes. Los rodillos de sujeción infinitamente variables permiten una adaptación universal a los recipientes.

Dimensiones (ancho x fondo x alto): 425 x 335 x 135 mm

Placa de sujeción: 320 x 320 mm

Capacidad del matraz: 28 x AS 2.1 (25 ml), 24 x AS 2.2 (50 ml), 15 x AS 2.3 (100 ml), 6 x AS 2.4 (250 ml), 6 x AS 2.5 (500 ml), 2 x AS 2.6 (1000 ml), 1 x AS 2.7 (2000 ml)

Peso: 1.600 g



designed for scientists

Los tornillos de apriete para la fijación del accesorio están incluidos en el volumen de suministro del agitador (HS / KS 260, KS 3000).





designed for scientists

Información Técnica

Dimensiones (An x Al x Pr) [mm]	425 x 135 x 335
Peso [kg]	1.618
Temperatura ambiental permitida [°C]	5 - 40

