



designed for scientists



AS 501.4 Soporte para pinzas de fijación

/// Hoja de Datos

Para agitar matraces, matraces Erlenmeyer y frascos con sección redonda (pinzas de fijación no incluidas en el suministro).

Dimensiones (An x Pr x Al): 475 x 460 x 95 mm

Capacidad del matraz: 55 x AS 2.1 (25 ml), 50 x AS 2.2 (50 ml), 30 x AS 2.3 (100 ml), 12 x AS 2.4 (250 ml), 11 x AS 2.5 (500 ml), 8 x AS 2.6 (1000 ml), 5 x AS 2.7 (2000 ml)

Peso: 2.640 g



designed for scientists

Información Técnica

Dimensiones (An x Al x Pr) [mm]	500 x 480 x 120
Peso [kg]	2.53
Temperatura ambiental permitida [°C]	5 - 40

