



designed for scientists



## H 100 Cubierta protectora

/// Hoja de Datos

For magnetic stirrers RET basic (36220XX) and RCT basic with production year before 2022 (33800XX + 38100XX).

Resistant to most acids, alkaline solutions and organic solvents.

Material: silicone

Max. temperature: 135 °C



designed for scientists

### Información Técnica

|   |                      |
|---|----------------------|
| Funcionamiento de temperatura max. [°C] | 200                  |
| Dimensiones (An x Al x Pr) [mm]         | 161.2 x 35.5 x 270.7 |
| Peso [kg]                               | 0.041                |
| Temperatura ambiental permitida [°C]    | 5 - 40               |

