

IKA

designed for scientists



H 102 Cover for RH basic

/// Hoja de Datos

Cubierta protectora para RH basic

www.ika.com

Sujeto a modificaciones técnicas

 IKAworldwide

 IKAworldwide /// #lookattheblue

 @IKAworldwide



designed for scientists

Información Técnica

| | |
|---|----------------------|
| Funcionamiento de temperatura max. [°C] | 200 |
| Dimensiones (An x Al x Pr) [mm] | 161.9 x 54.5 x 248.5 |
| Peso [kg] | 0.05 |
| Temperatura ambiental permitida [°C] | 5 - 40 |

