

IKA

designed for scientists



H 5000

/// Hoja de Datos

Beaker, acero inoxidable, 5 l,
Ø interior: 220 mm
Altura: 142 mm

www.ika.com

Sujeto a modificaciones técnicas



IKAworlwide



IKAworlwide /// #lookattheblue



@IKAworlwide



designed for scientists

Información Técnica

Diámetro interno [mm]	220.0
Forma de la parte inferior	cilíndrico hueco
Volumen utilizable max. [ml]	5000
Material	Acero inoxidable 1.4301
Alto [mm]	142.0
Peso [kg]	1.585
Temperatura ambiental permitida [°C]	5 - 40