

IKA

designed for scientists



IKAFLON 118 float Stirring drive (1 pc.)

/// Hoja de Datos

Floating magnetic stirring bar for I-MAG, 92 mm stirrer in 118 mm cage



designed for scientists

Información Técnica

Largo [mm]	NA
Diámetro [mm]	118
Forma	Redondo
Material	PPS
Material: Imán	NdFeB
Fuerza de atracción magnética	extrefuerte
Funcionamiento de temperatura max. [°C]	150
Peso [kg]	0.27
Temperatura ambiental permitida [°C]	5 - 150

