




designed for scientists



IKAFLON 140 slide (1 pc.)

/// Hoja de Datos

Magnetic stirring bar made of PTFE, with centric PEEK ring for reduction of friction and better concentricity with uneven vessel bottoms. For large vessels, turbulent mixing at low speeds.

Length: 140 mm, diameter: 27 mm, diameter with ring: 35 mm.

All materials are FDA compliant.



designed for scientists

Información Técnica

Largo [mm]	140
Diámetro [mm]	35
Forma	Redondo
Recubrimiento	PTFE
Material	PEEK + PTFE
Material: Imán	SmCo
Fuerza de atracción magnética	fuerte
Funcionamiento de temperatura max. [°C]	280
Peso [kg]	0.22
Temperatura ambiental permitida [°C]	5 - 280

