

Werksbescheinigung / Declaration of compliance with the order

gemäß / in accordance with DIN EN 10204 - 2.1

Mischgefäß mit Edelstahlkugeln / Stainless steel-Ball-Mill Tube:

Modell / Model:

BMT-20 S

Das Produkt wurde ausschließlich auf Eignung mit reinem **Acetonitril** getestet, dabei wurden keine Beeinträchtigungen festgestellt. Dieser Test gilt ausschließlich für den Einsatz der reinen Chemikalie, in wie weit Beimischungen anderer Substanzen negative Auswirkungen haben können, kann hier weder getestet noch ausgeschlossen werden und unterliegt damit der Verantwortung des Kunden. /

This product was only tested with pure **Acetonitril**. No impairment or damage to the product was detected. This test is valid only for the pure chemical. Admixture of other substances could have negative effects. We are not able to test the effects of mixtures between Acetonitril and other substances on this product. Responsibility for this is completely on customers side.

Bedingungen / Conditions:

Die Tubes BMT-20 S sind ausschließlich für den Einmalgebrauch ausgelegt. /
The tubes BMT-20 S are appropriate only for a single-use.

max. Anwendungszeit: 29 Minuten / 4000 U/min.
max. Application time: 29 minutes / 4000 rpm

Staufen, 03.02.2009

Andreas Frei
Corporate Director
Research & Development lab/ana

Jürgen Goebel
Normenstelle /
Standard Department

Dieses Dokument wurde elektronisch erstellt und ist ohne Unterschrift gültig. (DIN EN 10204, Abschnitt 5) /
This document was created electronically and is valid without signature. (DIN EN 10204, Chapter 5)



Designed
to work perfectly

Germany

IKA®-Werke
GmbH & Co.KG
Janke & Kunkel-Str. 10
79219 Staufen
Germany
Tel. +49 7633 831-0
Fax +49 7633 831-98
sales@ika.de

USA

IKA® Works, Inc.
2635 Northchase
Pkwy SE
Wilmington
NC 28405-7419,
USA
Tel. +1 910 452-7059
Fax +1 910 452-7693
usa@ika.net

China

IKA® Works Guangzhou
173-175 Friendship Road
Guangzhou
Economic & Technological
Development District
510730 Guangzhou, China
Tel. +86 20 8222-6771
Fax +86 20 8222-6776
sales@ikagz.com.cn

India

IKA® India Pvt. Ltd
814/475, Survey # 129/1
Mysore Road, Kengeri
Bangalore - 560 060,
Karnataka,
India
Tel. +91-80-26253900
Fax +91 80 26253901
info@ika.in

Asia / Australia

IKA® Works (Asia) Sdn Bhd
No. 17 & 19, Jalan PJU 3/50
Sunway Damansara
Technology Park
47810 Petaling Jaya
Selangor, Malaysia
Tel. +60 3 7804-3322
Fax +60 3 7804-8940
sales@ika.com.my