


IKA

designed for scientists



C 6030 Venting station

/// Hoja de Datos

To vent the combustion gases in a controlled manner into an absorption solution for further analysis.

www.ika.com

Sujeto a modificaciones técnicas

 IKAworldwide

 IKAworldwide /// #lookattheblue

 @IKAworldwide



designed for scientists

Información Técnica

Dimensiones (An x Al x Pr) [mm]	150 x 250 x 250
Peso [kg]	1.883
Temperatura ambiental permitida [°C]	5 - 40

