

IKA

designed for scientists



DB 1.4

/// Hoja de Datos

Utilizado para microtubos VWR/Eppendorf (2,0 ml)

Tamaño de poro: Ø 10,9 mm

N.º de orificios: 20

Profundidad de orificio: 38,1 mm

Dimensiones (An x Pr x Al): 95 x 76 x 51 mm

www.ika.com

Sujeto a modificaciones técnicas

 IKAworldwide

 IKAworldwide /// #lookattheblue

 @IKAworldwide



designed for scientists

Información Técnica

Dimensiones (An x Al x Pr) [mm]	95 x 51 x 76
Peso [kg]	0.823
Temperatura ambiental permitida [°C]	5 - 40

