

# IKA

designed for scientists



## DB 4.4

/// Hoja de Datos

Utilizado para tubos de fondo redondo (12/13 mm)

Tamaño de poro: Ø 13,9 mm

N.º de orificios: 20

Profundidad de orificio: 48,4 mm

Dimensiones (An x Pr x Al): 95 x 76 x 51 mm

[www.ika.com](http://www.ika.com)

Sujeto a modificaciones técnicas

 IKAworldwide

 IKAworldwide /// #lookattheblue

 @IKAworldwide



designed for scientists

### Información Técnica

Dimensiones (An x Al x Pr) [mm]	95 x 51 x 76
Peso [kg]	0.62
Temperatura ambiental permitida [°C]	5 - 40

