

IKA

designed for scientists



DB 4.8

/// Hoja de Datos

Utilizado para tubos de fondo redondo (35 mm)

Tamaño de poro: Ø 35,0 mm

N.º de orificios: 4

Profundidad de orificio: 47,6 mm

Dimensiones (An x Pr x Al): 95 x 76 x 51 mm

www.ika.com

Sujeto a modificaciones técnicas



IKAworlwide



IKAworlwide /// #lookattheblue



@IKAworlwide



designed for scientists

Información Técnica

Dimensiones (An x Al x Pr) [mm]	95 x 51 x 76
Peso [kg]	0.547
Temperatura ambiental permitida [°C]	5 - 40

