

# IKA

designed for scientists



## DB 5.4

/// Hoja de Datos

Utilizado para frascos de 19 mm

Tamaño de poro: Ø 19,7 mm

N.º de orificios: 12

Profundidad de orificio: 45,0 mm

Dimensiones (An x Pr x Al): 95 x 76 x 51 mm

[www.ika.com](http://www.ika.com)

Sujeto a modificaciones técnicas

 IKAworldwide

 IKAworldwide /// #lookattheblue

 @IKAworldwide



designed for scientists

### Información Técnica

|                                      |              |
|--------------------------------------|--------------|
| Dimensiones (An x Al x Pr) [mm]      | 95 x 51 x 76 |
| Peso [kg]                            | 0.525        |
| Temperatura ambiental permitida [°C] | 5 - 40       |

