

# IKA

designed for scientists



## DB 8.1

/// Hoja de Datos

Utilizado para cubetas 2 x 6 (12,5 mm)

Profundidad de orificio: 25,4 mm

Dimensiones (An x Pr x Al): 95 x 76 x 51 mm

[www.ika.com](http://www.ika.com)

Sujeto a modificaciones técnicas



IKAworlwide



IKAworlwide /// #lookattheblue



@IKAworlwide



designed for scientists

### Información Técnica

Dimensiones (An x Al x Pr) [mm]	95 x 51 x 76
Peso [kg]	0.8332
Temperatura ambiental permitida [°C]	5 - 40

