



designed for scientists



DIS-300-S-M.10 Tubo de disolución, 300ml, 10 unidades

/// Hoja de Datos

Volume: up to 300 ml

Application: Crushing and dissolution of tablets and dragées. This tube is equipped with a dissolver blade and opens up new application fields. It was especially designed for dissolving samples in different solvents.



designed for scientists

Información Técnica

Rango de volumen (H2O) [l]	.1 - .3
Diámetro del rotor [mm]	35.5
Diámetro del tubo [mm]	77
Velocidad circunferencial a 6000 rpm [m/s]	11.2
Volumen del tubo [ml]	300
Conformidad con la FDA	sí
Material: Membrana	PP
Material: Recipiente de pruebas y tapón	PP
Material: Varilla de giro	PEEK+PTFE+grafito+fibra de carbono
Material: Disolvente	Acero inoxidable
Peso [kg]	0.66
Temperatura ambiental permitida [°C]	5 - 40

