



designed for scientists



H.PUR.12.AD Tube set incl. M 16x1 connection

/// Hoja de Datos

Hose set DN 12 (2 x 1.5 m), incl. 4 hose clamps, DN 12 olives and M 16x1 lock nuts for connecting a double walled reactor vessel.

Material: PUR

Temperature range: -30 - 90 °C

Max. operating pressure: 3 bar

Connection: for M 16x1 thread (male)

inner: 12 mm; outer: 16 mm

Chemical Resistance: water; silicon oils; (ethylene glycol); mineral oil; ethanol

IKA

designed for scientists

