

# IKA

designed for scientists



## IB HF 24 eco

/// Hoja de Datos

Baño de plástico, forma alta, tamaño L, volumen 24 l

Temperatura máx.: 100 °C (solo agua)

Dimensiones interiores: 490 x 286 x 200 mm

Dimensiones exteriores: 584 x 338 x 205 mm

[www.ika.com](http://www.ika.com)

Sujeto a modificaciones técnicas

 IKAworldwide

 IKAworldwide /// #lookattheblue

 @IKAworldwide



designed for scientists

### Información Técnica

Peso [kg]	2.8828
Temperatura ambiental permitida [°C]	5 - 40

