



designed for scientists



IKA Tip 1000 µl bag (1.000 pcs)

/// Hoja de Datos

Universal fit pipette tips, non-sterile.

Suitable for all blue IKA pipettes.

Size: 1000 µl

- Autoclavable material (polypropylene)
- Excellent volume accuracy due to precise tip orifice
- Minimal residual wetting due to controlled quality of raw materials used
- Guaranteed dust-free due to special ionization treatment

www.ika.com

Sujeto a modificaciones técnicas



IKAworlwide



IKAworlwide /// #lookattheblue



@IKAworlwide



designed for scientists

Packaging unit: 1000 pcs/bag



IKA

designed for scientists

