



designed for scientists



IKA Tip 1000 µl box (600 pcs.)

/// Hoja de Datos

Universal fit pipette tips, 1000 µl size, non-sterile.
Fits to all blue IKA pipettes.

Made of polypropylene and can be autoclaved.
Excellent volume accuracy due to precise tip orifice.
Minimal residual wetting due to controlled quality of raw materials used.
Guaranteed dust-free due to special ionisation treatment.

www.ika.com

Sujeto a modificaciones técnicas



IKAworlwide



IKAworlwide /// #lookattheblue



@IKAworlwide



designed for scientists

60 pcs/box, 10 boxes included in scope of delivery.





designed for scientists

