



designed for scientists



## IKA Tip 200 µl bag (1.000 pcs)

/// Hoja de Datos

Universal fit pipette tips, non-sterile.

Suitable for all yellow IKA pipettes.

Size: 200 µl

- Autoclavable material (polypropylene)
- Excellent volume accuracy due to precise tip orifice
- Minimal residual wetting due to controlled quality of raw materials used
- Guaranteed dust-free due to special ionization treatment



designed for scientists

\*\* Recommended for use with IKA single channel pipettes \*\*

Packaging unit: 1000 pcs/bag



# IKA

designed for scientists

