

# IKA

designed for scientists



## IO T 1.00

/// Hoja de Datos

Bandeja de rejilla de alambre, cromada

Dimensiones (L x Pr):

548 mm x 421 mm

Capacidad máxima de carga: 30 kg

[www.ika.com](http://www.ika.com)

Sujeto a modificaciones técnicas



IKAworlwide



IKAworlwide /// #lookattheblue



@IKAworlwide



designed for scientists

### Información Técnica

Dimensiones (An x Al x Pr) [mm]	546 x 34 x 410
Peso [kg]	2.1

