



designed for scientists



## LR 2000.3 Recipiente de reacción

/// Hoja de Datos

Double-walled reactor vessel including accessories for the laboratory reactors LR-2.ST.

Volume: 2.0 l

Material in contact with medium: stainless steel 1.4571.

Including quick connection adapters with DN 8 Olive for connection of temperature hoses.

Scope of delivery:

Stainless steel reactor incl. quick connection adapters with DN 8 olive.



designed for scientists

### Información Técnica

Material en contacto con producto	AISI 316L / 1.4571
Volumen [l]	2
Diámetro interior del recipiente [mm]	153
Diámetro exterior del recipiente [mm]	194
Altura (dimensión de la sala) [mm]	230
Anchura (dimensión de la sala) [mm]	194
Profundidad (dimensión de la sala) [mm]	225
Dimensiones de la brida	DN150
Forma del fondo	Cabezal torisférico
Manta calefactora	sí
Diámetro interior de la entrada de la camisa de atemperación [mm]	10
Aplicación de vacío	sí
Conexión interior de la salida de la camisa de atemperación [mm]	10
Peso [kg]	3.3