




designed for scientists



## LT 5.24 Adaptador de tubo flexible

/// Hoja de Datos

Adaptador de manguera

Rango de temperatura: -60 a 180 °C

Longitud: R 1/8" x M 16 x 1

Material: metal



designed for scientists

### Información Técnica

|                                      |              |
|--------------------------------------|--------------|
| Dimensiones (An x Al x Pr) [mm]      | 8 x 86 x 134 |
| Peso [kg]                            | 0.085        |
| Temperatura ambiental permitida [°C] | 5 - 40       |

