

IKA

designed for scientists



MultiDrive TC 1 Protective housing

/// Hoja de Datos

Carcasa protectora para MultiDrive MT 150 y BT 250

www.ika.com

Sujeto a modificaciones técnicas



IKAworlwide



IKAworlwide /// #lookattheblue



@IKAworlwide



designed for scientists

Información Técnica

Apto para lavavajillas	sí
Apropiado para solventes	sí
Dimensiones (An x Al x Pr) [mm]	110 x 180 x 128
Peso [kg]	0.365
Temperatura ambiental permitida [°C]	5 - 40
Humedad relativa permitida [%]	80