

# IKA

designed for scientists



## RV 10.606 Distilling spider with 5 flasks 50 ml (NS 24/29)

/// Hoja de Datos

Para destilación simultánea en 5 matraces de evaporación, 50 ml.

[www.ika.com](http://www.ika.com)

Sujeto a modificaciones técnicas

 IKAworldwide

 IKAworldwide /// #lookattheblue

 @IKAworldwide



designed for scientists

### Información Técnica

Peso [kg]	0.4588
Temperatura ambiental permitida [°C]	5 - 40

