




designed for scientists



SY 100.2 Reactor lid

/// Hoja de Datos

Reactor lid for the EasySyn reactor system.

Flange size: DN 100.

Suitable for reactor vessels with volumes of 500 ml (SY 500 and SY 500 D) and 1000 ml (SY 1000 and SY 1000 D).

Openings: 4 x NS 29/32, 1 x NS 14/23, 2 x GL 18



designed for scientists

Información Técnica

Altura (dimensión de la sala) [mm]	95
Diámetro externo [mm]	138
Tipo de brida	acc. DIN 12214
Dimensiones de la brida	DN100
Aplicación de vacío	sí
Conexión central de la tapa	NS 29/32
Conexión laterales de la tapa	3x NS 29/32, 1x NS 14/23, 2x GL18
Material de la tapa del reactor	Vidrio de borosilicato 3,3 DIN ISO 3585
Funcionamiento de temperatura [°C]	-50 - 230
Presión Max.	atmospheric pressure
Vacío alcanzable [mbar]	3