



VOLS-1 flex

/// Hoja de Datos

The VOLS-1 flex adapter set for low volumes consists of

- Mounting VOLS 1.1
- Chamber cover VOLS 1.21
- Extension connector VOLS 1.11
- Water jacket VOLS 1.2

Compared to the standard VOLS-1 adapter set, VOLS-1 flex does not include a standard chamber (VOL-C-RTD-1) and spindle (VOL-SP-6.7).



designed for scientists

VOLS-1 flex can thus be flexibly combined with the wide range of spindles and chambers. VOLS-1 flex is particularly suitable for applications with the following spindle/chamber combinations:

- VOL-SP-2.1 with VOL-C-RTD-2.1
- VOL-SP-3.8 with VOL-C-RTD-3.8
- VOL-SP-4.2 with VOL-C-RTD-4.2

- VOL-SP-7.1 with VOL-C-RTD-1
- VOL-SP-9.0 with VOL-C-RTD-1
- VOL-SP-9.4 with VOL-C-RTD-1
- VOL-SP-10.4 with VOL-C-RTD-1
- VOL-SP-11.0 with VOL-C-RTD-1
- VOL-SP-13.5 with VOL-C-RTD-1
- VOL-SP-16.1 with VOL-C-RTD-1

The digits in the VOL-SP spindle name represent the usable volume. Spindles and chambers are available separately as accessories.



designed for scientists

