


THERMOCELL PACKAGE



Immersion PT100 probe
(Ref 000696)

FEATURES AND BENEFITS

- Temperature range from Room T°C to +300°C.
- Economic electric oven.
- Compatible with bayonet coupling device.
- Compatible with norms ASTM D4402 and D3236.

OPTIONS & ACCESSORIES

- Immersion PT100 probe with holder (Ref. 000696)
- Lid for cup MB-SV13R (Ref. 130185)

WHAT'S INCLUDED WITH THERMOCELL (Ref. N130150)

- A - Choice of one MK-SV spindle (p 83).
 - B - Cup MB-SV13R (Ref. 116213).
 - C - Holder H-SV13RD for disposable cup MB-SV13RD. (Ref. 116313).
 - D - 5 Disposable cups MB-SV13RD (Ref. 116413 for 100pcs).
 - E - Centring tool ST-R (Ref. 114436).
- RT-1 PLUS oven (Ref. N130130).

