

**Packed weight:** 11kg  
**Product weight:** 9.2kg

**Packed dimensions;**  
A 470mm  
B 470mm  
C 450mm

**Product dimensions;**  
A 360mm  
B 380mm  
C 275mm

<b>D</b>	298mm (at opening)
<b>D</b>	235mm (tray length)
<b>E</b>	323mm (at opening)
<b>E</b>	240mm (tray width)
<b>F</b>	118mm
<b>G</b>	136mm
<b>H</b>	227mm


## Universal tray (TU12)

Vessel	Capacity (bath w/out lid)	Capacity (bath with lid)
25ml flask	24	24
50ml flask	18	18
100ml flask	10	10
150ml flask	9	9
250ml flask	6	6
500ml flask	4	4
1000ml flask	2	0
125ml Duran® YOUTILITY	15	15
250ml Duran® YOUTILITY	9	9
500ml Duran® YOUTILITY	6	3
1000ml Duran® YOUTILITY	5	0
25ml Duran® 21-801	25	25
50ml Duran® 21-801	23	23
100ml Duran® 21-801	14	14
150ml Duran® 21-801	12	12
250ml Duran® 21-801	9	9
500ml Duran® 21-801	6	5
1000ml Duran® 21-801	5	2


## Flask tray (TF12)

Vessel	Capacity (bath w/out lid)	Capacity (bath with lid)
25ml flask	18	18
50ml flask	14	14
100ml flask	9	9
250ml flask	5	5
500ml flask	4	4
1000ml flask	2	2
SH-DWP	2	2


## Flask tray (TF12) - Optimal flask clamp positions




25ml




50ml




100ml



250ml



500ml



1000ml