

Packed weight: 13.5kg
Product weight: 11.2kg

Packed dimensions;
A 450mm
B 730mm
C 450mm

Product dimensions;
A 335mm
B 565mm
C 275mm

| | |
|----------|---------------------|
| D | 500mm (at opening) |
| D | 420mm (tray length) |
| E | 298mm (at opening) |
| E | 235mm (tray width) |
| F | 118mm |
| G | 136mm |
| H | 227mm |


Universal tray (TU18)

| Vessel | Capacity (bath w/out lid) | Capacity (bath with lid) |
|-------------------------|------------------------------|-----------------------------|
| 25ml flask | 44 | 44 |
| 50ml flask | 35 | 35 |
| 100ml flask | 18 | 18 |
| 150ml flask | 15 | 15 |
| 250ml flask | 11 | 11 |
| 500ml flask | 7 | 5 |
| 1000ml flask | 4 | 0 |
| 125ml Duran® YOUTILITY | 24 | 24 |
| 250ml Duran® YOUTILITY | 15 | 15 |
| 500ml Duran® YOUTILITY | 10 | 5 |
| 1000ml Duran® YOUTILITY | 8 | 0 |
| 25ml Duran® 21-801 | 45 | 45 |
| 50ml Duran® 21-801 | 28 | 28 |
| 100ml Duran® 21-801 | 18 | 18 |
| 150ml Duran® 21-801 | 15 | 15 |
| 250ml Duran® 21-801 | 15 | 15 |
| 500ml Duran® 21-801 | 8 | 8 |
| 1000ml Duran® 21-801 | 6 | 0 |


Flask tray (TF18)

| Vessel | Capacity (bath w/out lid) | Capacity (bath with lid) |
|--------------|------------------------------|-----------------------------|
| 25ml flask | 33 | 33 |
| 50ml flask | 26 | 26 |
| 100ml flask | 17 | 17 |
| 250ml flask | 14 | 14 |
| 500ml flask | 6 | 6 |
| 1000ml flask | 4 | 4 |
| SH-DWP | 4 | 4 |


Flask tray (TF18) - typical vessel configuration




25ml




50ml




100ml



250ml



500ml



1000ml