


Overall dimensions




Please note images are for representation purposes only

Packed weight: 36.3kg
Product weight: 29.1kg

Packed dimensions:
A: 465mm
B: 610mm
C: 755mm

Product dimensions:
A: 245mm
B: 430mm
C: 640mm

Working area



Working area:
D: 118mm
E: 154mm
F: 145mm
(10mm below top)

LT ecocool™ 100 minimum immersion volume depth: 85mm

LT ecocool™ 150 minimum immersion volume depth: 130mm

Please note images are for representation purposes only

Vessel capacity without lid

| | |
|-------------------------|----------------|
| Duran 21-801: 50ml | 9 ⁺ |
| Duran 21-801: 100ml | 4 |
| Duran 21-801: 150ml | 3 |
| Duran 21-801: 250ml | 2 |
| Duran 21-801: 500ml | 1 |
| Duran 21-801: 750ml | 1 |
| Duran 21-801: 1000ml | 1 |
| Duran YOUTILITY: 125ml | 4 |
| Duran YOUTILITY: 250ml | 2 |
| Duran YOUTILITY: 500ml | 2 |
| Duran YOUTILITY: 1000ml | 1 |

Vessel capacity with tank lid

| | |
|--------------------------|----|
| Erlenmeyer flask: 100ml | 4 |
| Erlenmeyer flask: 250ml | 2 |
| Erlenmeyer flask: 500ml | 1 |
| Erlenmeyer flask: 1000ml | 1* |

* Bridge plate released

+ May need height support in the tank depending on fill level

