

PP-20/3, platform

DESCRIPTION

Flat platform with non-slip rubber mat can accommodate various low profile containers.

SPECIFICATIONS

Working area	355 × 455 mm × 3x
Height between platforms	140 mm
Overall dimensions (W×D×H)	380 × 480 × 300 mm



CAT. NUMBER

BS-010126-DK

Platform PP-20-3



PSU-20i
Without platform
Multi-functional Orbital Shaker

Shaker PSU-20i provides three motion types: orbital, reciprocal and vibrating, which can be performed separately, pairwise and sequentially in repeated cycles. Shaker is designed for applications both in small specialized ...

[read more](#)