

Clifton Deep Chamber baths NE4-25D Stirred Gerber Digital Bath

Clifton Range®



This Deep Chamber Bath™ is suitable for Gerber testing or where deeper water volumes are required.

The stainless steel deep tank contained within a sturdy outer case and supplied complete with a stainless steel shelf all providing good water circulation and temperature stability. Ideal for applications such as Gerber milk testing.



NE4-25D

GENERAL	
Temp Range	Ambient +5°C - 99°C
Temp Alarm	Under/over temperature alarm is automatically set 4°C above and 4°C below set point
Temp Control	PID Digital, LED display, resolution 0.1°C
Sensitivity/Uniformity	±0.1°C / ±0.1°C
Timer	Visual indicator and end message - Variable 00:01 to 99:59 hh:mm
Heater Power (220-240V)	1250W
Voltage	230V 50-60Hz
Drain	Drain tap fitted
DIMENSIONS	
Tank Capacity	25 Litre
External (mm)	440w x 265d x 370h mm
Internal (mm)	262w x 235d x 330h mm
Water Level max	280 mm

Accessories

Product Code	Description
6880	12 Hole stainless steel Butyrometer Rack
CAL	Calibration Record