

**Clifton Deep Chamber baths™
NE1-25 Unstirred
Thermostatic Bath**

Clifton Range®



Unstirred Clifton DeepChamberBath™, with Thermostatic analogue control provides a deep volume of temperature controlled water. Comprising of a stainless steel tank housed in a durable paint finished outer case. Painted surfaces feature a tough anti-bacterial finish with good chemical resistance. Complete with a removable stainless steel perforated shelf.



NE1-25

GENERAL	
Temp Range	Ambient +5°C - 99°C
Temp Control	Analogue dial pointer, graduated in 5°C increments
Sensitivity/Uniformity	±0.25°C / ±0.2°C
Heater Power (220-240V)	1000W
Voltage	230V 50-60Hz
DIMENSIONS	
External (mm)	345w x 345d x 365h mm
Internal (mm)	300w x 300d x 250h mm

Accessories



Polypropylene Spheres (Pack of 200)

Bath	NE1-25
No. of packs	4 Packs
Accessory Code	BP0368