



Clifton digital ShallowBaths™ are perfect for heating small height samples such as micro tubes, micro plates and bottles. ShallowBaths™ require less water and less process chemicals, so are therefore economical to use. Controls feature auto alarm setting, illuminated on/off switch, heater and over temperature indications.

Comprising of a high quality corrosion resistant stainless steel tank, with a bright clean finish housed in a durable painted outer case. Painted surfaces feature anti-bacterial finish with a good chemical resistance.

	NE2-2D	NE2-5D	NE2-9D
GENERAL			
Chamber Capacity	2 Litres	5 Litres	9 Litres
Temperature Range	Ambient + 5°C - 99°C		
Temperature Control	PID Digital, LED display showing actual or set temperature, resolution 0.1°C		
Temperature Alarm	Over temperature alarm is automatically set 4°C above set point		
Sensitivity/Uniformity	±0.2°C / ±0.1°C		
Heater Power (220-240V)	200W	800W	
Voltage	120V or 230V models available 50-60Hz		
DIMENSIONS			
External (mm)	332w x 185d x 135h	361w x 332d x 135h	535w x 332d x 135h
Internal (mm)	300w x 150d x 65h	325w x 300d x 65h	500w x 300d x 65h
Water Level Max. (mm)	40		
Water Level Min. (mm)	20		

ShallowBaths™

Accessories



Stainless Steel Lift Off Gable Lid

Bath	NE2-2D	NE2-5D	NE2-9D
Accessory Code	SL1-4	SL1-14	SL1-22



Stainless Steel Flat Lid

Bath	NE2-2D	NE2-5D	NE2-9D
Accessory Code	LD-4	LD-14	LD-22



Clear Lift Off Gable Lid

Bath	NE2-2D	NE2-5D	NE2-9D
Accessory Code	GL1-4AF	GL1-14AF	GL1-22AF



Polypropylene Spheres (Pack of 200)

Bath	NE2-2D	NE2-5D	NE2-9D
No. of Packs	1 Pack	2 Packs	
Accessory Code	BP0368		



Stainless Steel Test Tube Racks

Bath	NE2-2D	NE2-5D	NE2-9D
No. of Racks	1	4	6
Accessory Code	6870 : 26 holes x 17mm DIA		
	6871 : 16 holes x 26mm DIA		
	6872 : 36 holes x 13mm DIA		
	6873 : 18 holes x 19mm DIA		
	6900 : 12 holes x 32mm DIA		

When using the racks on a 2 litre chamber, remove 1 shelf from the rack.



Bath Armor™

Bath	NE2-2D	NE2-5D	NE2-9D
No. of packs	1 x BA-5	1 x BA-5	1 x BA-12
Accessory Code	BA-5 : 4 Litre Pack		
	BA-12 : 12 Litre Pack		

Universal Accessories

CAL	Calibration Record
-----	--------------------