

# KERN 313-07

# KERN

Weight set, E2, 1 mg - 2 kg, Knob, stainless steel polished (OIML), in wooden case



Article number for guaranteed plus-deviation 921-306

Pictograms

STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML E2
Product family	313

## Measuring System

Test weight - nominal value	1 mg - 2 kg
Weight OIML class	E2

## Construction

Shape weight	Knob
Dimension (W×D×H)	305×190×140 mm
Number of weights	27
Case of weight set	in wooden case
Material	stainless steel polished (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	305×190×140 mm
Net weight	9,01 kg
Shipping method	Parcel service
Net weight approx.	10 kg
Gross weight approx.	10 kg
Shipping weight	9,01 kg

## Services (optional)

Article number for conformity assessment	952-306
Article number for DAkkS calibration	962-306