

KERN 313-09

KERN

Weight set, E2, 1 mg - 10 kg, Knob, stainless steel polished (OIML), in wooden case



Article number for guaranteed plus-deviation 921-308

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML E2
Product family	313

Measuring System

Test weight - nominal value	1 mg - 10 kg
Weight OIML class	E2

Construction

Shape weight	Knob
Dimension (W×D×H)	540×200×240 mm
Number of weights	29
Case of weight set	in wooden case
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	2 d
Dimensions packaging (W×D×H)	540×200×240 mm
Net weight	30,61 kg
Shipping method	Freight
Net weight approx.	32 kg
Gross weight approx.	34 kg
Shipping weight	32,5 kg

Services (optional)

Article number for conformity assessment	952-308
Article number for DAkkS calibration	962-308