



# Optional accessories & functions for Incubator/Drying Oven/Hot Air Sterilizer

## Shelves



## Specifications without shelves

Number of the shelves of standard accessories can be reduced upon request. Please let us know your preference upon estimating.

◆ Applicable: all models

## Load capacity of shelves


If you need to load more than the standard load capacity of shelves, we can produce the shelves of special specification which satisfy required load capacity. We propose specifications that match the weight and shape of the sample.

◆ Applicable: all models

## Recorder + floating sensor

A floating sensor can be inserted into the chamber from the standard measurement port to measure and record the temperature inside the chamber. As a incubators in pharmaceutical manufacturing and hot air sterilizers as GMP management device, it is a necessary function for daily use record. Both paper type and paperless type can be selected.

◆ Applicable: all models




| Applicable Model   |
|--|
| RT-120HM, RY-120HM, ST-120HM, SY-120HM, PW-120HM                                 |
| RT-120HL, RY-120HL, ST-120HL, SY-120HL, PW-120HL                                 |
| ITD-100, ILD-110HG, ILD-100HGC, RT-120HG, RY-120HG, ST-120HG, SY-120HG, PW-120HG |
| ITD-120, ILD-120, ILD-120C   |
| ITD-140, ITD-160, ILD-150, ILD-150C  |
| GTR-100  |
| GTR-120  |
| GTR-150, GTR-200   |
| SB-120HM/HL  |
| SB-120HG   |
| DCH-120HM  |
| DCH-120HL/HG   |
| KMD-100HM  |
| KMD-120HL  |
| KMD-100, KMD-100F, KMD-120, KMD-120F   |

## Indicator light

An indicator light can be mounted on the top of the device so that you can check the operating status of the device from a distance.

◆ Applicable: GTR series  
In case of other models, consultation are required



The optional accessories above are available both when you place an order for the machine and after the order.

Please note that the optional functions above can be fitted to the machine only when you place an order and cannot be retrofitted to the machine after the order.